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Using new instruments for the first time

Please check products for possible damage in transit. All instruments must be cleaned, disinfected and sterilised **before use**.

Field of application

The products must be used exclusively for the intended purpose by appropriately trained and qualified personnel. The doctor concerned is responsible for selecting the instrumentation for specific applications, for appropriate training and for adequately informing the operating personnel, as well as for operational use.

Checks

Instruments must be checked to ensure that they are in working order each time they are used. For reliable application, the functionality of the product must be particularly carefully checked in the areas of the blades, points, locks, latches, safety catches and all moving parts with regard to cracks, breakages and other possible damage. Discard damaged products immediately.

Repair

The guarantee and warranty will become void if products are repaired by companies or persons who have not been approved by GEOMED Medizin-Technik,

Disposal

If instruments can no longer be repaired and conditioned, then they should be sent to the usual hospital disposal facility.

Materials

The materials used are stainless steels to DIN/EN/ISO 7153-1, pure titanium or titanium alloys to DIN/ISO 5832-2.

Please also refer to the recommendations of the Instrument Preparation Working Group <http://www.a-k-i.org>.

Reconditioning instructions

	<p>This processing instruction was taking into account under the requirements of EN ISO 17664 / ANSI AAMI ST 81 sterilization of medical devices - information to be provided by the manufacturer for the reprocessing of medical devices.</p> <p>Instruments are supplied unsterilised and must be cleaned and if necessary disinfected and sterilised before use.</p> <p>Instruments may only be prepared by personnel who have the necessary specialist knowledge and training, and are able to assess the risks involved together with the associated effects.</p> <p>Additional aspects must be taken into account for certain instruments; for more details, please see section entitled "Special information".</p>
General preparation principles	<p>All instruments must be cleaned, disinfected and sterilised each time they are used; this also applies particularly to the first time they are used after delivery, as all instruments are supplied unsterilised (cleaning and disinfection after removal of protective transit packaging; sterilisation after unpacking). Effective cleaning and disinfection is an indispensable prerequisite for effective sterilisation.</p> <p>When using instruments, please ensure that, as part of your responsibility for the sterility of instruments, only methods for cleaning/disinfection and sterilisation, which have been adequately validated for the equipment and products, are implemented, that the equipment used (disinfector, steriliser) is regularly maintained and tested, and that the validated parameters are maintained for each cycle.</p> <p>Please also observe the legal regulations applicable in your country as well as the hygiene regulations of the medical practice or hospital. This applies particularly to the different requirements relating to effective prion inactivation.</p>

File name / Revision	Produced by:	Date produced:	Date modified:	Page
GA-010_11_Reusable_surgical_instruments	ai / hh	25.01.12	08.02.2017	1 of 5

Instructions for use

Reusable surgical instruments

Cleaning and disinfection - Basic principles	<p>An automated process (disinfector) should be used wherever possible for cleaning and disinfecting instruments. Due to its significantly lower effectiveness and reproducibility, a manual process - even using an ultrasonic bath - should only be used if an automated process is not available.</p> <p>The pre-treatment must be carried out in both cases.</p>
Pre-treatment	<p>Coarse contamination must be removed from the product directly after use (within a maximum of 2 hrs).</p> <p>For this purpose, use running water or a disinfectant solution; the disinfectant solution should be aldehyde-free (otherwise solidification of blood contamination may occur), have proven effectiveness (e.g. DGHM or FDA approval or CE marking), be suitable for the disinfection of instruments, and be compatible with the instruments (see section "Material compatibility").</p> <p>For the manual removal of contamination, use only a soft brush or a clean soft cloth, which is used only for this purpose, and never metal brushes or steel wool.</p> <p>Please note that disinfectants used for pre-treatment serve only to protect personnel, and cannot replace the later disinfection steps, which are to be carried out when cleaning is complete.</p> <p>If instruments can be taken apart, then take them apart before preparation.</p>
Manual cleaning and disinfection	<p>When selecting the cleaning agents and disinfectants to be used, it must be ensured</p> <ul style="list-style-type: none"> • that these are basically suitable for the cleaning and disinfection of instruments, • that the cleaning agent - if applicable - is suitable for ultrasonic cleaning (no foaming), • that a disinfectant with proven effectiveness is used (e.g. VAH/DGHM or FDA approval or CE marking) and that this is compatible with the cleaning agent used, and • that the chemicals used are compatible with the instruments (see section "Material compatibility"). <p>Combined cleaning agents/disinfectants should only be used in the case of extremely low prior contamination (no visible contamination) of the instruments.</p> <p>It is essential that the concentrations and reaction times specified by the manufacturer of the cleaning agents and disinfectants be observed. Use only freshly prepared solutions, sterile water or water with low germ content (max. 10 germs/ml) or water with low endotoxin content (max. 0.25 endotoxin units/ml) (e.g. purified water (PW)/highly purified water (HPW)). For drying, use only clean, lint-free cloths, which are used only for this purpose, or oil-free, filtered air.</p> <p>Cleaning</p> <ul style="list-style-type: none"> • Place the instruments in the cleaning bath for the specified reaction time so that the instruments are adequately covered (if necessary, ultrasonic assistance or careful brushing with a soft brush). When doing so, make sure that the instruments do not touch one another. • Then remove the instruments from the cleaning bath and rinse them thoroughly for at least 1 minute using running water if possible. • Check the instruments (see section "Checking" and "Maintenance"). <p>Disinfection</p> <ul style="list-style-type: none"> • Place the cleaned and checked instruments in the disinfection bath for the specified reaction time so that the instruments are adequately covered. When doing so, make sure that the instruments do not touch one another. • Then remove the instruments from the disinfection bath and rinse them thoroughly for at least 1 minute using running water if possible. • Dry and pack the instruments immediately after removal if possible (see section "Packing", if necessary after further drying in a clean place). <p>Evidence of the basic suitability of the instruments for effective manual cleaning and disinfection has been provided by an independent accredited test laboratory using the cleaning agent Bodedex forte and the disinfectant Korsorex plus (Bode Chemie GmbH & Co., Hamburg). This was carried out using the procedure described above.</p>
Automated cleaning and disinfection	<p>When selecting the cleaning/disinfection unit (CDU, disinfector), it must be ensured</p> <ul style="list-style-type: none"> • that the CDU basically has proven effectiveness (e.g. DGHM or FDA approval or CE marking according to DIN EN ISO 15883), • that, if possible, a proven programme for thermal disinfection is used (A_0 value > 3000 or - in the case of older units - at least 10 min at 93 °C) (in the case of chemical disinfection there is a risk of disinfectant residue on the instruments), • that the programme used is suitable for the instruments and has sufficient rinsing cycles,

11/11/11

11/11/11

Instructions for use

Reusable surgical instruments

	<p>IQ/OQ (commissioning) and product-specific performance qualification (PQ))</p> <ul style="list-style-type: none"> Maximum sterilisation temperature 138 °C (280 °F; plus tolerance according to DIN EN ISO 17665 (formerly: DIN EN 554/ANSI AAMI ISO 11134) Sterilisation time (exposure time at sterilisation temperature) at least 20 min (at 121 °C (250 °F)) or 5 min at 132 °C (270 °F)/134 °C <ul style="list-style-type: none"> ¹ The use of the less effective gravity displacement method is only permitted if the vacuum method is not available. <p>Evidence of the basic suitability of the instruments for effective steam sterilisation has been provided by an independent accredited test laboratory using the EuroSelectomat steam steriliser (MMM Münchener Medizin Mechanik GmbH, Planegg) in conjunction with the fractionated vacuum method, and also using a steriliser from H+P Labortechnik GmbH, Oberschleißheim in conjunction with the gravity displacement method; in addition, the Systec V-150 steam steriliser (Systec GmbH Labor-Systemtechnik, Weltenberg) was also used for both methods. In doing so, typical conditions in clinics and doctors' surgeries as well as the procedure described above were taken into account.</p> <p>The flash sterilisation method is basically not allowed.</p> <p>Also, hot-air sterilisation, radiation sterilisation, formaldehyde or ethylene oxide sterilisation and plasma sterilisation must not be used.</p>
Storage	Sterilised instruments must be stored in the sterilisation pack, dry and free from dust, in a closed cupboard.
Material compatibility	<p>When selecting cleaning agents and disinfectants, please ensure that they do not contain the following constituents:</p> <ul style="list-style-type: none"> Acids (< pH 5)/oxidising acids Alkalis (> pH 10) Organic solvents Benzine, phenol or ammonia Halogens, halogenated hydrocarbons, sodium chloride (in higher concentration), Oxidants / peroxides / hypochlorite <p>No instruments must be subjected to temperatures higher than 141 °C (286 °F)!</p> <p>All liability is excluded if this is disregarded!</p>
Special information	<p>The aspects listed below must also be taken into account with the following products/product groups:</p> <ul style="list-style-type: none"> Instruments with joints (e.g. forceps, clamps, scissors, needle holders, abdominal retractors etc.) <ul style="list-style-type: none"> Open and close the instrument several times during pre-treatment as well as during cleaning and disinfection Open and close the instrument several times during pre-treatment and place in the disinfectant in the half-open position when carrying out automatic cleaning and disinfection If required, oil the joint (but not other surfaces of the product) using as little oil as possible Instruments with channels/cavities (e.g. suction tubes, suction curettes etc.) <ul style="list-style-type: none"> Rinse through at least 5 times during manual cleaning and disinfection and also at every rinsing step using a single-use syringe (minimum volume 5-50 ml depending on the size of the lumen) Connect directly to the rinsing connector of the CDU using a suitable adapter if necessary during automatic cleaning and disinfection Do not continue to use the instrument if the patency is impaired Instruments with teeth (e.g. saws, rasps, ...) <ul style="list-style-type: none"> Ensure that residues are completely removed from the teeth
Additional information	<p>Further information on the preparation of medical products:</p> <ul style="list-style-type: none"> Internet: http://www.rki.de Internet: http://www.a-k-i.org Hygiene requirements in the preparation of medical products. Recommendation by the Commission for Hospital Hygiene and the Prevention of Infection at the Robert Koch Institute (RKI) and the Federal Institute for Drugs and Medical Devices (BfArM) with regard to the "Hygiene requirements in the preparation of medical products"

Instructions for use

Reusable surgical instruments

	<p>DANGER! Danger to life!</p> <p>If this product is used on a patient with Creutzfeldt-Jakob Disease (CJD) or even if vCJD is only suspected, the product must not be re-used and must be destroyed.</p>
---	--

Reusability

Provided that appropriate care is taken and as long as the products are undamaged and not contaminated, the preparation process has little effect on the instruments. The end of the product's life is normally determined by the wear and usage of the instrument.

Warranty

As the manufacturer and distributor of the products, GEOMED Medizin-Technik GmbH & Co. KG accepts **no** liability for immediate damage or consequential damage caused by improper use, handling or by improper preparation, sterilisation and maintenance.

	<p>HERSTELLER und SERVICE Adresse:</p> <p>GEOMED Medizin-Technik GmbH & Co. KG Ludwigstalerstrasse 27 D-78532 Tuttlingen / Germany</p> <p>Tel: #49 (0)7461 – 9355-0 Fax: #49 (0)7461 – 9355-93 Mail: info@geomed.de</p>
	<p>CE-Kennzeichnung gemäss Richtlinie 93/42/EWG CE marking according to directive 93/42/EEC Marquage CE conforme à la directive 93/42/CEE Identificación CE en conformidad con la directriz 93/42/CEE</p>
	<p>CE-Kennzeichnung gemäss Richtlinie 93/42/EWG CE marking according to directive 93/42/EEC Marquage CE conforme à la directive 93/42/CEE Identificación CE en conformidad con la directriz 93/42/CEE</p>