Nội dung

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
1 Product Liability and Warranty

1.1 General Information

We thank you for having decided to buy a KLS Martin product. This product carries the CE mark, which means that it satisfies the essential requirements laid down in the EC Directive concerning medical devices.

We are the manufacturer of this product:

Gebrüder Martin GmbH & Co. KG

A company of the KLS Martin Group
KLS Martin Platz 1 · D-78532 Tuttlingen · Germany
Postfach 60 · D-78501 Tuttlingen · Germany
Tel. +49 7461 706-0 · Fax +49 7461 706-193
info@klsmartin.com · www.klsmartin.com

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1.2 Intended Use

The instruments may only be used for their intended purpose in the specified medical fields, with use being restricted to adequately trained and qualified personnel. The treating physician or user is responsible for selecting the right instrument(s) for the surgical task / application at hand as well as for their safe handling. This includes ensuring an adequate level of training, knowledge and experience.

1.3 Warranty

NOTICE

Risk of damage to the instruments due to improper use!

The responsibility for proper instrument cleaning, disinfection und sterilization rests with the operator / product user. Be sure to observe your national / local regulations, including potential restrictions.

Gebrüder Martin, as the manufacturer of the products, accepts no liability for direct or consequential damage caused by improper use, handling, processing, sterilization or maintenance. Unauthorized instrument repair (by firms or persons not specifically authorized by Gebrüder Martin to perform such work) shall void the warranty given!

Non-observance of these notices, as well as improper handling or use of products supplied by us, will void your rights under the warranty. Consequently, Gebrüder Martin shall not be liable for any resulting damage in such cases.

1.4 User's Inspection

Immediately upon receipt, the goods must be checked for completeness and potential damage in transit. Notice of any such damage must be given immediately.

Before using the instruments for the first time, as well as before each subsequent use and before returning them for repair or servicing, they must be cleaned, disinfected and sterilized in accordance with our instrument processing instructions, which are available on the Internet at www.klsmartin.com/phnet/ but can also be obtained from KLS Martin on request.

1.5 Hotline

- Should you have any questions on how to handle the unit / product or question on its clinical application, please do not hesitate to contact the Product Management:
Tel: +49 7461 706-234
Fax: +49 7461 706-190
- Should you have any technical questions, please do not hesitate to turn to our Martin Service Center:
Tel: +49 7461 706-343
Fax: +49 7461 706-408
E-Mail: service@klsmartin.com
- Should you have any questions concerning maintenance contracts or training courses, please contact our Technical Service Manager:
Tel: +49 7461 706-332
E-Mail: service@klsmartin.com

NOTICE

Each product and its packaging are marked with a specific lot number (LOT). Please always indicate this number together with the item number in case of complaint.

2 Notices concerning this Document

Improper handling and care as well as non-intended use can lead to premature wear and / or pose a risk to patients and users!

Be sure to read, understand and follow the instructions given below!

- Every user is required to read this document completely and follow them carefully.
- In particular, be sure to heed all cautions, warnings and danger notices.
- Keep this document accessible to users at all times.

2.1 Symbols Used in this Document

Throughout this document, important information (such as general or safety-related notices) is marked with the following symbols and signal words:

▲ WARNING

Danger of death or serious injury!

Indicates a situation which, if not avoided, **could** result in death or serious injury!

▲ CAUTION

Danger of minor injury!

Indicates a situation which, if not avoided, could result in minor or moderate injury!

NOTICE

Risk of material damage!

Indicates a situation which, if not avoided, could lead to material damage (loss of time, data loss, device / machine failure, etc.)!

2.2 Safety Checks

If the device is not functionally reliable and / or safe to operate, it must be marked as such and withdrawn from service. A technical check is mandatory in any such case.

3 Use

3.1 Prior to First Use

⚠ WARNING

Danger of infection due to non-sterile handling!

Improper sterilization and non-sterile handling of the instruments can lead to serious health hazards for patients.

All instruments must be cleaned and sterilized before using them or the first time, as well as prior to each subsequent use. See section 4 "Processing, Cleaning, Care, Disinfection and Sterilization", page 25.

In addition, we recommend sending the products through the cleaning and disinfecting cycle three times prior to first use. Current studies have shown that such pre-treatment supports the natural formation of a passivation layer right from the start.

4 Processing, Cleaning, Care, Disinfection and Sterilization

4.1 Inspection & Testing Prior to Reuse

Before each use, the instruments must be thoroughly inspected for damage such as fractures, cracks or deformation, as well as for functional reliability. Special attention must be paid to cutting edges, tips, joints, box locks, ratchets and all movable parts. If wear, corrosion, deformation, porosity or other damage is detected, the instrument must be immediately withdrawn from service.

Due to their alloy, stainless steel instruments typically develop a passive film in the form of a protective layer. However, this film does not protect them well against chemical attack by chloride ions and aggressive media and liquids!

Therefore, in addition to the instrument manufacturer's endeavors to select the right materials and process them carefully, the user must make an important contribution by ensuring proper instrument processing along with adequate and regular care.

NOTICE

Risk of damage to the instruments due to improper use!

The responsibility for proper instrument cleaning, disinfection and sterilization rests with the operator / product user. Be sure to observe your national / local regulations, including potential restrictions.

We recommend using the procedures specified in sections 4.2 to 4.8 for reprocessing our reusable surgical instruments.

4.2 Machine Cleaning

4.2.1 General Notices

It is always preferable to use machines (washer-disinfectors) for instrument cleaning and disinfection because, unlike manual procedures, machine processes can be easily standardized.

Be sure to observe and follow the operating and loading instructions provided by the machine manufacturer. In addition, only the cleaning agents recommended by the manufacturer should be used for the application at hand.

- Jointed instruments should always be processed in open condition. Be sure to arrange the items so that the water can easily flow out of cannulations, blind (non-through) holes and hollow bodies.
- Complex instruments must be totally taken apart before cleaning.
- For instruments with long or narrow lumens, standard procedures should be used only if the hot disinfectant can easily flow through the lumens and safe rinsing is guaranteed.
- The instrument trays used for cleaning must always be loaded correctly to ensure proper cleaning.
- Store the instruments correctly in the tray. Be sure to prevent "rinsing shadows"!
- For ergonomic reasons, a total tray weight of 8 – 10kg should not be exceeded.
- When removing the instruments from the machine after cleaning, be sure to check them for cleanness (visible dirt). This especially applies to cannulated instruments or those with blind holes. If necessary, repeat the cleaning cycle or pre-clean manually as specified in section 4.3 "Manual Cleaning", page 28.
- The final rinse must guarantee that any residues left after the cleaning stage are completely removed. Be sure to use only fully demineralized water for the final rinse!

4.2.2 Wet and Dry Disposal

“Wet disposal” and “dry disposal” are basically different cleaning procedures.

- **Dry disposal:**
Following use in the OR, the instruments are deposited for return to the CSSD in dry condition, without applying a disinfectant or immersing them in a disinfecting solution. So for dry disposal, the instruments are not to be cleaned with or placed in a physiological saline solution.
If pre-cleaning is necessary after use, the instruments must be rinsed or wiped clean immediately.
We recommend using the “dry disposal” method for machine cleaning.
- **Wet disposal:**
Immediately after use, the instruments are to be placed into a non-fixating cleaning and disinfecting solution.
Please observe the instructions provided by the manufacturer of your disinfectant.

4.2.3 Detergents and Disinfectants & Temperature Ranges

- Be sure to use commercially available detergents and disinfectants that have been approved for use in each case.
- Always observe the recommended concentrations, exposure times and temperatures.
- Make sure that no residues are left on the instruments after cleaning.
- Use demineralized water for the final rinse.

	Stainless Steel Instruments	Aluminum Instruments
Cleaning/disinfection	acid / neutral / alkaline with / without the addition of tensides	neutral / mildly alkaline with / without the addition of tensides
Allowable temperature for chemical disinfection	max. 60 °C / 140 °F	max. 60 °C / 140 °F
Allowable temperature for thermal methods	max. 93 °C / 199 °F	max. 93 °C / 199 °F
Drying	at max. 125 °C / 257 °F	

NOTICE

Use of neutral-pH treating agents recommended!

Alkaline detergents and/or acid neutralizers are capable of causing discoloration – usually grayish brown – on metal-coated (e. g. TiNi) surfaces even after a short time. According to the current state of research, this does **not** compromise the proper functioning of the instruments. Nonetheless, using neutral-pH agents is the best option to ensure gentle treatment of your instrument surfaces.

NOTICE**Media containing chlorides attack surfaces!**

Water with a chloride content exceeding 120 mg/l can attack the surfaces of your instruments during cleaning. If you use water with a chloride content of less than 120 mg/l, you still need to take the concentration effect during drying into account as well.

NOTICE**Ensure reliable drying!**

Reliable drying is an essential factor for successful sterilization!

4.3 Manual Cleaning

Instruments should be disinfected and cleaned immediately after use, following the instructions given under "Cleaning and Sterilization of Surgical Instruments", see section 1.4, "User's Inspection", page 23. Be sure to prevent residues from drying on the instruments as this would make proper cleaning and disinfection more difficult than necessary.

The following should be observed in manual cleaning processes:

- The solutions used for manual cleaning must always be prepared in accordance with the manufacturer's instructions.
- Use a suitable brush for cleaning lumens, cannulations, blind holes and cavities, making sure that every part of the inner surface can be properly accessed.
- Use a soft brush and a neutral or mildly alkaline detergent for removing blood and other residues.
- Never use metal brushes or metal sponges for manual cleaning.
- To ensure proper instrument functioning, verify that all movable parts have been thoroughly cleaned.
- Clean jointed instruments in closed as well as open condition.
- Take instruments fully apart where applicable.
- Pay special attention to slots, ratchets, joints and box locks, narrow lumens, blind holes and other areas that are hard to access.
- Suitable trays or baskets (e. g. sterilization trays or wire baskets) must be used for storing surgical instruments properly during the cleaning process.

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4.4 Ultrasonic Treatment

Effective ultrasonic cleaning requires placing the surgical instruments into suitable sterilization trays or wire baskets in an open condition. As the use of plain warm water alone cannot produce satisfactory results, it is necessary to add a suitable cleaning agent. Be sure to observe and follow the manufacturer's instructions with regard to concentration as well as temperature. Since an excessive dirt content of the cleaning solution has an adverse effect on the cleaning results, the solution must be replaced at regular intervals in accordance with the manufacturer's instructions. In the same manner, prescribed immersion or ultrasonic treatment times must be strictly observed. As a rule, ultrasonic cleaning must always be followed by a rinsing cycle. Be sure to check the instruments (where applicable) for loosened components after the ultrasonic bath. To prevent water spots ("spotting"), fully demineralized or distilled water should be used for the final rinse.

4.5 Chemical Disinfection

- The solutions employed for chemical disinfection must always be used in accordance with the manufacturer's instructions.
Never use hypochlorite-based disinfectants containing active chlorine (e. g. chlorine bleaching agent (sodium hypochlorite) or Javel water) or disinfectants containing iodine.
- Pure water must be used for preparing the dilutions specified for the chemical disinfectants. The addition of detergents is not permitted. The manufacturer's instructions regarding exposure times and concentration must be duly observed in each case.
- Disinfecting solutions must be prepared afresh every day. Extended / multiple use can easily lead to the following problems:
 - Increased concentration due to evaporation (corrosion risk)
 - Excessive dirt load (corrosion risk plus lower effectiveness)
- Following disinfection, it is important to rinse all items sufficiently under clear running water. To prevent the formation of water spots, the use of fully demineralized water is highly recommended.
- Immediately after completion of the cleaning and rinsing cycles, surgical instruments must be sufficiently dried.

4.6 Inspection after Cleaning

- Following cleaning, the instruments must be macroscopically clean, i.e. free from visible dirt or deposits.
- Instruments with stains or spots must be withdrawn from service at once and given special treatment.
- All movable parts, working tips and (scissor) blades should be inspected with particular care.
- If damage or malfunction is detected, the instrument must also be withdrawn from service immediately.

4.7 Instrument Care

"Care" means treating the instruments with instrument oil or milk (white oil-in-water emulsion). Instruments with joints or box locks (scissors, forceps, clamps, etc.) or with metal sliding surfaces (rib shears, punches, etc.) must be treated with steam-sterilizable care agents based on paraffin oil. The paraffin oil must comply with the pharmacopoeia in force at the time and must be physiologically safe as specified in the German Pharmacopoeia ("Deutsches Arzneibuch", 10th edition (DAB 10)", "European Pharmacopoeia (Ph. Eur.)" or "United States Pharmacopoeia (USP)".

Care agents prevent metal-on-metal friction, thus ensuring the easy movement of your instruments. Laser-marked products can be adversely affected when using basic cleaning agents containing phosphoric or hydrofluoric acid because the marking may fade away and the coding function get impaired or lost as a result.

As a rule, surgical instruments must be subjected to regular care, which means each time before a functional test is carried out. At the same time, it is important to prevent "gumming" of the joints due to an accumulative effect, especially in instruments that are continuously in use.

4.8 Sterilization

4.8.1 Prior to Sterilization

- Prior to sterilization, the instruments must be adequately packaged, e. g. using containers that meet EN 868-8 requirements.
- The packaging method used must comply with the relevant standards.
- Check the instruments for cleanness and integrity.
- Clean and disinfect the instruments and rinse them with distilled water, then dry them carefully, see sections 4.2 to 4.5.

4.8.2 Steam Sterilization

Danger of infection due to non-sterile handling!

Improper sterilization and non-sterile handling of the instruments can lead to serious health hazards for patients.

Sterilization must be carried out according to a validated steam sterilization process, for example in a sterilizer satisfying EN 285:2009, ANSI / AAMI ST 79 requirements and validated in accordance with ISO 17665-1:2006 requirements.

It is essential to keep the steam free from any foreign bodies such as rust particles and other impurities. This helps to prevent instrument corrosion or surface damage caused by deposits. The steam used for sterilization must comply with EN 285:2009. The user instructions provided by the steam sterilizer manufacturer must be duly observed. Instruments incorporating locks or ratchets must be sterilized in an open condition or with the ratchet set to the first notch.

The following program variants can be used:

- Validated steam sterilization process, "134 °C (273 °F) / 2 bar" program
- Validated steam sterilization process, "121 °C (250 °F) / 1 bar" program

The sterilization and holding times are subject to national provisions and regulations and therefore cannot be defined in general. It is the operator's responsibility to ensure that the desired results will be achieved with the implemented processing and sterilization procedure, including the equipment and materials used and the staff employed in the processing and sterilization department (CSSD). This requires validation and routine monitoring of the process used.

In specific, the following sources should be observed for cleaning, instrument care, disinfection and sterilization:

- ISO 17664: "Sterilization of medical devices/Information to be provided by the manufacturer for the processing of resterilizable medical devices"
- EN 285: "Sterilization - Steam sterilizers - Large sterilizers"
- ISO 17665-1: Sterilization of health care products - Moist heat - Requirements for the development, validation and routine control of a sterilization process for medical devices
- EN 556-1: Sterilization of medical devices - Requirements for medical devices to be designated "STERILE" - Part 1: Requirements for terminally sterilized medical devices
- Recommendation of the Instrument Preparation Working Group (Arbeitskreis Instrumenten-Aufbereitung / AKI), <http://www.a-k-i.org>, "Proper Maintenance of Instruments"
- Recommendation of the German Society for Hospital Hygiene (Deutsche Gesellschaft für Krankenhaushygiene / DGKH), <http://www.dgkh.de>, "Hygiene Requirements for Reprocessing Medical Devices"
- Recommendation of the Commission for Hospital Hygiene and Infection Prevention at the Robert Koch Institute (RKI) and the Federal Institute for Drugs and Medical Devices (BfArM), <http://www.rki.de>, "Hygiene Requirements for Reprocessing Medical Devices"

4.8.3 Information Concerning Alternative Sterilization Methods

Steam sterilization has established itself as a very safe and reliable sterilization method all over the world and therefore is considered the method of choice for goods to be sterilized that are insensitive to temperature and humidity. To this end, reference is usually made to steam sterilization using a validated steam sterilization process (see also EN 554). Consequently, there is no necessity to use other sterilization methods – e. g. low-temperature plasma (LTP) sterilization or formaldehyde or ethylene oxide – for sterilizing steam-sterilizable medical devices. Since plasma sterilization is the subject of controversial discussion among experts with regard to its effectiveness in cavities and lumina, Gebrüder Martin does not perform validation of sterilization processes using gas/plasma sterilizers for steam-sterilizable medical devices. However, sterilization plant operators are free to validate their own alternative sterilization procedures for the medical devices to be sterilized.

5 Storage and Transportation

- Store the instruments in a clean, cool and dry place.
- Protect them against mechanical damage.
- Use adequate containers / packaging for safe storage and transportation.
- Handle with utmost care; never throw these products or allow them to fall down.
- Use approved sterilization packaging (complying with EN 868 / ISO 11607 requirements, for example) for sterilization and subsequent transportation and storage.
- When returning products, be sure to clean and disinfect all items and use sterile packaging.

6 Processing Restrictions, Disposal

Frequent reprocessing has only a minor influence on the service life of surgical instruments, which is mainly a function of use-related wear and tear or damage. Please dispose of your instruments in accordance with relevant local regulations, or have them properly recycled, once they have reached the end of their life cycle.

The national regulations for waste disposal must be complied with!

7 Supplementary Instructions on Processing marCore[®] MIS Instruments

These additional instructions are intended only for marCore[®] MIS instruments with item numbers ranging from 24-599-00-04 to 24-599-99-04.

The instructions provided in the foregoing sections of this document remain fully valid for marCore[®] MIS instruments as well. They are merely supplemented as follows.

7.1 Product Overview


Fig. 7-1: Product overview


Fig. 7-2: Instrument – example

- | | |
|---------------|--------------------------------|
| 1 Working end | 4 Nut |
| 2 Shank | 5 Cleaning connector Luer-Lock |
| 3 Handle | |

The handle (3) and the nut (4) are made of titanium and have been color-coded. Gentle cleaning of these parts requires the use of neutral-pH or mildly alkaline detergents. The shank (2) and the working ends (1) of these instruments are stainless steel with a protective metal coating that increases the service life of the instruments a great deal as long as they are used as intended. The cleaning connector (5) of the instruments is also stainless steel.


Fig. 7-3: Color coding

7.2 Supplements to the Previous Sections

Supplement to Section 3 Use

The working ends of the instruments including the joint should be cleaned intraoperatively after each instrument exchange between the operating personnel and the surgeon using a physiological saline solution (Ringer's solution).

Following use, the instruments should be placed back, in open condition, into the storage tray intended for this purpose, and returned to the CSSD for reprocessing.


Fig. 7-4: Instruments in the mesh tray

4 Joint at the working end

Supplement to Section 4 Processing, Cleaning, Care, Disinfection and Sterilization

The instruments can be machine-processed in the same way as other instruments provided they are fixed in place securely, e. g. in special mesh trays or MIS mobile units, and the washing/flushing system used is adequate for processing MIS instruments.

Supplement to Section 4.2 Machine Cleaning, Section 4.2.1 General Notices

- The instrument trays used for cleaning must be loaded so that the desired cleaning result will not be compromised (see e.g. Fig. 7-1).
- The needle holders must be set to neutral position so that the ratchet is open.
- The cleaning and disinfecting agents used must be compatible with the instrument materials and provide the required hygienic and microbiological effectiveness. A thorough rinse must ensure that no cleaning or care agent residues are left on the instrument surfaces.
- Special mesh trays or suitable fixation means must be used for the instruments. If necessary, use a mobile unit (slide-in load carrier) providing special flushing technology with Luer-lock cleaning connectors (siehe Fig. 7-5) to ensure reliable cleaning and internal flushing of the instruments.


Fig. 7-5: Cleaning connector

5 Luer-lock cleaning connector for connection to the flushing system of a MIS mobile unit

Instruments made of different materials

Metal-coated (e. g. TiNi) surfaces and all titanium parts (see Fig. 7-6) can be adversely affected even after a short time when using alkaline cleaning agents and / or acid neutralizers. This usually takes the form of a slight, grayish or brownish discoloration of the coated surfaces or fading of the color-coded handles and nuts. Judging from the current state of research, this does not in any way compromise the proper functioning of the instruments. Nonetheless, using neutral-pH agents is the best option to ensure gentle treatment of your instrument surfaces.

Only neutral or mildly alkaline cleaning agents suitable for the titanium handles are permitted for use.

Titanium handles

Shank with titanium nitride coating

Fig. 7-6: Titanium and titanium nitride coating

Supplement to Section 4.3 Manual Cleaning

In this context, please note the following:

- All instruments must be stored on instrument trays allowing adequate cleaning, e. g. sterilization trays / mesh trays.
- All instruments must be flushed with a cleaning gun via the dedicated (Luer-lock, see Fig. 7-5) cleaning connector and must subsequently be dried with medical compressed air.
- The joints at the distal end (working end) must also be flushed several times with a cleaning gun via the dedicated (Luer-lock, see Fig. 7-5) cleaning connector. They too must then be dried with medical compressed air.

Supplement to Section 4.4 Ultrasonic Treatment

To avoid damage, the instruments must be fixed in place using special holding devices (see Fig. 7-4).

Supplement to Section 4.6 Inspection After Cleaning

To prevent transport damage, put the instruments back into the dedicated mesh trays after inspection and secure them using suitable fixation devices.

Supplement to Section 4.7 Instrument Care

Following cleaning, all instruments must be oiled in open condition, both at the joints of the working ends (4), and at the handle joints (6), applying oil from the front and from the side at the points shown below. For care agent suitability, see section 4.7 "Instrument Care", page 30.


Fig. 7-7: Instruments in the mesh tray

4 Joint at the working end 6 Handle joint