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Legend relating to illustrations on fold-out page

No. Description

- 1 Container cover
- 2 Paper filter
- 3 Filter holder
- 4 Button
- 5 Cover seal
- 6 Container box
- 7 Handle plate
- 8 Carrying handles
- 9 Coding label
- 10 Seal slot
- 11 Color field
- 12 Slot section for log label
- 13 Locking plate

1 Product Liability and Warranty

1.1 General Information

We thank you for having decided to buy a Gebrüder Martin product. This product carries the CE mark, which means that it satisfies the essential requirements laid down in the EC Directive concerning medical devices.

We are the manufacturer of this product:
 Gebrüder Martin GmbH & Co. KG
 A company of the KLS Martin Group
 KLS Martin Platz 1 • D-78532 Tuttlingen • Germany
 Postfach 60 • D-78501 Tuttlingen • Germany
 Tel. +49 7461 706-0 • Fax +49 7461 706-193
 info@klsmartin.com • www.klsmartin.com

1.2 Warranty

Our Standard Terms and Conditions of Sale effective at the time shall apply. Agreements diverging from these Standard Terms and Conditions do not restrict the legal rights of the buyer. Any warranty exceeding the above provisions shall require a contractual form and shall exclude component-related vandalism and consumables. Improper interventions or alterations performed by third parties during the period of limitation shall void any and all warranty claims. Unauthorized actions performed on the product shall invalidate any liability claims against Gebrüder Martin.

2 Product Overview

The product is illustrated on the fold-out page (cover page) of these Instructions for Use. For a description of the various features, see opposite page.

3 Intended Use

The marSafe® Sterilization Container is to be used for storing, sterilizing and transporting sterile supplies and placing them ready for use under sterile conditions. All models and their components fully meet the requirements of the relevant international and national standards, specifically ISO 11607-1:2009 and ISO 11607-2:2006, EN 868-8:2009 and DIN 58953-9:2010.

3.1 Sterilization

The containers are suitable for steam sterilization. It is essential, however, to use a validated steam sterilization process (e.g. a sterilizer complying with EN 285:2009, validated in accordance with ISO 17665-1:2006 requirements).

- Sterilization containers with perforated cover and closed container box are validated for steam sterilization in a sterilizer in accordance with EN 285 by **fractionated vacuum processes**.
- Sterilization containers with perforated cover and perforated container box are also suitable for steam sterilization in a sterilizer in accordance with EN 285 by **gravitation processes**.

4 Safe Usage

Before use, check device for functionality and proper condition. Ensure that only qualified personnel with the required training or experience use the sterile containers.

- Keep Instructions for Use accessible to personnel.
- Observe general guidelines concerning contaminated and sterilized goods as well as goods to be sterilized.

5 Structure and Function of the Container

5.1 Paper filter

The microbial barrier of the marSafe® sterilization container consists of a paper filter (2) and the corresponding filter holder (3). The paper filter is intended for single use only. Remove the filter holder by pressing the button (4). Insert the paper filter. Press the filter holder in until you hear it lock in place.

5.2 Cover Seal

An all-around silicone seal (5) is provided on the inside of the container cover to ensure the microbially tight fit of the container cover on the container box (6). The seal is easily accessible for cleaning.

5.3 Container/Cover Lock

Closing the container

Place the container cover (1) onto the container box (6). Lift up the locking plate (13) and place it over the edge of the cover. Finally, push the locking plate of the toggle lever lock against the container box. Place the container box on a solid surface to facilitate the closing process.

Opening the container

The container is opened by tilting the locking plate (13) of the toggle lever lock upwards. Then lift up the locking plate and fold it outward. Finally, lift off the container cover (1). Seals should be removed before opening the container cover.

5.4 Protective Cover

The optional protective cover (55-290-99-04) protects the microbe retention system from potential external influences. The protective cover is fixed using two bolted retainers. The protective cover can also be sterilized and does not have to be removed.

5.5 Seal Slot

The seal slot (10) of the handle plate (7) that can be fitted with two seal types protrudes from the center of the locking plate (13) when the container is closed. The seals (item no. 55-996-84-04 and 55-996-85-04) are inserted using the lateral eyelets and then closed. The seals are inserted into the slot from above and lock automatically.

5.6 Carrying Handles

The front sides of the container each feature a carrying handle (8). The carrying capacity of the carrying handles complies with the requirements set out in EN 868-8:2009, Appendix C.

5.7 Coding Labels and Color Fields

A color field (11) is integrated into the locking plate (13) on both front sides of the container; on the left it serves for attaching a coding label (9). Laser-lettered coding labels with customized information concerning tray contents and destination are available from Gebrüder Martin upon request. The coding label is inserted on the left behind the color field of the locking plate. Color fields are available in 6 different colors to allocate labels to specific medical fields. The color fields can easily be inserted and removed by the user.

NOTICE

A color field must be inserted to ensure the correct sealing of the containers.

5.8 Log Cards or Labels

Both front sides of the container feature a slot section (12) with a retaining spring in the color field (11) for accommodating and fixing universal log labels (e. g. item no. 55-291-11-04, 55-291-12-04, 55-291-13-04). A defective retaining spring can be replaced by the user.

6 Initial Commissioning

The container must be thoroughly cleaned before loading and the first sterilization cycle (see section 11 "Disinfecting, Cleaning and Care", page 8).

To assemble the container, follow the instructions in section 5 "Structure and Function of the Container", page 7.

Then carry out a functional check in accordance with section 7.1 "Functional Check Prior to Loading the Container", page 8.

For the validation of the packaging process, the international DGSV (German Association for Sterile Services) guideline "Validation of Packaging Processes according to DIN EN ISO 11607-2" can be applied.

7 Loading the Container

When loading the container, a clearance of about 1 cm must be maintained below the filter holder (3) to ensure unhindered flow of the steam into the container box.

7.1 Functional Check Prior to Loading the Container

Prior to using the container, all container components must be visually inspected for potential defects. Any damaged components found must be repaired or replaced by original parts. The cover seal must be free from visible defects, or it must be replaced. Only cleaned container assemblies may be used.

7.2 Loading the Container with Linen or other Textiles

When loading the container with textiles, be sure to pile the folded items (max. 6 kg) upright into the container. Do not overfill! It should always be easily possible to slide the stretched hand into the pile. In accordance with the German standard DIN 58953-9:2010, we also recommend wrapping the items to be sterilized with a suitable, adequately sized cloth inside the container. This ensures easy handling of the sterile supplies and contamination-free provision for use later. Use of a gravitation process when loading the container with textiles must be approved by the manufacturer.

7.3 Loading the Container with Instruments

We recommend using trays for sterilizing instruments inside the container in accordance with DIN 58952-2:2012 and DIN 58952-3:2012, fixing the instruments in place in the tray with storage aids such as silicone strips.

According to DIN EN 868-8:2009 and DIN 58953-9:2010, as well as for ergonomic reasons, the following maximum loads are recommended:

Size [cm ²]	Weight [kg]
60 x 30	10
47 x 30	7
30 x 30	5

Due to their special design, container boxes allow the tray baskets to be taken out of the box under sterile conditions after the container has been opened in the operating room, even if no additional internal packaging is used.

8 Loading and Unloading the Sterilizer

The loading instructions provided by the manufacturer of the sterilizer must always be duly observed. Heavy containers should always be placed at the bottom of the sterilizing chamber. *marSafe*[®] containers can be stacked both during sterilization and storage and are also stackable with the MicroStop[®] system. Always use the carrying handles and never the container cover for lifting or carrying the container when loading or unloading it. The heaviest containers should always be placed at the bottom of the sterilizing chamber.

NOTICE

The container must not be additionally packed or covered in the sterilizer. The perforation fields in the container box or container cover must not be covered. Do not place film packaging or silicone mats directly on the container.

9 Storing Sterile-Supply Containers

The containers must be stored so that they are protected against dust, moisture and damage. The storage life depends on local storage conditions and cannot be specified universally. Therefore, the Hygiene Board must determine the acceptable storage period. The responsibility for storage conditions and periods always lies with the hospital management or the medical director. Sterile containers must be kept separate from any goods that represent a potential danger. Loss of the quality of seal on the sterile packaging is generally event-related and not time-related.

10 Placing the Sterile Supplies Ready for Use

Prior to opening the containers, verify that the seals are correctly in place and fully intact. Check in particular whether the indicators show the correct color change. Before opening the container, destroy and remove the seals according to the instructions. The seals 55-996-84-04 and 55-996-85-04 can easily be removed using the predetermined breaking point. The seals 55-996-82-04 and 55-996-83-04 must be bent at the tab and turned clockwise and then be removed in downward direction.

11 Disinfecting, Cleaning and Care

11.1 General Information

The parameters for cleaning and disinfecting the container are to be determined by the operator / user, e. g. as part of the hygiene plan. The treatment agents must be selected to be compatible with the container component materials and also effective hygienically and microbiologically. A thorough final rinse, preferably with demineralized water, must ensure that no cleaning and care agent residues are left on any of the components of the container.

NOTICE

Observe all national legal requirements, national and international standards and guidelines and the in-house hygiene regulations for reprocessing.

- Observe all applicable national regulations in regard to reprocessing of the products in case of patients suffering from Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of the disease.
- Note that successful reprocessing of this medical product can be guaranteed only after prior validation of the reprocessing process. The responsibility for this lies with the operator/reprocessor.
- Due to process tolerance, the data provided by the manufacturer serve only as reference values for the evaluation of reprocessing processes carried out by the operator/reprocessor.

- The container cover (1) and container box (6) must be separated for cleaning.

The following sections provide information on manual and mechanical cleaning and disinfecting options, focusing on material compatibility and value conservation.

11.2 Materials Used

The container box (6), the handle plate (7), and the locking plate (13) are made of aluminum and feature a protective anodized coating that makes the surface more abrasion-resistant, thus enhancing its resistance to wear. However, while anodized surfaces are resistant to pH-neutral or alkaline cleaners up to a pH of 10.5, they are susceptible to destruction by chemical attack, particularly when alkaline agents with a pH of more than 10.5 and acidic neutralizers are used in machine-based processing. Therefore, the container box may be treated only with pH-neutral or alkaline agents up to a pH of 10.5 or with cleaning agents explicitly approved for cleaning anodized aluminum surfaces. Be sure to observe the manufacturers' instructions relating to both your automatic cleaning device and the cleaning agents used.

The color field (11) is made of a special high-performance plastic material suitable for sterilization. Both parts can be treated with any program suitable for cleaning aluminum products. Thermal disinfection up to 95°C (203°F) is possible as well, posing no problem.

NOTICE

Color fields must never be treated with ACETONE/PETROLEUM ETHER! In addition, be sure to use cleaning agents only in concentrations that do not exceed the maximum specified by the manufacturer (this applies to both manual and machine cleaning processes). Note that ACETONE/PETROLEUM ETHER as well as excessive cleaner concentrations cause stress cracking, thus destroying the material. For the same reason, no disinfectants containing glucoptamine should be used for disinfecting the *marSafe*[®] container cover because such agents would corrode the silicone gasket, thus reducing its life or even destroying it completely.

11.3 Manual Cleaning and Disinfection

All accessible components can be cleaned with neutral, surfactant-based cleaning agents (mild dishwashing liquids; neutral cleaning agents). For disinfection, aldehyde-free surface disinfectants may be used. The final rinsing may only be carried out using demineralized water. Do not use metal brushes or scouring agents for cleaning.

After this procedure, all parts should be checked for visible dirt. If there is any residual dirt, repeat the process.

Recommendations for manual cleaning and disinfection:

Phase	Step	Temperature	Water quality	Chemicals		Note
1	Cleaning	Room temperature (cold)	Drinking water	--		Remove gross contamination
2	Drying	Room temperature	--	--		Lint-free cloth or medical-grade compressed air
3	Wiping Disinfection	--	--	Alcohol 70%	Aldehyde-free surface disinfectant	
4	Final Rinsing	Room temperature (cold)	Fully demineralized water	Not required	Fully remove all chemicals	Remove chemicals using demineralized water
5	Drying	Room temperature	--	--		Lint-free cloth

11.4 Mechanical Cleaning and Disinfection

For machine-based cleaning processes, use only pH-neutral or alkaline cleaning agents up to a pH of 10.5 which the respective manufacturer has expressly declared suitable for cleaning anodized aluminum components. To protect the aluminum surfaces from damage, no acidic neutralizers should be added. The final rinsing may be carried out using demineralized water only. Be sure to observe also the instructions provided by the manufacturer of your automatic cleaning machine as well as those provided by the manufacturer of the products used. Moreover, it is essential to use only cleaning machines that are suitable for cleaning containers, allowing safe placement inside the tray baskets and providing adequately arranged spray nozzles or arms. For more information visit <http://kls-martin.com/phnet/>

11.5 Removing Residues

Stubborn residues or dirt on the container box that cannot be removed during the usual cleaning process (such as adhesive labels, indicator strips, lettering) can be removed with commercial cleaners for anodized surfaces. Following this special treatment, both the container boxes and covers should be either rinsed for a second time or subjected to routine cleaning.

12 Maintenance and Repair

The manufacturer will assume responsibility for the safety and reliability of the Sterile Containers only if repairs or modifications are carried out by authorized persons and if the containers are used as intended and in accordance with the present Instructions for Use. Whenever parts affecting microbial impermeability are repaired, a careful visual inspection must be carried out before putting the container back in service. Only genuine Gebrüder Martin spare parts may be used for repair. Use of third-party spare parts shall void the warranty.

NOTICE

The container may be repaired only by Gebrüder Martin or by a person or firm that has been specially authorized by Gebrüder Martin to carry out such work.

If the repair is carried out by a person or firm duly authorized by Gebrüder Martin, the operator of the container is required to obtain from the repairer a certificate with details on the nature and scope of the work done. This certificate must show the name and address of the firm and must be dated and signed. Moreover, the repaired container or container parts must be marked with the repairer's ID label in such a case (i.e. whenever the repair has not been performed by the manufacturer itself).

12.1 Maintenance Schedule

Depending on the conditions of use, we recommend users to perform the following maintenance steps in compliance with EN 868-8:2009 after 500 usage cycles or once a year. All maintenance is, however, limited to checks of daily routine.

Component	Check for	How?	Remedy
Seal	Integrity	Visual inspection for damage	Replacing the container cover
Container cover	Integrity	Visual inspection for damage and deformation	Replacing the container cover
Filter holder	Integrity	Visual inspection of contour and locking	Replacing the filter holder
Container tray	Integrity	Visual inspection of contour and locking with the container cover	Repair or replacement
Sealing surface	Integrity	Visual inspection for planarity and damage	Replacement

13 Dimensions

The container dimensions fully comply with EN 868-8:2009 recommendations. For dimensional details, please refer to Gebrüder Martin's special brochure.

14 Safety Instructions

The contents of a Sterile Container can be considered sterile only if the container has been sterilized correctly, all locks are intact and all container seals are undamaged.

⚠ CAUTION

Danger of burns!

The container is hot after sterilization. Ensure adequate protective clothing is worn!

To avoid risks and damage, container stacks should not exceed a height limit of 55 cm (which corresponds to 2 STUs). Always use suitable facilities for transporting and storing the containers, taking care that they cannot topple over or drop down. Be sure to load the containers only with resterilizable items, thereby observing the product manufacturers' instructions in each case. The container is not suitable for sterilizing medicines or biological products.

15 Spare Parts

Spare part	Item no.
Filter holder	04-250-00-04
Retaining bracket for log labels	04-275-00-04

For further accessory equipment, please refer to the special brochure by Gebrüder Martin.

NOTICE

Only genuine and approved Gebrüder Martin parts may be used!

16 General Information

It is mandatory to observe all existing general guidelines and hygiene principles when handling contaminated items, items to be sterilized and sterile supplies. To avoid potential recontamination risks, carefully read and follow these Instructions for Use. Whenever deviations from the "to-be" status described above are detected when opening a container, its contents must be considered non-sterile and the entire sterilization process must be repeated. The life of the container and its seals complies with the requirements set out in the European standard EN 868-8:2009.

NOTICE

Only trained personnel may use the containers.

17 Reference Literature and Standards

EN 285:2013, ISO 17665-1:2006, ISO 11607-1:2009 and ISO 11607-2:2006, EN 868-8:2009, DIN 58952-3:2012, DIN 58953-9:2010.

18 Manufacturer's Information acc. to EN 868-8:2009, 5

Specification of major components	Section 5 "Structure and Function of the Container", page 7 and section 2 "Product Overview", page 7
Checks, maintenance and/or replacements	Section 6, 12, 14 and 15
Container life (sterilization cycles)	Acc. EN 868-8:2009, 4.3
Seal life	Acc. EN 868-8:2009, 4.3
Cleaning procedures	Section 11 "Disinfecting, Cleaning and Care", page 8