EMSCULPT

18 Jul 2018

FIRST & ONLY FOR MUSCLE & FAT

EMSCULPT®

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More than 53 offices
around the globe

1,500 employees
worldwide


More than
300 engineers

2 focused divisions:
medical and aesthetics

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*The American Society for Aesthetic Plastic Surgery

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
Robert Weiss MD FAAD FACPh
Former ASLMS President
Maryland Laser Skin & Vein Institute, Hunt Valley MD

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Carolyn Jacob MD, FAAD
Board-Certified Dermatologist
Chicago Cosmetic Surgery and Dermatology, Chicago IL

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
Brian M. Kinney, MD FACS
Plastic Surgery Excellence, Beverly Hills CA

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
Another patient dies at a Miami cosmetic surgery clinic

BY DANIEL CHANG

dchang@miamiherald.com

December 15, 2017 07:58 PM

Updated December 16, 2017 12:02 AM

A 40-year-old woman from Louisiana on Thursday became at least the fourth patient to die since 2013 from complications of cosmetic surgery at a South Florida clinic — a trend that has raised alarm among public health officials about the dangers of popular procedures such as liposuction and the so-called Brazilian butt lift.

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EFFECT ON FAT

Concise MoA Paper

EMSCULPT® TECHNOLOGY FOR NON-THERMAL INDUCTION OF MUSCLE GROWTH AND FAT APOPTOSIS

MECHANISM OF ACTION

The role of fat and muscles in aesthetic appearance. The majority of human body composition comprises of fat (approx. 25%) and muscle (42% male/36% female) tissues. Many procedures address subcutaneous fat which correlates with the overall body shape, and its reduction can deliver a slimmer look. The underlying muscles are however equally important as an increased muscle tone better defines the body contour by reducing the localized prolapse/fatly and adds to a healthier aesthetic appearance.

HIFEM® technology

EMSCULPT is based on High-Intensity Focused Electro-Magnetic (HIFEM) field technology which has the ability to induce supramaximal muscle contractions. The rapidly changing magnetic field induces electric currents in the tissue where it depolarizes neural membranes and governs motor units in the target muscles, causing concentric contractions.^{1,2} The effects are highly selective due to its physiological characteristics only motor neurons are activated, while other neurons or tissues are not responsive to the current, and therefore stay unaffected.

Nerve depolarization for induction of supramaximal contractions

During normal voluntary muscle contractions, the muscle fibers relax between each nervous stimulus due to the central nervous system's inability to signal another impulse while the previous one is still in action. EMSCULPT generates impulses that are independent of brain function, and at such a rapid frequency that

doesn't allow such a relaxation phase. Under normal conditions, the highest amount of tension that could be developed and held physiologically is called maximal voluntary contraction (MVC). Usually, it lasts only for a split second. Contractions with a tension higher than MVC are defined as supramaximal. EMSCULPT has the ability to generate supramaximal contractions and hold them for multiple seconds, which significantly increases the physiologic stress/workload needed to allow muscles to adapt.

The effects on muscle tissue

When exposed to supramaximal contractions, the muscle tissue is forced to adapt to such extreme conditions and responds with a deep remodeling of its inner structure, i.e., the growth of myofibrils (muscle hypertrophy) and creation of new protein strands and muscle fibers (muscle hyperplasia).^{3,4} Increased muscle density and volume lead to a better definition and muscle tone.

BEFORE AFTER



Figure 2: Illustration of muscle structure changes after exposure to extreme load.

The effects on fat tissue

During physical activity, muscles need energy to produce contractions. The energy is derived primarily from adenosine triphosphate (ATP) and secondarily from creatine phosphate and glycogen to fuel the muscles. When these are depleted, the body's catabolic processes take place in the form of lipolysis - i.e. the breakdown of lipids (triglycerides) into free fatty acids (FFA) and glycerol.^{5,6} These

Figure 1: Muscle response to EMSCULPT treatment compared to common exercise.

EMSCULPT.COM | INFO@BTNET.COM

EMSCULPT®

EMSCULPT®

Scientific MoA Review

Supramaximal deep muscle remodeling and ER stress-induced fat disruption using the Emsculpt High Intensity Focused Electro-magnetic field device

Robert Weiss MD¹, Brian Kinney MD²

¹Maryland Dermatology Laser Skin & Vein Institute, Hunt Valley, MD, USA

²Plastic Surgery Excellence, Beverly Hills, CA, USA

Abstract

Emsculpt is a high intensity focused electromagnetic-magnetic field device, causing a unique combination of non-invasive subcutaneous fat reduction and deep muscle remodeling that can achieve strong aesthetic improvement of the treated area. The device induces supramaximal muscle contractions via its High Intensity Focused Electro-Magnetic (HIFEM) field. This results in cell membrane porosity changes, electroporation and an outflow of free fatty acids (FFAs) from adipocytes. The fat cells cannot process such amount of FFAs released during the contractions, and they thus induce apoptosis via a chemical reaction called Endoplasmic Reticulum stress ("ER stress"). At the same time, the muscle is forced to adapt to the supramaximal muscle contractions which results in a highly effective deep muscle remodeling through muscle hypertrophy and muscle fiber hyperplasia. The introduced mechanisms lead to bulking and firm toning of the treated area.

Keywords: apoptosis; deep muscle stimulation; HIFEM; ER stress; supramaximal contractions; hypertrophy, hyperplasia

1. Introduction

1.1. Induction of muscle contractions via nerve pathways

The unit induces supramaximal muscle contractions by generating a High Intensity Focused Electro-Magnetic field (HIFEM). This field has an ability to selectively interact with peripheral motor neurons.¹⁻⁴ The interaction between the HIFEM fields and excitable nerve tissue creates electric currents in the neurons. Electric currents depolarize neural membranes in order to initiate action potentials spreading down the peripheral nerve, and governing motor units in the target muscles. Contractions mediated through the nerve allow the impulse to activate the whole muscle, encompassing the deep and

Due to its physiological characteristics, other tissue is less responsive to the current, and therefore stays relatively unaffected. Emsculpt's field has the ability to penetrate deep into the tissue. This allows for fully effective activation of all motor nerves in the treated area in a wide range of patient profiles with different physiological characteristics.

1.2. Unique focused technology

The magnetic field is targeting directly the fibers of peripheral motoric nerves⁷ in the stimulated area, thus leads to a contraction of the whole muscle group innervated by the specific nerve or nerve plexus. Also, motor nerve cells respond much more selectively to the HIFEM. Importantly, it selectively activates larger fibers, more than the smaller fibers that mediate pain.

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

EMSCULPT®

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WIKIPEDIA
The Free Encyclopedia

- Main page
- Contents
- Featured content
- Current events
- Random article
- Donate to Wikipedia
- Wikipedia store

Interaction

- Help
- About Wikipedia
- Community portal
- Recent changes
- Contact page

Tools

- What links here
- Related changes
- Upload file
- Special pages
- Permanent link
- Page information

Not logged in | Talk | Contributions | Create account | Log in

Article | Talk | Read | Edit | View history |

High-intensity focused electromagnetic field

From Wikipedia, the free encyclopedia
(Redirected from HIFEM)

This article is an **orphan**, as no other articles link to it. Please introduce links to this page from related articles; try the Find link tool for suggestions.
(July 2018)

High-Intensity focused electromagnetic field (HIFEM®) is a **non-invasive** medical technology that is used for strengthening^{[1][2]} and re-education^{[3][4]} of **muscles** via interaction of the **magnetic field** with the tissue of the patient. In **aesthetic medicine** it is also used for non-invasive reduction of **abdominal fat**^{[1][2][5][6]}.

HIFEM is the label of the technology that is used in **aesthetic medicine**, **urology** and **gynecology**. It uses **focused electromagnetic field** with intensity measured in **Tesla**^[7] and is based on **Faraday's principle** of **electromagnetic induction**. **Electromagnetic field** passes non-invasively through the body and interacts with **motor neurons** which subsequently trigger **supramaximal muscle contractions** due to the **action potential**^[1]. The exposure of muscles to these contractions leads to muscle strengthening.

In **aesthetic medicine**, HIFEM is used as a **non-invasive body contouring** treatment for **abdomen** and **buttocks**^{[1][8][9]}. The same technology is used in **urology** and **gynecology** to treat **urinary incontinence** through strengthening of the **pelvic floor muscles**^{[4][10][11]}.

Contents [hide]

- 1 Mechanism of action
- 2 Use
 - 2.1 Effect on muscles
 - 2.2 Effect on fat
 - 2.3 Effect on pelvic floor
- 3 Regulatory approval

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FIRST & ONLY FOR MUSCLE & FAT

EMSCULPT®

CLINICAL TRAINING VIDEOS AVAILABLE

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*Study in progress, displayed results are for the first 16 patients (ASLMS initial abstract)

FIRST & ONLY FOR MUSCLE & FAT

EMSCULPT®

RELATED RESEARCH AVAILABLE ONLINE

American Society for Laser Medicine and Surgery Abstracts

369

AMERICAN SOCIETY FOR LASER MEDICINE AND SURGERY

2018 ELECTRONIC POSTERS (ePOSTERS) TOWN HALL AND ePOSTERS

ePOSTER TOWN HALL: *BODY CONTOURING* CARBOXYTHERAPY FOR SUBCUTANEOUS

A NOVEL NON-INVASIVE
ON SIMULTANEOUS INDU
IN ADIPOSE AND MUSCLE
AND EFFICACY OF A HIGH
FOCUSED ELECTRO-MAG
DEVICE USED FOR ABDOM
SHAPING

Carolyn Jacob, Katya Pasko
*Chicago Cosmetic Surgery and Der
Derma Vita Clinic, Sofia, Bulgaria*
Background: Current approaches
contouring are represented by fat
inducing thermal effects (radiofreq
These are ideal for patients with fi
Our study investigates effects of a
both the subcutaneous adipocytes
structure in a non-thermal manne
lower-to-medium BMI patients.

Study Design/Materials and Me
23.8 kg·m²) received four 30-min t

NIH U.S. National Library of Medicine

ClinicalTrials.gov

[Find Studies](#) ▾

[About Studies](#) ▾

[Submit Studies](#) ▾

[Resources](#) ▾

[About Site](#) ▾

[Home](#) > [Search Results](#) > [Study Record Detail](#)

[Save this study](#)

Trial record **6 of 14** for: BTL [Lead]

[Previous Study](#)

[Return to List](#)

[Next Study](#)

Changes in Abdominal Tissues Induced by Magnetic Device - MRI Evaluation

The safety and scientific validity of this study is the responsibility of the study sponsor and investigators. Listing a study does not mean it has been evaluated by the U.S. Federal Government. Read our [disclaimer](#) for details.

ClinicalTrials.gov Identifier: NCT03545633

Recruitment Status ⓘ: Active, not recruiting

First Posted ⓘ: June 4, 2018

Last Update Posted ⓘ: June 6, 2018

Sponsor:

BTL Industries Ltd.

Information provided by (Responsible Party):

BTL Industries Ltd.

FIRST & ONLY FOR MUSCLE & FAT

EMSCULPT®



HISTOLOGICAL IN VIVO STUDY: THE MECHANISM OF ACTION

INDUCTION OF FAT APOPTOSIS BY NON-THERMAL DEVICE: SAFETY AND MECHANISM OF ACTION OF NON-INVASIVE HIFEM TECHNOLOGY EVALUATED IN A HISTOLOGICAL PORCINE MODEL. R. Weiss, M.D., J. Bernardy, MVDr.

92%

Increase in apoptotic index
8 hours post-treatment

Microscopic analysis of the fat tissue confirmed that the amount of apoptotic cells increased significantly after the treatments.



FIRST & ONLY FOR MUSCLE & FAT

EMSCULPT®

ULTRASONOGRAPHY STUDY

CHANGES IN SUBCUTANEOUS ABDOMINAL FAT THICKNESS FOLLOWING HIGH-INTENSITY FOCUSED ELECTRO-MAGNETIC (HIFEM) FIELD TREATMENTS: A MULTI CENTER ULTRASOUND STUDY.

Bruce Katz MD¹, Robert Bard MD², Richard Goldfarb MD³, Aaron Shiloh MD⁴, Dilyana Kenolova MD⁵

FIRST & ONLY FOR MUSCLE & FAT

EMSCULPT®

MAGNETIC RESONANCE STUDY

HIGH INTENSITY FOCUSED ELECTRO-MAGNETIC THERAPY (HIFEM) EVALUATED BY MAGNETIC RESONANCE IMAGING (MRI): SAFETY AND EFFICACY STUDY OF A DUAL TISSUE EFFECT BASED ON NON-INVASIVE ABDOMINAL BODY SHAPING. Brian M Kinney MD FACS¹, Paula Lozanova MD²

3.6_{cm}

Average
circumferential
reduction

18.6%

Average fat
reduction

15.4%

Average increase
in Muscle mass




10.4%

Average reduction
in Diastasis Recti

BEFORE

2-MONTH FU

FIRST & ONLY FOR MUSCLE & FAT

EMSCULPT®

COMPUTED TOMOGRAPHY STUDY

COMPUTED TOMOGRAPHY BASED EVIDENCE OF SIMULTANEOUS CHANGES IN HUMAN ADIPOSE AND MUSCLE TISSUES FOLLOWING A HIGH INTENSITY FOCUSED ELECTRO-MAGNETIC APPLICATION: A NEW METHOD FOR NON-INVASIVE BODY SCULPTING. David E. Kent MD¹, Carolyn I. Jacob MD²

4.17 cm

Average
circumferential
reduction

19%

Average fat
reduction


16%

Average increase
in Muscle thickness


11%

Average reduction
in Diastasis Recti

BEFORE

1-MONTH FU

SAME PATIENT AS ABOVE

FIRST & ONLY FOR MUSCLE & FAT

EMSCULPT®

WAIST CIRCUMFERENCE REDUCTION TESTED IN A MULTICENTER STUDY

A NOVEL NON-INVASIVE TECHNOLOGY BASED ON SIMULTANEOUS INDUCTION OF CHANGES IN ADIPOSE AND MUSCLE TISSUES: SAFETY AND EFFICACY OF A HIGH INTENSITY FOCUSED ELECTRO-MAGNETIC FIELD DEVICE USED FOR ABDOMINAL BODY SHAPING. Carolyn I. Jacob MD¹, Katya Paskova MD²

4.4_{cm}

Average waist reduction on
avg. BMI 23.8

95%

Satisfaction with treatment
results

89%

BA recognition
by blinded evaluators

HIGHER-BMI PATIENT RESULTS

LOWER-BMI PATIENT RESULTS

BEFORE

3-MONTH FU

FIRST & ONLY FOR MUSCLE & FAT

EMSCULPT®

A LARGE-SCALE MULTICENTRIC STUDY: NON-INVASIVE BUTTLIFTING EFFECTS

HIGH INTENSITY FOCUSED ELECTRO-MAGNETIC TECHNOLOGY (HIFEM) FOR NON-INVASIVE BUTTOCKS
LIFTING AND TONING OF GLUTEAL MUSCLES: EFFICACY AND SAFETY STUDY.

C. Jacob MD¹, B. Kinney MD², M. Busso MD³, S. Chilukuri MD⁴, JD McCoy NMD⁵, C. Bailey⁶, R. Denkova MD⁷

BEFORE

1-MONTH FU

BEFORE

1-MONTH FU

75

Patients multi-center
study

85%

Reported a significant
improvement

FIRST & ONLY FOR MUSCLE & FAT

EMSCULPT®

FIRST & ONLY FOR MUSCLE & FAT

EMSCULPT®

SUMMARY

FIRST & ONLY FOR MUSCLE & FAT


EMSCULPT®

NON-INVASIVE

**7 PEERREVIEWED
STUDIES**

**MUSCLE & FAT
(NEW PATIENTS)**

NON-THERMAL

NO CONSUMABLES

**OPERATOR
INDEPENDENT**

NEW PATIENTS

**NON-INVASIVE BUTTOCK
PROCEDURE**

NO DOWNTIME

FIRST & ONLY FOR MUSCLE & FAT

EMSCLPT®

Q&A

FIRST & ONLY FOR MUSCLE & FAT

EMSCULPT®