Tnl

Troponin I (CLIA)

■ Order Information

Catalog No.	Package Size
105-005659-00	2×50 tests
105-005676-00	2×100 tests

■ Intended Use

The Mindray CL-series Tnl assay is a Chemiluminescent Immunoassay (CLIA) for the quantitative determination of Tnl in human serum or plasma. This assay can be used as an aid in the differential diagnosis of acute coronary syndrome (ACS), especially for the timely management of acute myocardial infarction (AMI).

■ Summary

Troponin-I (Tnl) is a regulatory subunit of the troponin complex which consists of three subunits, i.e. Troponin I, Troponin T, and Troponin C. Troponin I can form troponin-I-C together with Troponin-C, and it can also form troponin-I-T with Troponin-T. The main physiological action of Troponin I is inhibiting actomyosin ATPase to impede muscle contraction due to lack of calcium.

Following myocardial infarction (MI) or ischemic damage, the integrity of myocardium cytomembrane will be destroyed and many structure proteins with other large molecules are released into cardiac muscle mesenchyme within hours¹. These biomarkers of necrosis include cardiac Troponin I, Troponin T, CK-MB, myoglobin, and so on. Cardiac Troponin I clearly differs from skeletal muscle Troponin I, which can be used to differentiate between skeletal muscle lesions and myocardial injury. Compared with other biomarkers, Cardiac Troponin I is the first choice in biomarkers for myocardial damage due to its high sensitivity and high tissue-specificity². Thus, Cardiac Troponin I plays an important role in the supplementary diagnosis of MI and identifying patients with acute coronary syndromes (ACS) who are at greater risk for cardiac events³.

According to redefinition by European Society of Cardiology (ESC) and American College of Cardiology (ACC), MI is diagnosed when blood levels of cardiac troponin are above the 99th percentile of the healthy population reference limit in the clinical setting of acute ischemia⁴. The coefficient of variation at the 99th percentile for troponin assays must be less than or equal to 10%. Clinical studies for patients with acute myocardial infarction (AMI) have demonstrated elevated levels of cardiac Troponin I are detectable in serum within 3 to 6 hours after the onset of chest pain, they will reach peak concentrations in approximately 12 to 16 hours, and then remain elevated for 4-9 days^{5,6}.

Elevated cardiac Troponin I levels have also been reported in cases of unstable angina pectoris (UAP) and congestive heart failure (CHF). Other studies have shown that detectable levels of cardiac Troponin I might correlate with higher incidence of mortality in patients with UAP and non-ST segment elevation MI (NSTEMI)⁷. Thus, the measurement of cardiac Troponin I can be useful in the risk stratification of patients with UAP and NSTEMI.

■ Assay Principle

The Mindray CL-series Tnl assay is a two-site immunoenzymatic assay to determine the level of Tnl.

In the first step, sample, pretreatment solution, paramagnetic microparticles coated with monoclonal anti-Tnl antibody (mouse) and monoclonal anti-Tnl antibody (mouse)-alkaline phosphatase conjugate are added into a reaction cuvette. After incubation, Tnl present in the sample binds to both anti-Tnl antibody coated microparticles and anti-Tnl antibody alkaline phosphatase-labeled conjugate to form a sandwich complex. Microparticles are magnetically captured while other unbound substances are removed by washing.

In the second step, the substrate solution is added to the reaction Cuvette. It is catalyzed by anti-Tnl antibody (mouse)-alkaline phosphatase conjugate in the immunocomplex retained on the microparticles. The resulting chemiluminescent reaction is measured as Relative light units (RLUs) by a

photomultiplier built inside the system. The amount of Tnl present in the sample is proportional to the Relative light units (RLUs) generated during the reaction. The Tnl concentration can be determined via a calibration curve.

■ Reagent Components

Ra	Paramagnetic microparticles coated with monoclonal anti-Tnl antibody (mouse). Minimum concentration: 1.2 g/L solid. TRIS [®] buffer: 50 mmol/L. Preservatives: 0.05% ProClin 300 and 0.09% sodium azide.
Rb	Anti-Tnl antibody (mouse)-alkaline phosphatase conjugate. Minimum concentration: 2.08 µg/mL. MES [®] buffer: 50 mmol/L. Preservatives: 0.05% ProClin 300 and 0.09% sodium azide.
Rc	Pretreatment solution. TRIS buffer: 50 mmol/L. Preservatives: 0.05% ProClin 300 and 0.09% sodium azide.

a) TRIS = Tris(hydroxymethyl)-aminomethane

b) MES = 2-(N-morpholino)ethanesulfonic acid buffered saline

■ Storage and Stability

The unopened Tnl (CLIA) reagent kit is stable up to the stated expiration date when stored at 2-8°C.

The Tnl (CLIA) reagent kit can be stored onboard and used for a maximum of 28 days after opening at 2-8°C.

■ Reagent Preparation

The reagents in the kit have been assembled into a ready-for-use unit that cannot be separated.

■ Materials Required but not Provided

Mindray CL-series Chemiluminescence Immunoassay Analyzer.
Cat.No.105-005910-00: Troponin I Calibrators, 1×2.0 mL for each of calibrator C0, C1 and C2.

Cat.No.105-005941-00: Cardiac Marker Multi Control (L), 3×2.0 mL.

Cat.No.105-005942-00: Cardiac Marker Multi Control (H), 3×2.0 mL.

Cat.No.105-005927-00: Cardiac Marker Multi Control (L), 6×2.0 mL.

Cat.No.105-005928-00: Cardiac Marker Multi Control (H), 6×2.0 mL.

Cat.No.105-004452-00: Wash Buffer, 1×10 L.

Cat.No.105-009044-00: Substrate Solution, 4×75 mL.

Cat.No.105-004274-00: Substrate Solution, 4×115 mL.

Reaction Cuvette.

■ Applicable Instrument

Mindray CL-series Chemiluminescence Immunoassay Analyzer

■ Specimen Collection and Preparation

Specimen Types

• Human serum and plasma (K₂EDTA, K₃EDTA, sodium heparin and lithium heparin) specimens are recommended for this assay.

• Blood collection tubes from various manufacturers may contain additives which could affect the test results in some cases. Not all available tubes of on the market were tested by Mindray. Each laboratory should determine the acceptability of different blood collection tubes and serum/plasma separation products.

• When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Specimen conditions

• Do not use:
– heat-inactivated specimens
– grossly hemolyzed specimens
– specimens with obvious microbial contamination

• For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.

• Avoid more than five freeze cycles.

• Follow blood collection tube manufacturer's recommendations for centrifugation. Centrifuge the specimens after clot formation is complete. Ensure residual fibrin and cellular matter has been removed prior to analysis.

• For optimal results, inspect all samples for bubbles. Remove bubbles with a pipette tip prior to analysis. Specimens must be mixed thoroughly after thawing. Thawed samples should be centrifuged prior to use.

• If the sample was covered with lipid layer after centrifugation, the sample should be transferred to a clean tube and centrifuged before testing. Do not transfer the lipid layer. Handle carefully to prevent cross contamination.

Specimen Storage

• Specimens should be tested timely after sample collection. If the assay cannot be completed within 8 hours, refrigerate the samples at 2 to 8°C. If testing will be delayed for more than 36 hours, specimens should be frozen at -20°C or colder. The specimens can be stored at -20°C for up to 30 days.

• Avoid more than five freeze cycles.

■ Assay Procedure

For optimal performance of this assay, operators should read the related system operation manual carefully, to get sufficient information such as operation instructions, sample preservation and management, safety precaution, and maintenance. Prepare all required materials for the assay as well.

Before loading the Tnl (CLIA) reagent kit on the machine for the first time, unopened reagent bottle should be inverted gently for at least 30 times to resuspend the microparticles that have settled during shipment or storage. Visually inspect the bottle to ensure the microparticles have been completely resuspended. If the microparticles remain adhered to the bottle, continue inverting until the microparticles have been completely resuspended. If the microparticles cannot be resuspended, it is recommended not to use this bottle of reagent. Contact Mindray Customer Service for help. Do not invert opened reagent bottle.

This assay requires 75 µL of sample for a single test. This volume does not include the dead volume of the sample container. Additional volume is required when performing additional tests from the same sample. Operators should refer to the system operation manual and specific requirement of the assay to determine the minimum sample volume.

■ Calibration

The Mindray CL-series Tnl (CLIA) has been standardized against a commercial Tnl test (CLIA).

The specific information of master calibration curve of Tnl (CLIA) reagent kit is stored in the two-dimensional barcode on the reagent pack, which is used in combination with Tnl calibrators for calibration of the specific reagent lot. Before initiating calibration on each new lot of reagent, load the assay master curve by scanning the two-dimensional barcode on the reagent pack. When performing the calibration, scan the two-dimensional barcode on the calibrator pack, and then test the Tnl calibrators at three levels. Valid calibration curve is required before any Tnl test. Recalibration is recommended every 4 weeks, or when a new reagent lot is used, or the quality controls are out of specified range. For detailed instruction of calibration, refer to the system operation manual.

■ Quality Control

It is recommended that quality controls should be run once every 24 hours if the tests are in use, or after every calibration. The quality control frequency should be adapted to each laboratory's individual requirements. The recommended two levels of quality controls for this assay are Cardiac Marker Multi Control (L) and Cardiac Marker Multi Control (H).

Quality control results should be within the acceptable ranges. If a control is out of its specified range, the associated test results are invalid and the samples must be retested. Recalibration may be required. Examine the assay system referring to the system operation manual. If the quality control results are still out of the specified range, please contact Mindray Customer Service for help.

■ Calculation

The analyzer automatically calculates the analyte concentration of each sample on the master calibration curve read from the barcode, and a 4-Parameter Logistic Curve Fitting (4PLC) with the relative light units (RLUs) generated from Tnl calibrators of defined concentration values. The results are shown in the unit of ng/mL.

■ Dilution

Samples with Tnl concentrations above the upper limit can be diluted with Mindray Sample Diluent. The recommended dilution is 1:10 (either automatically by the analyzer or manually). The concentration of the diluted sample must be >1 ng/mL. After manual dilution, multiply the result by the dilution factor. After automated dilution by the analyzers, the system automatically multiply the result by the dilution factor when calculating the sample concentration.

■ Expected values

The 99th percentile upper reference limit

A study on a cohort of 280 healthy individuals using serum samples has determined the reference range of Mindray CL-series Tnl assay.

Number of Samples	Age Range	99 th Percentile
280	19-80 years	0.04 ng/mL

Cutoff value for AMI

297 patients who have symptoms of chest pain were tested with Mindray CL-series Tnl assay. Based on the World Health Organization (WHO) criteria for the definition of AMI from the 1970's, 66 individuals were diagnosed with AMI. Using the clinical results combined several literature results⁴⁻⁶, we established Receiver Operating Characteristics (ROC) curve to determine the diagnostic cutoff and the Area under Curve (AUC) was greater than 0.90. The most appropriate decision threshold for determining AMI from these distributions is 0.5 ng/mL. We recommend the cutoff of CL-series Tnl assay for determining AMI is 0.5 ng/mL.

Due to the variation in geography, race, sex, and age, it is highly recommended that each laboratory should establish its own reference range.

■ Limitation

The upper limit of this assay is 50 ng/mL. A specimen with Tnl concentration lower than the upper limit can be quantitatively determined, while the specimen with a concentration higher than the upper limit will be reported as >50 ng/mL or diluting the samples with Mindray Sample Diluent.

The concentration of Tnl in a given specimen, determined with assays from different manufacturers, can vary due to differences in assay methods, calibration, and reagent specificity. The assay results should be used in conjunction with other evidences to make clinical decisions, such as symptoms, results of other tests.

Specimen from individuals who have been exposed to mouse monoclonal antibodies may contain human anti-mouse antibodies (HAMA)⁸. Such specimens may show either falsely elevated or falsely depressed values with assay kits employing mouse monoclonal antibodies^{9,10}. However, no obvious interference of HAMA has been observed in this assay.

■ Performance Characteristics

■ Lower limits of measurement

Limit of Blank, Limit of Detection and Limit of Quantitation

Limit of Blank (LoB) = 0.006 ng/mL

Limit of Detection (LoD) = 0.010 ng/mL

Limit of Quantitation (LoQ) = 0.015 ng/mL

The Limit of Blank, Limit of Detection and Limit of Quantitation were determined in accordance with the CLSI (Clinical and Laboratory Standards Institute) EP17-A2 requirements.

The Limit of Blank is the 95th percentile value from n ≥ 60 measurements of analyte-free samples over several independent series. The Limit of Blank corresponds to the concentration below which analyte-free samples are found with a probability of 95 %.

The Limit of Detection is determined based on the Limit of Blank and the standard deviation of low concentration samples.

The Limit of Detection corresponds to the lowest analyte concentration which can be detected (value above the Limit of

Blank with a probability of 95 %).

The Limit of Quantitation is the lowest analyte concentration with a total allowable relative error of $\leq 30\%$.

■ Functional Sensitivity

The Tnl (CLIA) reagent assay has a functional sensitivity of ≤ 0.04 ng/mL, which meets the requirements of Tnl assay. Functional sensitivity is defined as the concentration of Tnl that can be measured with an interassay CV of 10%.⁵

■ Measuring Range

The measuring range is defined as the range of values which meets the limits of acceptable for linearity. The measuring range of Tnl (CLIA) reagent assay is 0.006-50 ng/mL. Values below the analytical sensitivity are reported as <0.006 ng/mL. Values above the measuring range are reported as >50 ng/mL (or up to 500 ng/mL for 10-fold diluted samples).

■ Analytical Specificity

Hemoglobin up to 500 mg/dL, bilirubin up to 20 mg/dL, triglycerides up to 1000 mg/dL, total protein up to 10 g/dL, rheumatoid factors (RF) up to 1500 IU/mL and antinuclear antibody will not interfere with the CL-series Tnl assay. Criterion: Recovery within $\pm 10\%$ (RF: $\pm 15\%$) of initial value.

In vitro tests were performed on 24 commonly used pharmaceuticals. No interference was observed from these substances. Criterion: Recovery within $\pm 10\%$ of initial value.

Tnl Calibrator C0 was spiked with Skeletal Troponin I, Cardiac Troponin C, Recombinant human Cardiac Troponin T, Actin, Tropomyosin, Myoglobin, Myosin Light Chain, Human CK-MB. No obvious cross reactivity was observed as all the results were $\leq 0.05\%$. The results are summarized in the table below.

Substance	Cross-reactant Concentration (ng/mL)	Cross Reactivity	Acceptance Criteria
Skeletal Troponin I	1000	0.00%	Reported Cross Reactivity $\leq 0.05\%$
Cardiac Troponin C	1000	0.00%	
Recombinant human Cardiac Troponin T	1000	0.00%	
Actin	1000	0.00%	
Tropomyosin	1000	0.00%	
Myoglobin	1000	0.00%	
Myosin Light Chain	1000	0.00%	
Human CK-MB	1000	0.00%	

*Representative data; results in individual laboratories may vary.

■ High Dose Hook

For the Mindray CL-series Tnl assay, there is no high dose hook effect at Tnl concentration up to approximately 1000 ng/mL.

■ Accuracy

Two controls with traceable and predefined values were used to verify the accuracy of this assay. The results showed that the relative deviations were within $\pm 10\%$. The results are summarized in the table below.

Sample	Measured Tnl Value (ng/mL)	Defined Tnl Value (ng/mL)	Relative Deviation
Level 1	0.49	0.51	-4.70%
Level 2	34.36	35.95	-4.41%

*Representative data; results in individual laboratories may vary.

■ Precision

Precision was determined by following Clinical and Laboratory Standards Institute (CLSI) Protocol EP5-A2¹¹. Two levels of quality controls and three levels of human serum were tested in duplicate in two separate runs per day, for a total of 20 days, using a single lot of reagents and a single calibration curve. The precision data are summarized in the table below.

Sample	Mean Tnl (ng/mL)	Within-run CV	Between-run CV	Within-Device CV
Control L	0.32	2.01%	0.90%	2.59%
Control H	14.38	0.90%	1.19%	1.80%

Sample	Mean Tnl (ng/mL)	Within-run CV	Between-run CV	Within-Device CV
HS-1	0.50	1.21%	1.11%	2.00%
HS-2	9.99	1.06%	1.70%	2.75%
HS-3	34.81	1.40%	1.29%	2.04%

*Representative data; results in individual laboratories may vary.

■ Linearity

A study was performed based on guidance from CLSI EP06-A¹². A high concentration Tnl sample (approximately 50 ng/mL) was mixed with a low concentration sample (<0.006 ng/mL) at different ratios, generating a series of dilutions. The Tnl of each dilution was determined using the Mindray CL-Series Tnl Assay. Linearity was demonstrated in the range of 0.006 ng/mL to 50 ng/mL, the correlation coefficient r is ≥ 0.9900 . The linearity data are summarized in the table below.

Concentration (ng/mL)	1	2	3	4	5	6
Expected Tnl	0.00	11.41	22.83	34.24	45.65	57.06
Measured Tnl	0.01	12.28	24.13	35.55	46.43	57.06

*Representative data; results in individual laboratories may vary.

■ Method Comparison

The Mindray CL-Series Tnl assay was compared to a commercially available diagnostic kit in a correlation study with about 212 serum specimens. The statistical data obtained by Deming computing mode are summarized in the table below.

Concentration Range	Slope	Intercept	Correlation Coefficient
0.006 - 50 ng/mL	1.0032	0.00005	0.9978

■ Warnings and Precautions

- For in vitro diagnostic use only. For laboratory professional use.
- Follow all the rules in handling laboratory reagents and take necessary safety precautions.
- The concentration of Tnl in a given specimen determined with different manufacturers can vary due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the Tnl assay used. Values obtained with different assay methods cannot be used interchangeably.
- Do not use reagent kits beyond the expiration date.
- Do not use reagents mixed from different reagent lots.
- Always keep the reagent pack in the upright position to ensure no microparticle has been lost prior to use.
- Reagent pack opened for more than 28 days is not recommended for use.
- Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
- All the specimen and reaction wastes should be considered potentially biohazard. The handling of specimens and reaction wastes should be in accordance with the local regulations and guidelines.
- The Material Safety Data Sheet (MSDS) is available upon request.
- Please confirm the integrity of the package before use. Do not use the reagents with damaged packages.
- If the reagents are opened unintentionally before use, they shall be used as soon as possible.
- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the local competent authority.
- Instability or deterioration should be suspected if there are visible signs of leakage, turbidity, precipitates or microbial growth.
- Do not freeze. The results can't be assured when the reagents are stored at inappropriate condition.
- This kit contains components classified as follows in accordance with the Regulation (EC) No. 1272/2008:

Warning

H317 May cause an allergic skin reaction.

H412 Harmful to aquatic life with long lasting effects.

Prevention:

P261 Avoid breathing dust /fume/gas/mist/vapours/spray.
P272 Contaminated work clothing should not be allowed out of the workplace.

P280 Wear protective gloves/protective clothing/eye protection/face protection.
P273 Avoid release to the environment.

Response:

P302 + P352 IF ON SKIN: Wash with plenty of water.
P333 + P313 If skin irritation or rash occurs: Get medical advice/attention.
P362 + P364 Take off contaminated clothing and wash it before reuse.

Disposal:

P501 Dispose of contents/container in accordance with local regulation.

■ Graphical Symbols

■ References

- Mair J, Wagner I, et al. Cardiac troponin I to diagnose myocardial injury. Lancet 1993; 341: 838-839.
- Bodor GS. Cardiac troponin-I: a highly specific biochemical marker for myocardial infarction. J Clin Immunoassay 1994; 17:40-44.
- Antman EM, Tanasijevic MJ, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. N Engl J Med 1996; 335:1342-1349.
- The Joint European Society of Cardiology/American College of Cardiology Committee. Myocardial Infarction redefined—a consensus document of the Joint European Society of Cardiology/American College of Cardiology Committee for the Redefinition of Myocardial Infarction. J Am Coll Cardiol 2000; 36: 959-969.
- Antman EM, Tanasijevic MJ, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. N Engl J Med 1996; 335:1342-1349.
- Tanasijevic MJ, Cannon CP, et al. The role of cardiac troponin I in risk stratification of patients with unstable coronary artery disease. Clin Cardiol 1999; 22:13-16.
- Antman EM, Fox KM. Guidelines for the diagnosis and management of unstable angina and non-Q-wave myocardial infarction: Proposed revisions. Am Heart J 2000; 139: 461-75.
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
- Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
- CLSI. EP5-A2: Vol. 24, No. 25, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline – Second Edition.
- CLSI. EP06-A: Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline. 2003.

© 2015-2022 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights Reserved

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R.China.

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-81888998

Fax: +86-755-26582680

EC-Representative: Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffeustraße 80, Hamburg 20537, Germany

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- Mair J, Wagner I, et al. Cardiac troponin I to diagnose myocardial injury. *Lancet* 1993; 341: 838-839.
- Bodor GS. Cardiac troponin-I: a highly specific biochemical marker for myocardial infarction. *J Clin Immunoassay* 1994; 17:40-44.
- Antman EM, Tanasijevic MJ, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. *N Engl J Med* 1996; 335:1342-1349.
- The Joint European Society of Cardiology/American College of Cardiology Committee. Myocardial infarction redefined—a consensus document of the Joint European Society of Cardiology/American College of Cardiology Committee for the Redefinition of Myocardial Infarction. *J Am Coll Cardiol* 2000; 36: 959-969.
- Antman EM, Tanasijevic MJ, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. *N Engl J Med* 1996; 335:1342-1349.
- Tanasijevic MJ, Cannon CP, et al. The role of cardiac troponin I in risk stratification of patients with unstable coronary artery disease. *Clin Cardiol* 1999; 22:13-16.
- Antman EM, Fox KM. Guidelines for the diagnosis and management of unstable angina and non-Q-wave myocardial infarction: Proposed revisions. *Am Heart J* 2000; 139: 461-75.
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988; 34:27-33.
- Kricka L. Interferences in immunoassays - still a threat. *Clin Chem* 2000; 46: 1037-1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. *Clin Chem* 2002; 48: 613-621.
- CLSI. EP5-A2: Vol. 24, No. 25, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline - Second Edition.
- CLSI. EP06-A: Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline. 2003.

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Identificador exclusivo do dispositivo

■ Referências

- Mair J, Wagner I, et al. Cardiac troponin I to diagnose myocardial injury. Lancet 1993; 341; 838-839.
- Bodor GS. Cardiac troponin-I: a highly specific biochemical marker for myocardial infarction. J Clin Immunoassay 1994; 17:40-44.
- Antman EM, Tanasiewicz MJ, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. N Engl J Med 1996; 335:1342-1349.
- The Joint European Society of Cardiology/American College of Cardiology Committee. Myocardial Infarction redefined—a consensus document of the Joint European Society of Cardiology/American College of Cardiology Committee for the Redefinition of Myocardial Infarction. J Am Coll Cardiol 2000; 36: 959-969.
- Antman EM, Tanasiewicz MJ, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. N Engl J Med 1996; 335:1342-1349.
- Tanasiewicz MJ, Cannon CP, et al. The role of cardiac troponin I in risk stratification of patients with unstable coronary artery disease. Clin Cardiol 1999; 22:13-16.
- Antman EM, Fox KM. Guidelines for the diagnosis and management of unstable angina and non-Q-wave myocardial infarction: Proposed revisions. Am Heart J 2000; 139: 461-75.
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
- Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
- CLSI. EP5-A2: Vol. 24, No. 25. Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline - Second Edition.
- CLSI. EP06-A: Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline. 2003.

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Fabricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Endereço: Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R. China.

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■ Referencias

- Mair J, Wagner I, et al. Cardiac troponin I to diagnose myocardial injury. Lancet 1993; 341; 838-839.
- Bodor GS. Cardiac troponin-I: a highly specific biochemical marker for myocardial infarction. J Clin Immunoassay 1994; 17:40-44.
- Antman EM, Tanasiewicz MJ, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. N Engl J Med 1996; 335:1342-1349.
- The Joint European Society of Cardiology/American College of Cardiology Committee. Myocardial Infarction redefined—a consensus document of the Joint European Society of Cardiology/American College of Cardiology Committee for the Redefinition of Myocardial Infarction. J Am Coll Cardiol 2000; 36: 959-969.
- Antman EM, Tanasiewicz MJ, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. N Engl J Med 1996; 335:1342-1349.
- Tanasiewicz MJ, Cannon CP, et al. The role of cardiac troponin I in risk stratification of patients with unstable coronary artery disease. Clin Cardiol 1999; 22:13-16.
- Antman EM, Fox KM. Guidelines for the diagnosis and management of unstable angina and non-Q-wave myocardial infarction: Proposed revisions. Am Heart J 2000; 139: 461-75.
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
- Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
- CLSI. EP5-A2: Vol. 24, No. 25, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline - Second Edition.
- CLSI. EP06-A: Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline, 2003.

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
■ Simboli grafici

■ Riferimenti

- Mair J, Wagner I, et al. Cardiac troponin I to diagnose myocardial injury. Lancet 1993; 341; 838-839.
- Bodor GS. Cardiac troponin-I: a highly specific biochemical marker for myocardial infarction. J Clin Immunology 1994; 17:40-44.
- Antman EM, Tanasiewicz MJ, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. N Engl J Med 1996; 335:1342-1349.
- The Joint European Society of Cardiology/American College of Cardiology Committee. Myocardial Infarction redefined—a consensus document of the Joint European Society of Cardiology/American College of Cardiology Committee for the Redefinition of Myocardial Infarction. J Am Coll Cardiol 2000; 36: 959-969.
- Antman EM, Tanasiewicz MJ, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. N Engl J Med 1996; 335:1342-1349.
- Tanasiewicz MJ, Cannon CP, et al. The role of cardiac troponin I in risk stratification of patients with unstable coronary artery disease. Clin Cardiol 1999; 22:13-16.
- Antman EM, Fox KM. Guidelines for the diagnosis and management of unstable angina and non-Q-wave myocardial infarction: Proposed revisions. Am Heart J 2000; 139: 461-75.
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
- Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.
- CLSI. EP5-A2: Vol. 20, No. 25. Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline – Second Edition.
- CLSI. EP06-A: Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline. 2003.

© 2015-2022 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Tutti i diritti riservati

Produttore: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Indirizzo: Mindray Building, Keji 12th Road South, High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R.China.

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Seri kodu

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■ Referanslar

- Mair J, Wagner I, et al. Cardiac troponin I to diagnose myocardial injury. Lancet 1993; 341: 838-839.
- Bodor GS. Cardiac troponin-I: a highly specific biochemical marker for myocardial infarction. J Clin Immunoassay 1994; 17:40-44.
- Antman EM, Tanasijevic MJ, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. N Engl J Med 1996; 335:1342-1349.
- The Joint European Society of Cardiology/American College of Cardiology Committee. Myocardial Infarction redefined—a consensus document of the Joint European Society of Cardiology/American College of Cardiology Committee for the Redefinition of Myocardial Infarction. J Am Coll Cardiol 2000; 36: 959-969.
- Antman EM, Tanasijevic MJ, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. N Engl J Med 1996; 335:1342-1349.
- Tanasijevic MJ, Cannon CP, et al. The role of cardiac troponin I in risk stratification of patients with unstable coronary artery disease. Clin Cardiol 1999; 22:13-16.
- Antman EM, Fox KM. Guidelines for the diagnosis and management of unstable angina and non-Q-wave myocardial infarction: Proposed revisions. Am Heart J 2000; 139: 461-75.
- Boscatto LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988; 34:27-33.
- Kricka L. Interferences in immunoassays - still a threat. Clin Chem 2000; 46: 1037-1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. Clin Chem 2002; 48: 613-621.

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
- Mair J, Wagner I, et al. Cardiac troponin I to diagnose myocardial injury. *Lancet* 1993; 341; 838-839.
- Bodor GS. Cardiac troponin-I: a highly specific biochemical marker for myocardial infarction. *J Clin Immunoassay* 1994; 17:40-44.
- Antman EM, Tanasijevic MJ, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. *N Engl J Med* 1996; 335:1342-1349.
- The Joint European Society of Cardiology/American College of Cardiology Committee. Myocardial Infarction redefined—a consensus document of the Joint European Society of Cardiology/American College of Cardiology Committee for the Redefinition of Myocardial Infarction. *J Am Coll Cardiol* 2000; 36: 959-969.
- Antman EM, Tanasijevic MJ, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. *N Engl J Med* 1996; 335:1342-1349.
- Tanasijevic MJ, Cannon CP, et al. The role of cardiac troponin I in risk stratification of patients with unstable coronary artery disease. *Clin Cardiol* 1999; 22:13-16.
- Antman EM, Fox KM. Guidelines for the diagnosis and management of unstable angina and non-Q-wave myocardial infarction: Proposed revisions. *Am Heart J* 2000; 139: 461-75.
- Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988; 34:27-33.
- Kricka L. Interferences in immunoassays - still a threat. *Clin Chem* 2000; 46: 1037-1038.
- Bjerner J, et al. Immunometric assay interference: incidence and prevention. *Clin Chem* 2002; 48: 613-621.
- CLSI. EP5-A2: Vol. 24, No. 25, Evaluation of Precision Performance of Quantitative Measurement Method; Approved Guideline - Second Edition.
- CLSI. EP06-A: Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline. 2003.

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© 2015-2022 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. Tous droits réservés.
Fabricant : Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Représentant en Europe : Shanghai International Holding Corp. GmbH (Europe)