



AST Reagent Kit (Asparaginic Acid Substrate Method)

Instructions for Use

REF CC1012

PRODUCT NAME

AST Reagent Kit (Asparaginic Acid Substrate Method)

PACKAGE SPECIFICATION

R1:1×20 mL	R2:1×10 mL	R1:1×40 mL	R2:1×20 mL
R1:1×60 mL	R2:1×30 mL	R1:2×60 mL	R2:2×30 mL
R1:2×30 mL	R2:2×15 mL	R1:2×40 mL	R2:2×20 mL
R1:2×50 mL	R2:1×50 mL	R1:2×60 mL	R2:1×60 mL
R1:2×60 mL	R2:3×20 mL	R1:2×65 mL	R2:1×70 mL
R1:2×80 mL	R2:1×80 mL	R1:2×120 mL	R2:2×60 mL
R1:3×20 mL	R2:3×10 mL	R1:3×40 mL	R2:3×20 mL
R1:4×40 mL	R2:4×20 mL	R1:4×50 mL	R2:2×50 mL
R1:4×55 mL	R2:2×55 mL	R1:4×60 mL	R2:2×60 mL
R1:4×60 mL	R2:4×30 mL	R1:4×60 mL	R2:6×20 mL
R1:4×65 mL	R2:2×65 mL	R1:4×65 mL	R2:2×70 mL
R1:4×100 mL	R2:2×100 mL	R1:4×653 mL	R2:4×283 mL
R1:2×65 mL	R2:1×65 mL	R1:4×65 mL	R2:2×65 mL
12×72 T (R1:12×16.8 mL R2:12×8.4 mL)			

INTENDED USE

Used for the *in vitro* quantitative determination of aspartate aminotransferase in human serum and plasma. Mainly used clinically to assist in the auxiliary diagnosis of viral hepatitis, obstructive jaundice and myocardial infarction.

For professional and laboratory use only.

TEST PRINCIPLE

$$\begin{array}{l} \text{Aspartic acid} + \alpha\text{-ketoglutaric acid} \xrightarrow{\text{AST}} \text{Glutamic acid} + \text{Oxaloacetic acid} \\ \text{Oxaloacetic acid} + \text{NADH} + \text{H}^+ \xrightarrow{\text{AST}} \text{Malic acid} + \text{NAD} \end{array}$$
 The rate of NADH reduction was measured at a wavelength of 340 nm and was used to calculate AST viability.

MAIN COMPONENTS

Kit composition	Reagent components	Content
Reagent 1	Tris(hydroxymethyl)aminomethane buffer (pH 7.6)	12.114 g/L
	α -ketoglutaric	6 g/L
	Reduced Coenzyme I	1 g/L
	Sodium azide	2 g/L
	Lactate dehydrogenase	1500 U/L
Reagent 2	Malate dehydrogenase	2000 U/L
	L-aspartic acid	50 g/L
	Sodium azide	2 g/L

The components in different batches of a multi-component kit are not interchangeable.

STORAGE AND SHELF LIFE

Unopened reagents should be stored at 2°C-8°C away from light, with a shelf life of 18 months. Opened reagents are stable for 42 days when stored at 2°C-8°C.

Please refer to the label on the reagent kit for the production date and expiration date.

APPLICABLE INSTRUMENTS

The kit is applicable to the following instruments: fully automatic biochemistry analyzers from Hitachi High-Tech (Shanghai) International Trading Co., Ltd., models: 7100, 7170, 7180, 7600, LABOSPECT 008 AS, 3100, 3500; fully automatic biochemistry analyzers from Beckman Coulter Commercial Enterprise (China) Co., Ltd., models: DXC800, AU480, AU680, AU5800; fully automatic biochemistry analyzers from Canon Medical Systems (China) Co., Ltd., models: TBA-120FR, TBA-2000FR, TBA-FX8; fully automatic biochemistry analyzers from Shenzhen Mindray Bio-Medical Electronics Co., Ltd., models: BS-420, BS-490, BS-600, BS-800, BS-820, BS-2000; fully automatic biochemistry analyzers from Dirui Industrial Co., Ltd., models: CS-400, CS-600B, CS-1200; fully automatic biochemistry analyzers from Siemens Healthineers Diagnostics (Shanghai) Co. Ltd., models: 1800, 2400, ADVIA Chemistry XPT; fully automatic biochemistry analyzers from Roche Diagnostics (shanghai) Co., Ltd., models: cobas 6000 c 501, cobas 8000 c 502, 701, 702; clinical chemistry analyzers from Getein Biotech, Inc, models: CM-400, CM-430, CM-480, CM-600, CM-630, CM-680, CM-800, CM-830, CM-880, CM-2000, CM-1600, CM-1200, CM-1000; automatic biochemical analyzers from Changchun Blaser Medical Technology Co., LTD, models: BBA-400, BBA-300, BBA-480.

If you need the application parameters of the fully automatic biochemistry analyzers, please contact our company.

SAMPLE REQUIREMENTS

Serum and heparinized plasma samples should be separated promptly after collection to avoid hemolysis. Serum and plasma samples can maintain stable enzyme activity for 4 days at 15-25 °C, for 7 days at 2-8 °C and for 3 months at -20 °C.

TEST PROCEDURE

1. Reagent preparation: Use directly.
2. Test conditions: (Different load parameters can be requested based on different testing instruments)

Primary/secondary wavelength	340 nm/415 nm	Calibration type	Linearity
Sample/R1/R2	15-30/200/100 μ L	Time of mixture of serum + R1	3 min
Method	Rate method	Reaction time after addition of R2	3 min
Calibration method	Two-point calibration	Direction of reaction	Downward

Operating procedures:

Dual Reagent Operation

Substances added	Blank tubes	Test tubes
Reagent 1	200 μ L	200 μ L
Distilled water	15-30 μ L	-
Sample	-	15-30 μ L
Mix well, incubate at 37°C for 3 min		
Reagent 2	100 μ L	100 μ L
Mix well, incubate at 37°C for 60-90s, continuously monitor the absorbance change at the measurement wavelength for 1-3min, and calculate Δ A/min.		

3. Calibration procedure: Randox calibration serum is recommended.
4. Quality control procedure: Select quality control serum from Randox, and its measured value should be within the range of its label claim. If the result deviates from the range, find out the reason according to the steps below:
 - 4.1 Check whether the parameter settings and light source are correct.
 - 4.2 Check whether the cuvettes and sampling probes are clean.
 - 4.3 Check whether water is contaminated, and bacterial growth will cause incorrect results.
 - 4.4 Check reaction temperature.
 - 4.5 Check the expiration date of the kit.

5. Result calculation

$$\text{AST viability (U/L)} = \frac{(\Delta A \text{ test/min} - \Delta A \text{ blank/min}) \times K}{\text{Total reaction volume (mL)} \times 1000}$$

$$K = \frac{\text{Sample volume (mL)} \times \text{millimolar extinction coefficient} \times 1.0}{\text{Note: } 1000 = \text{conversion factor from U/mL to U/L; } 1.0 = \text{cuvette optical diameter; Millimolar extinction coefficient} = 6.22}$$

Note: 1000 = conversion factor from U/mL to U/L; 1.0 = cuvette optical diameter; Millimolar extinction coefficient = 6.22

REFERENCE RANGE

The reference range for adults is 0-45 U/L;

The reference range is for reference only, it is recommended that each laboratory should establish its own reference range.

RESULT INTERPRETATION

Since hemolysis interferes with determination, it should be avoided as much as possible during operation.

LIMITATIONS

There is no interference with measurement when hemoglobin is ≤ 500 mg/dL, ascorbic acid ≤ 50 mg/dL, bilirubin ≤ 40 mg/dL, and triglycerides ≤ 1000 mg/dL.

PERFORMANCE CHARACTERISTICS

1. Appearance

Reagent 1 in the kit is a colorless clear liquid, which may contain a small number of subvisible particles that do not affect determination. Reagent 2 is a colorless or slightly yellow clear liquid, which may contain a small number of subvisible particles that do not affect determination.

2. Reagent blank absorbance

2.1 Reagent blank absorbance

Reagent blank absorbance $A_{340\text{nm}} \geq 1.000$.

2.2 Rate of change in absorbance of reagent blanks

Reagent blank absorbance change rate $|\Delta A_{340\text{nm}}|/\text{min} \leq 0.004$.

3. Accuracy

The relative deviation should not fall outside the range of $\pm 15.0\%$.

4. Linear range

For serum sample testing within the reagent linear range of [5, 1000] (37°C):

a) The linear correlation coefficient (r) should not be less than 0.9900;

b) The deviation from linearity should not fall outside the range of $\pm 10\%$ U/L for testing within the linear range of [5, 100] U/L;

the deviation from linearity should not fall outside the range of $\pm 10\%$ for testing within the linear range of (100, 1000] U/L;

5. Analytical sensitivity

When a sample has a concentration of 95.1 U/L, its absorbance difference should be ≥ -0.050 .

6. Precision

6.1 Within-run precision

Within-run precision should not be more than 5.0%.

6.2 Between-run precision

Between-run precision should not be more than 10.0%.

PRECAUTIONS

1. General precautions

1.1 This product is for *in vitro* diagnosis only.

1.2 For clinical diagnosis, please make a comprehensive judgment based on the measurements, clinical symptoms and other findings.

1.3 Please use this product according to the IFU.

1.4 The results of the kit are only used as a basis for clinical diagnosis of various diseases, and should be considered in conjunction with the patient's symptoms/signs, medical history, other laboratory tests and response to treatment.

1.5 Reagents from different manufacturers used to test the same sample may produce different results, which should be considered in conjunction with clinical findings.

2. Precautions for operation

2.1 Please treat the specimens as dangerous substances that may be infected with HIV, HBV, HCV, etc. Please use disposable gloves to avoid or reduce the associated risk for infection.

2.2 If the reagents get into the eyes or mouth, or come into contact with the skin, rinse them quickly and thoroughly with water, and receive medical treatment from a doctor when necessary.

2.3 Avoid direct sunlight during operation.

3. Precautions for use

3.1 Please store the reagents according to the storage method, and avoid freezing. Please do not use frozen reagents whose quality may change.

3.2 Please do not use expired reagents whose test results may be inaccurate.

3.3 Please avoid adding reagents halfway during a test.

3.4 Please avoid direct sunlight during operation.

3.5 Reagents should not be used if they are turbid or if the absorbance of the reagent blank is less than 1.000.

4. Precautions for waste disposal

Samples, waste liquids, etc. are potentially biologically contaminated. Operators should comply with the SOP for laboratory safety and dispose of waste liquids in accordance with local regulations for medical waste, infectious waste, industrial waste, etc.

5. Other precautions

5.1 On a fully automatic biochemistry analyzer, the linearity range is related to the ratio of the amount of a sample to the amount of a reagent and the time of measurement.

5.2 The amounts of the reagent and sample can be changed proportionally according to the requirements of different instruments.

5.3 Please do not use the reagent bottles for other purposes.

5.4 A result calculated with the k value is not as reliable as that obtained using the SRM (calibrator).
















5.5 Please do not mix reagents in different batches.

REFERENCE

Shang Hong, et al. National Standard Operating Procedure for Clinical Testing (4th Edition). People's Medical Publishing House, 2015: 281.

DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on AST Reagent Kit (Asparaginic Acid Substrate Method) are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European Standard EN ISO 15223-1:2021.

Key to symbols used					
	Manufacturer		Use-by date		Catalogue number
	Date of manufacture		Batch code		Temperature limit
	<i>In vitro</i> diagnostic medical device		Keep away from sunlight		Biological risks
	Consult <i>instructions for use</i> or consult <i>electronic instructions for use</i>		Do not use if package is damaged and consult <i>instructions for use</i>		Authorized representative
	CE mark		This way up		Do not re-use



Getein Biotech, Inc.

Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China

Tel: +86-25-68568508

Fax: +86-25-68568500

E-mail: tech@getein.com.cn

overseas@getein.com.cn

Website: www.getein.com



CMC Medical Devices & Drugs S.L.

Add: C/ Horacio Lengo N° 18, CP 29006, Málaga, Spain.

Tel: +34951214054