



STORAGE AND SHELF LIFE

Unopened reagents should be stored at 2°C-8°C away from light, with a shelf life of 18 months. Opened reagents are stable for 42 days when stored at 2°C-8°C.

Please refer to the label on the reagent kit for the production date and expiration date.

Fe Reagent Kit (Ferrozine Chromogenic Method)

Instructions for Use

REF CC1034

PRODUCT NAME

Fe Reagent Kit (Ferrozine Chromogenic Method)


PACKAGE SPECIFICATION

R1: 1x20 mL	R2: 1x10 mL	R1: 1x40 mL	R2: 1x20 mL
R1: 1x60 mL	R2: 1x30 mL	R1: 2x20 mL	R2: 2x10 mL
R1: 2x30 mL	R2: 2x15 mL	R1: 2x40 mL	R2: 2x20 mL
R1: 2x50 mL	R2: 1x50 mL	R1: 2x60 mL	R2: 1x60 mL
R1: 2x60 mL	R2: 2x30 mL	R1: 2x60 mL	R2: 3x20 mL
R1: 2x65 mL	R2: 1x65 mL	R1: 2x80 mL	R2: 1x80 mL
R1: 2x120 mL	R2: 2x60 mL	R1: 3x20 mL	R2: 3x10 mL
R1: 3x40 mL	R2: 3x20 mL	R1: 3x60 mL	R2: 3x30 mL
R1: 4x40 mL	R2: 4x20 mL	R1: 4x50 mL	R2: 2x50 mL
R1: 4x55 mL	R2: 2x55 mL	R1: 4x60 mL	R2: 2x60 mL
R1: 4x60 mL	R2: 4x30 mL	R1: 4x60 mL	R2: 6x20 mL
R1: 6x63 mL	R2: 6x31 mL	R1: 4x65 mL	R2: 2x65 mL
R1: 4x65 mL	R2: 2x65 mL	R1: 4x100 mL	R2: 2x100 mL
R1: 4x639 mL	R2: 4x294 mL	R1: 4x641 mL	R2: 4x278 mL
R1: 4x653 mL	R2: 4x283 mL	R1: 24x4.2 mL	R2: 24x2.1 mL
1x2000 mL (R1: 1x1800 mL R2: 1x200 mL)			
2x2000 mL (R1: 2x1800 mL R2: 2x200 mL)			
12x72 T (R1: 12x17.2 mL R2: 12x8.6 mL)			
Calibrator (optional): 1x1 mL			

INTENDED USE

This reagent kit is intended for the *in vitro* quantitative determination of ferrum concentration in human serum and plasma. Clinically, it is mainly used to aid in the diagnosis of anemia. For professional and laboratory use only.

TEST PRINCIPLE



Within a certain range, the absorbance (A) is directly proportional to the concentration.

MAIN COMPONENTS

Kit composition	Reagent components	Content
Reagent 1	Hydroxylamine hydrochloride	2.0 g/L
	Thiourea	5.0 g/L
	Sodium acetate	18.0 g/L
Reagent 2	Ferrozine	2.0 g/L
Calibrator (optional)	Ferrous sulfate, aqueous matrix	25-40 μmol/L

The components in different batches of a multi-component kit are not interchangeable.
Calibrator traceability: Traceable to international Standard Reference Material (SRM) 909c.

APPLICABLE INSTRUMENTS

The kit is applicable to the following instruments: fully automatic biochemistry analyzers from Hitachi High-Tech (Shanghai) International Trading Co., Ltd., models: 7100, 7170, 7180, 7600, LABOSPECT 008 AS, 3100, 3500; fully automatic biochemistry analyzers from Beckman Coulter Commercial Enterprise (China) Co., Ltd., models: DXC800, AU480, AU680, AU5800; fully automatic biochemistry analyzers from Canon Medical Systems (China) Co., Ltd., models: TBA-120FR, TBA-2000FR, TBA-FX8; fully automatic biochemistry analyzers from Shenzhen Mindray Bio-Medical Electronics Co., Ltd., models: BS-420, BS-490, BS-600, BS-800, BS-820, BS-2000; fully automatic biochemistry analyzers from Dirui Industrial Co., Ltd., models: CS-400, CS-600B, CS-1200; fully automatic biochemistry analyzers from Siemens Healthineers Diagnostics (Shanghai) Co. Ltd., models: 1800, 2400, ADVIA Chemistry XPT; fully automatic biochemistry analyzers from Roche Diagnostics (shanghai) Co., Ltd., models: cobas 6000 c 501, cobas 8000 c 502, 701, 702; clinical chemistry analyzers from Getein Biotech, Inc, models: CM-400, CM-430, CM-480, CM-600, CM-630, CM-680, CM-800, CM-830, CM-880, CM-2000, CM-1600, CM-1200, CM-1000; automatic biochemical analyzers from Changchun Blaser Medical Technology Co., LTD, models: BBA-400, BBA-300, BBA-480. If you need the application parameters of the fully automatic biochemistry analyzers, please contact our company.

SAMPLE REQUIREMENTS

Fasting is required prior to sample collection. Serum and heparin anticoagulated plasma should be separated as soon as possible after sample collection to avoid hemolysis. Once separated, the sample is stable for 7 days at 15-25°C and for 3 weeks at 2-8°C.

TEST PROCEDURE

1. Reagent preparation: The reagent is ready for use directly
2. Test conditions:

Wavelength	570 nm/700 nm	Calibration Type	Linearity
Sample/R1/R2	30/200/100 μL	Time of Mixture of Serum and R1	1-5 min
Method	Two-point end assay	Reaction time after addition of R2	5 min
Calibration Method	Two-point calibration	Direction	Upward

(Absorbance (A) read by the instrument = $A_{\text{Primary Wavelength}} - A_{\text{Secondary Wavelength}}$)

Operating procedures

Sample (calibrator)	30 μL
R1	200 μL
Mix well, incubate at 37°C for 1-5 min, and read the absorbance A_1	
R2	100 μL
Mix well, incubate at 37°C for 5 min, and read the absorbance A_2 , $\Delta A = A_1 - A_2$	

3. Calibration procedure: A calibrator from Getein is recommended, and a Randox calibration serum can also be used.
4. Quality control procedure: Use the quality control serum from Randox, and its measured value should be within the range of its label claim. If the result deviates from the range, find out the reason according to the steps below:
 - 4.1 Check whether the parameter settings and light source are correct.
 - 4.2 Check whether the cuvettes and sample probes are clean.
 - 4.3 Check whether water is contaminated, and bacterial growth can lead to incorrect results.
 - 4.4 Check reaction temperature.
 - 4.5 Check the expiration date of the kit.
5. Result calculation

$$\text{Fe concentration } (\mu\text{mol/L}) = \text{Fe standard concentration} \times \frac{\Delta A_{\text{tested}}}{\Delta A_{\text{calibrated}}}$$

REFERENCE RANGE

Male: 11-30 μmol/L. Female: 9-27 μmol/L.

The reference range is established based on "National Standard Operating Procedure for Clinical Testing (2nd Edition) (1997)" and is for reference purpose only. It is recommended that each laboratory establish its own reference range.

RESULT INTERPRETATION

The concentration of ferrum in serum is extremely low; therefore, use cuvettes and reagents that are free from ferrum contamination.

LIMITATIONS

There is no interference with the measurement when hemoglobin \leq 150 mg/dL, ascorbic acid \leq 500 mg/dL, bilirubin \leq 200 mg/dL and triglycerides \leq 300 mg/dL.

PERFORMANCE CHARACTERISTICS

1. Appearance
Reagent 1 in the kit is a colorless clear liquid, which may contain a small number of insoluble particles that do not affect determination. Reagent 2 is a slightly yellow clear liquid, which may contain a small number of insoluble particles that do not affect determination.
2. Reagent blank absorbance
Reagent blank absorbance $A_{570nm} \leq 0.200$.
3. Accuracy
The relative deviation of test reference materials should not exceed $\pm 20\%$.
4. Linear range
4.1 Linear correlation coefficient
Linear correlation coefficient (r) should be ≥ 0.990 in the range of [2, 179] $\mu\text{mol/L}$.
4.1 Linear deviation
Within the range of [2, 15] $\mu\text{mol/L}$, the linear absolute deviation should not exceed $\pm 3 \mu\text{mol/L}$; within the range of (15, 179] $\mu\text{mol/L}$, the linear relative deviation should not exceed $\pm 10\%$.
5. Analytical sensitivity
When testing a 35.8 $\mu\text{mol/L}$ sample, the absorbance difference should be ≤ 0.140 .
6. Precision
6.1 Repeatability
The coefficient of variation (CV) should not be greater than 6.0%.
6.2 Between-run precision
Between-run precision should not be greater than 10.0%.

PRECAUTIONS

1. General precautions
1.1 This product is for *in vitro* diagnostic use only.
1.2 For clinical diagnosis, please make a comprehensive judgment based on the measurements, clinical symptoms and other findings.
1.3 Please use this product according to the IFU.
2. Precautions for operation
2.1 Treat the specimens as dangerous materials that may cause infection with HIV, HBV, HCV, etc. Please use disposable gloves to avoid or reduce the associated risk for infection.
2.2 If the reagents get into the eyes or mouth, or touch the skin, rinse them quickly and thoroughly with water, and receive medical treatment from a doctor when necessary.
3. Precautions for use
3.1 Please store the reagents according to the storage method, and avoid freezing. Do not use frozen reagents whose quality may change.
3.2 Do not use expired reagents whose test results may be inaccurate.
3.3 Avoid adding reagents halfway during a test.
3.4 Avoid direct sunlight during operation.
3.5 Do not use the reagent if it displays any signs of turbidity.
4. Precautions for waste disposal
Samples, waste liquids, etc. are potentially biologically hazardous. Operators should comply with the SOP for laboratory safety and dispose of waste liquids in accordance with local regulations for medical waste, infectious waste, industrial waste, etc.
5. Other precautions
5.1 On a fully automatic biochemistry analyzer, the linearity range is related to the ratio of the amount of a sample to the amount of a reagent and the time of measurement.
5.2 The amounts of the reagent and sample can be changed proportionally according to the requirements of different instruments.
5.3 Please do not use the reagent bottles for other purposes.
















- 5.4 A result calculated with the k value is not as reliable as that obtained using the calibration result.
- 5.5 Do not mix reagents in different batches.


REFERENCE


Hong Shang, et al. National Standard Operating Procedure for Clinical Testing (4th Edition). People's Medical Publishing House, 2015: 265-269.

DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on Fe Reagent Kit (Ferrozine Chromogenic Method) are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European Standard EN ISO 15223-1:2021.

Key to symbols used					
	Manufacturer		Use-by date		Catalogue number
	Date of manufacture		Batch code		Temperature limit
	<i>In vitro</i> diagnostic medical device		Keep away from sunlight		Biological risks
	Consult <i>instructions for use</i> or consult electronic <i>instructions for use</i>		Do not use if package is damaged and consult <i>instructions for use</i>		Authorized representative
	CE mark		This way up		Do not re-use

 Getein Biotech, Inc.
Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China
Tel: +86-25-68568508
Fax: +86-25-68568500
E-mail: tech@getein.com.cn
overseas@getein.com.cn
Website: www.getein.com

 CMC Medical Devices & Drugs S.L.
Add: C/ Horacio Lengo N° 18, CP 29006, Málaga, Spain.
Tel: +34951214054