Các model liên quan: BFB4, IP1, USF1

2. BFB4

BFB4 Double Channel Electromyographic Biofeedback

Fisiocomputer's BFB4 is a two-channel surface electromyographic biofeedback, that is an electromedical device designed to detect and present in graphic form superficial potentials transmitted by nerve fibers, that is essential for rehabilitation exercises with visual feedback.

The BFB4 consists of a single electronic unit that includes both a PC-Tablet, where the software is pre-installed, and the patient interface, where the two signal detection cables have to be inserted.

Those are connected to the patient adhesive electrodes by special clips. The power is exclusively supplied by the battery inside the tablet.

From a very detailed analysis of the available literature, we can see that countless scientific articles on the use of biofeedback in orthopedic rehabilitation have been written all over the world.

Given the statistical significance of these works it's clear that today it's possible to define a state of art in the field of rehabilitation protocols with biofeedback, so that even the Nomenclator National rate list includes at least two voices: "Individual motor re-education in patients with injured motor skills...including biofeedback" among the outpatient specialist national services recognized by LEA.

Given that, based on the bibliographic references analyzed during the Clinical Validation, regarding the therapeutic indications, these can be defined as:

exercises assisted by feedback in motor rehabilitation, with additional possibility of comparison in real time between two muscles (for example "agonist" / "antagonist").

Strengths:

- Two separate channels for the detection and presentation in graphic form of the superficial potentials transmitted by nerve fibers
- Motor rehabilitation exercises with real-time comparison between two muscles (agonist / antagonist)
 - Strongly indicated for pelvic floor rehabilitation, with virtual obstacles of variable size
 - Extreme sensitivity ($2\mu\text{V} / \text{div}$) for exercises below threshold
 - Data storage in graphical and numerical form
 - Patient registry and planned exercises
 - Glossary of images, movements and associated muscles for positioning electrodes

Physical Features:

The device consists of an 8" tablet inserted in a protection section, the box containing the interface electronics with cables is positioned behind it.

The BFB4 Fisiocomputer comes complete with the following elements:

- a) Main unit, which includes both the tablet 8" and the interface electronics. The identification plate is put behind the patient interface box. The battery that supplies all the power is contained in the PC Tablet. The micro-USB cable that connects the two units are accessible from the outside.
- b) (2 pieces) Patient-cables, with 2.5mm jack plug for insertion into the interface and at the other side three colored cables (black-red-green) ending with a female clip, for the connection to the electrodes. The three cables come from one box on which the identification plate is put.
- c) (Envelope of 30 pieces) Single-skinned adhesive, circular electrodes (diameter 26mm.) with male clip connection.
- d) Mains adapter for battery recharging. Italian plug (US plug adapter included) and USB Female lateral output. Connection cable to the interface from USB-A Male to micro-USB male. Input 100-240V ~ / 50 - 60Hz / 0,16A. Output 5V = / 1A.
- e) (optional) pelvic floor rehabilitation kit.
- f) (optional) Articulated and adjustable floor rod, complete with clamp to support the apparatus.

Functional Features:

- Exercises with virtual obstacles that constitute a feedback tool for the patient;
- Signal acquired in oscillographic mode;
- Patient Card with different settable values;
- Printing with a device connected via WiFi with preview possibility;
- Sensitivity in $\mu\text{V} / \text{div}$;
- Rest time in Sec. (Time between the end of one cycle and the next);
- Target intensity in μV (obstacle height);
- Target start time from the beginning of the cycle to the obstacle
- Target time in Sec. (Obstacles);
- Decontrol time in Sec. (From the last obstacle to the end of the cycle);
- Number of obstacles (from 0 to 3);
- Channel color (chosen from a table);
- $\mu\text{V} / \text{D}$ alarm (can be inserted / removed).

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3. IP1

IP1 HYPERTHERMIA The Focalized Thermotherapy

When we talk about Hyperthermia, we refer to the possibility to focus on a well defined area inside the body, at a variable depth (from 1 cm to 8 cm), enough energy to carry the relevant part (target volume) at a very specific and controlled high temperature.

Hyperthermia is indicated for pathologies of different nature. The therapeutic indications are as follows: bruising outcomes, myalgia, tendinitis, compartment syndromes, post-traumatic and post-surgical articular stiffness, vertebral pain, arthrosis and calcific periarthritis of the shoulder.

The great peculiarity that makes this thermotherapy unique is the excellent reproducibility of thermal parameters, thanks to the use of sophisticated systems for controlling and implementing of the radiant energy addressed to the target.

All subsystems that form the IP1 have been designed and realized with particular care to the aspects of security, reliability and functionality of the product compared to the simplicity of use.

This thermal induction method was found to be so far the most suitable one for treating and resolving situations of the aforementioned pathologies, that appeared tough and with an unsteady and chronic development; all this makes us understand how hyperthermia is far more effective than many others therapies that rely on heat as a mean to treat the patient.

Strengths:

- High precision in achieving the target temperature
- Extremely fast in achieving the temperature of the target
- Easy to use thanks to the smart mode
- Innovative touch screen display

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Physical Features:

- Body Features: Anti-corrosion tubular frame and steel, galvanized with lacquered aluminum frame;

- External dimensions:

Heights from the ground: 985 mm on the top shelf;

1130 mm above the console;

1540 mm at the center of the arm;

Width and depth to the ground: 710x525 mm;

Weight: 75 Kg;

Description: The apparatus is contained in a four-wheeled metal cabinet (two are hidden inside the body, two under the wooden shelves). The control console is located on the top, facing the operator side (the one that has above the handle, visible in the picture), and contains the control electronics, with a practical 7" TFT screen of the "touch-screen" type, which can be used by the operator for all the operations of both setting and control and verification of the supply. The applicator (antenna with bolus, in the foreground on the picture) is placed on a hinged arm fixed to the chassis rear (to the right of the console watcher). The mains plug, fuse holders and light switch are located below on the right side, between the grid and the wooden shelf visible in the picture.

Functional Features:

Touch Screen Panel:

- Smart mode with 3 parameters to set: target depth, temperature and power;
- Parametric Mode: more settable parameters;

RF Generator / Booster:

- Emission frequency: 433.92 MHz;
- Maximum power output: 100 Watt;
- Power setting: 0-100% step 5%;
- Disadvantage protection: total (infinity ROS);
- Direct power measurement: constant in real time;
- Reflected power measurement: constant in real time;
- security controls: absorbed current sense, operating temperature;

Main Applicator:

- Used technology: teflon microstrip;
- Pairing with the patient: silicon bag with thermoregulated liquid;
- Rated frequency: 433.92 MHz;
- Maximum working power: 200W;
- Transfer efficiency: 90% min. (From RF to useful heat);
- Maximum therapeutic penetration: 80mm;
- Maximum footprint: 186x196mm;
- Actual field (50% S.A.R.): 80x120mm (96cmq);

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USF1

CRYOSONIC

The therapy with cold ultrasound (also called cryoultrasound therapy), is a new therapy, of great effectiveness for the treatment of any kind of inflammation and pain. Together with the anti-inflammatory action of ultrasound, in fact, it is associated the analgesic, anti-edema and anti-pain action of cryotherapy (cold therapy).

How does cryoultrasound therapy work?

The anti-inflammatory, pain-relieving and anti-edema action is expressed due to various factors: the mechanical factor, represented by the "set in vibration" by ultrasounds of the treated tissue cells, thus determining a sort of "micro-massage" in depth; the thermal factor: the rise in temperature caused by ultrasound promotes vascularization and elimination of the molecules responsible for inflammation (toxic metabolites); however, the temperature increase occurs only at the surface level, subtracting the blood to the muscles and thereby obtaining hemostasis, (i.e stop of the leakage of blood), very useful after trauma; on the other hand, the action of cold, (which first causes vasoconstriction immediate, but then vasodilation by "vase paralysis), is essentially represented by an analgesic and relaxant effect; the chemical factor: ultrasounds accelerate the demolition of toxic metabolites responsible for inflammation; also, applying an ultrasound anti-inflammatory ointment, ultrasounds facilitate its penetration in the tissues, (thus increasing its effectiveness), with a technique called sonophoresis. The ultrasounds emission from the head may be continuous or intermittent; similarly, the operator may use the handpiece in a circular motion or simply holding it.

The Fisiocomputer USF1 is a medical device for physical therapy designed to match the delivery of ultrasound therapy at 1 MHz and cryotherapy (up to $-9,9^{\circ}\text{C}$).

The synergy of the two therapeutic techniques eliminates the contraindication of the ultrasound in acute and sub-acute traumas and in acute and sub-acute inflammatory states. The head of the cooled ultrasound, placed in correspondence of the lesion to be treated, causes an initial dermal vasoconstriction followed by a strong vasodilation, while in depth a circulatory decrease (hemostasis process) occurs. The decrease in temperature in the deep tissue also creates an increase in the density of the fabric itself, with an increase in the mechanical effect of ultrasound and with a completely canceled thermal effect. In the lesions and inflammation of the musculoskeletal system of the musculoskeletal system, the healing process is accelerated, with a faster recovery of the movement due to the rupture of the pain-inactivity spasm cycle. Therefore this technique is particularly indicated in sports rehabilitation, where the rapid recovery of an injured athlete is made possible because the USF1 allows to treat immediately both a recent trauma and an acute and subacute inflammatory state.

Strengths:

- Epicyclic function to simulate the circumference of the handpiece on the area to be treated
- Ability to treat multiple patients simultaneously
- 1 MHz ultrasound therapy and cryotherapy up to $-9,9^{\circ}\text{C}$

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Physical Features:

- Dimensions and weight:
 - External dimensions:
 - Ground print, arm excluded, 380 x 450 mm. (width x depth);
 - Height of the console at the top 1140 mm.;
 - Support shelf height 815 mm.;
 - Total weight: 73 Kg.;
 - Head diameter US: 62 mm.
- Materials used:
 - Ground base, upper shelf and console support: steel sheet 30/10;
 - Supporting frame: anticorrosive corner profiles;
 - Control console: 12/10 steel sheet;
 - Side panels: 30/10 anticorrosive sheet;
 - Articulated arm: aluminum profiles and light alloy castings;
 - Treatment head: AISI 304 steel applicator, medical PVC supports;
- Finishes:
 - Support frame: natural color anodization;
 - All other metal parts: RAL9006 fire painting;
 - Other non-metallic components: natural;
- Supply:
 - From urban power grid with IEC plug to rear panel;
 - Supplied with detachable mains cable with plug (Italian 19mm pitch.) 10A + T;
 - Mains supply: 230V ~ / 50Hz, maximum absorption: 660VA;
 - General mains fuses: 2xT6.3A (5x20);

Description:

The Fisiocomputer USF1 is a physical therapy device designed to deliver ultrasound combined with cryotherapy, combining in this way both the intrinsic advantage of directly applied cold to the tissues and that of the much greater tolerability of ultrasounds at low temperatures (therefore with the possibility of greater doses and consequent durations of the treatments much inferior). There is a floor unit (with four articulated wheels, two of them with brakes) on which the raised control console is mounted so as to leave the upper support surface free. On the left is fixed the support-base of the articulated arm (with three joints), which supports the Ultrasound treatment head. All connections to the head itself are contained in a flexible tube (of a special type for robotic applications) that protrudes from the body of the device immediately to the left of the support for the control console. Above the base frame, a shelf with folded up edges acts as a tray and / or shelf, and also contains, on the opposite side to the console of control, a sturdy bridge handle to facilitate the operations of movement / positioning of the apparatus. The power cord socket must be inserted into the IEC plug located on the rear panel (bottom). Next to the mains plug there are both the ignition light switch and the two mains fuse holders. In addition to the (detachable) mains cable and the treatment head cable (not detachable and contained in the aforementioned flexible sheath) there are not other cables and / or other connections (fixed or removable).

Technical Features:

- Emission Type: Continuous and pulsed (10/20/50/100 Hz);
- Emission frequency: 1 MHz;
- Transducer: Piezoelectric ceramics;
- Power density adjustable from 0.1 to 3.0 W / cm²;
- Cryotherapy: Adjustable up to -9.9°C;
- Delivery timer adjustable up to 60 'in steps of 1'.

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