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Semi-Compliant Rapid Exchange (RX) Coronary Dilatation Catheters

Abbott Medical
3200 Lakeside Drive
Santa Clara, CA 95054 USA
TEL: (800) 227-9902
FAX: (800) 601-8874
Outside USA TEL: +1 951-914-4669
Outside USA FAX: +1 951-914-2531

**Abbott Vascular
Netherlands BV**
Argonstraat 1, 6422 PH
Heerlen
The Netherlands
TEL: +31-(0)45-5435700

Abbott Vascular International
Park Lane, Culliganlaan 2B
1831 Diegem, BELGIUM
TEL: + 32 2 714 14 11
FAX: + 32 2 714 14 12

Abbott Medical UK Limited
Elder House
Central Boulevard
Blythe Valley Business Park
Solihull B90 8AJ
United Kingdom

English

TREK™ Coronary Dilatation Catheter RX and MINI TREK™ Coronary Dilatation Catheter RX

Semi-Compliant Rapid Exchange (RX) Coronary Dilatation Catheters

* Denotes statements specifically addressing EU Medical Device Regulation 2017/745 requirements

This Instructions for Use (IFU) can also be found on the Abbott eLabeling website at vascular.eifu.abbott.*

CAUTION

This medical device should be used only by physicians trained in angiography and percutaneous transluminal coronary angioplasty (PTCA) in a facility with appropriate equipment and environmental conditions for performing such procedures.

Refer to the instructions supplied with any interventional devices to be used in conjunction with these semi-compliant RX coronary balloon dilatation catheters, for their intended uses, contraindications, and potential complications.

CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE. OBSERVE ALL WARNINGS AND PRECAUTIONS NOTED THROUGHOUT THESE INSTRUCTIONS. FAILURE TO DO SO MAY RESULT IN COMPLICATIONS.

DESCRIPTION

This Instructions for Use (IFU) is applicable to the TREK™ Coronary Dilatation Catheter and MINI TREK™ Coronary Dilatation Catheter, semi-compliant RX Coronary Dilatation Catheters (CDC) which are medical devices intended for use in percutaneous transluminal coronary angioplasty (PTCA) to treat patients with coronary artery disease. These medical devices have an integrated shaft system and a balloon near the distal tip. The shaft has a combination of single lumen and dual lumen tubing. One lumen is used for inflation of the balloon with contrast medium. The second lumen, in the distal shaft, permits the use of a guide wire to facilitate advancement of the dilatation catheter to and through the stenosis to be dilated. The dilatation catheter is coated with HYDROCOAT™ Hydrophilic Coating, which is activated when wet. The total coating length is approximately 42 cm, which includes the distal catheter length measuring approximately 25 cm. The guide wire lumen is coated with a hydrophobic coating. These coatings help improve lubricity and ensure the device performs as intended.

This device has several markers. The balloon has radiopaque marker(s) to aid in positioning the balloon in the stenosis and is designed to provide an expandable segment of known diameter and length at a specific pressure. The proximal shaft has proximal markers, which aid in gauging dilatation catheter position relative to the guiding catheter tip (marker located closest to the dilatation catheter adapter is for femoral guiding catheters and the other marker is for brachial guiding catheters).

The RX design of these dilatation catheters does not incorporate a lumen for distal dye injections and distal pressure measurements. This device is designed for use with 0.014" (0.36 mm) guide wires and has a maximum rated burst pressure (RBP) of 206 psi (14 atm / 1419 kPa).

The maximum crossing profile for these balloons is as follows:

Balloon Diameter	Balloon Crossing Profile Maximum
1.20 mm	0.031" (0.79 mm)
1.50 mm – 2.00 mm	0.032" (0.81 mm)
2.25 mm – 2.50 mm	0.035" (0.89 mm)
2.75 mm – 3.00 mm	0.036" (0.91 mm)
3.25 mm	0.039" (0.99 mm)
3.50 mm	0.040" (1.02 mm)
3.75 mm	0.041" (1.04 mm)
4.00 mm	0.045" (1.14 mm)
4.50 mm	0.046" (1.17 mm)
5.00 mm	0.048" (1.22 mm)

HOW SUPPLIED

Sterile – Sterilized using ethylene oxide sterilization. Non-pyrogenic. Do not use if the package is open or damaged.

This single-use device cannot be reused on another patient as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and / or resterilization may compromise the integrity of the design and / or materials, leading to contamination due to narrow gaps and / or spaces and diminished safety

and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user.

Storage – Although these devices have been tested for proper function after exposure to excursions of extreme temperatures (55°C to -20°C), the following storage conditions are recommended:

- Keep dry
- Keep away from sunlight

Refer to device label for more information.

Contents – One (1) semi-compliant RX Coronary Dilatation Catheter, one (1) protective sheath, one (1) compliance card.

INTENDED PURPOSE*

The semi-compliant RX Coronary Dilatation Catheter is intended to dilate the stenotic structures of a coronary artery or bypass graft for the purpose of improving myocardial blood flow / perfusion.

INDICATIONS

The semi-compliant RX Coronary Dilatation Catheters are indicated for:

- Balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis, for the purpose of improving myocardial perfusion
- Balloon dilatation of a coronary artery occlusion, for the purpose of restoring coronary flow in patients with ST-segment elevation myocardial infarction
- Balloon dilatation of a stent after implantation (balloon models 2.00 mm – 5.00 mm only)
- Balloon dilatation of *de novo* chronic total coronary occlusions (CTO) (MINI TREK™ RX Coronary Dilatation Catheter only)

CONTRAINDICATION

The semi-compliant RX Coronary Dilatation Catheter is not indicated to be used to treat patients with a coronary artery spasm in the absence of a significant stenosis.

INTENDED CLINICAL BENEFIT*

The semi-compliant RX Coronary Dilatation Catheter (CDC) is a medical device that is intended for use in percutaneous transluminal coronary angioplasty (PTCA) to treat patients with coronary artery disease by providing a clinical benefit of improving myocardial blood flow / perfusion.

LIMITATIONS*

These coronary dilatation catheters are designed to be used on one individual during a single procedure, with a device lifetime normally intended for continuous use for less than 60 minutes.

DEVICE SHELF LIFE

The semi-compliant RX Coronary Dilatation Catheter device shelf life is 3 years from the date of manufacture which is specified on the label as the "Use by date". Refer to symbols glossary in this IFU.

COUNSELING AND PROVIDING INFORMATION TO THE PATIENT*

Physicians should review warnings, precautions, and adverse events for residual risks and possible side effects or hypersensitivity risks that may be relevant when counseling patients about use of this device and should consider the following:

- Discussion of the risk / benefit issues for this particular patient
- Discussion of alteration to current lifestyle immediately following the procedure and over the long term
- Discussion of any necessary follow up

WARNINGS

The distal section of the dilatation catheter is coated with HYDROCOAT™ Hydrophilic Coating, which is activated when wet, and the guide wire lumen is coated with MICROGLIDE™ Hydrophobic Coating. Please refer to **PREPARATIONS FOR USE** section of these instructions for further information on how to prepare and use these devices to ensure they perform as intended. Failure to abide by the warnings in this labeling might result in damage to the device coating, which may necessitate intervention or result in serious adverse events.

For single use only. DO NOT sterilize and / or reuse, as this can compromise device performance and increase the risk of cross contamination due to inappropriate reprocessing.*

Note the "Use-by-date" specified on the package.

PTCA should only be performed at centers where emergency coronary artery bypass graft surgery (CABG) is available.

PTCA in patients who are not acceptable candidates for coronary artery bypass graft surgery requires careful consideration, including possible hemodynamic support during PTCA, as treatment of this patient population carries special risk.

Persons with known history of allergies to any of the components of this device listed below may suffer an allergic reaction to this coronary dilatation catheter. Prior to its use on the patient, the patient should be counseled on the materials contained in the device, and a thorough history of allergies must be discussed. This device contains: silicone coating, polyethylene oxide coating, polyamide, polyether block amide (PEBA), polyethylene and stainless steel.*

Use only the appropriate balloon inflation media. Do not use air or any gaseous medium to inflate the balloon. If gaseous medium is used and balloon rupture occurs there is the potential of causing air embolism and / or vessel injury.

Balloon pressure should not exceed the rated burst pressure (RBP). Use of a pressure-monitoring device is recommended to prevent overpressurization.

If an unexpected loss of pressure is to occur, the user should place the inflation device into a neutral (i.e., zero) pressure state and, using fluoroscopy, slowly remove the balloon catheter from the patient. If, at any time during the withdrawal process, there is resistance noted, the cause of the resistance should be determined. If the resistance involves the guiding catheter / introducer sheath both the guiding catheter / introducer sheath and the balloon catheter are to be removed as a single unit. Once the balloon catheter has been completely removed from the patient determine the root cause of the loss of pressure (e.g., inflation device, balloon catheter, connections, etc.) and replace that component. Do not reuse the balloon catheter if it is determined to be the root cause of the loss of pressure.*

To reduce the potential for vessel damage, the inflated diameter of the balloon should approximate the diameter of the normal or undiseased vessel segment, just proximal and distal to the stenosis.

Do not use, or attempt to straighten, a catheter if the shaft has become bent or kinked; this may result in the shaft breaking. Instead, prepare a new catheter.

Treatment of moderately or heavily calcified lesions is considered to be moderate risk, with increase in the risk of acute closure, vessel trauma, balloon burst, balloon entrapment, and associated complications. If resistance is felt, determine the cause before proceeding. Continuing to advance or retract the catheter while under resistance may result in damage to the vessels and / or damage / separation of the catheter.

The safety and effectiveness of these devices has not been established in patients with lesions located in an unprotected left main coronary artery.

PRECAUTIONS

When the catheter is exposed to the vascular system, it should be manipulated while under high quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding.

Do not torque the catheter more than one (1) full turn.

In the event of catheter damage / separation, retrieval methods (use of additional wires, snares, and / or forceps) may result in additional trauma to the coronary vasculature and / or the vascular access site. Complications may include bleeding, hematoma, or pseudoaneurysm.

To confirm sterility has been maintained, ensure that the package sterile barrier has not been opened or damaged prior to use. Inspect all product and ensure that the device is not damaged.

Care must be taken to properly size the balloon prior to use.

During the procedure, appropriate anticoagulant and coronary vasodilator therapy must be provided to the patient as needed. Anticoagulant therapy should be continued for a period of time as determined by the physician after the procedure.

If the surface of the coronary dilatation catheter becomes dry, wet with heparinized normal saline to reactivate the coating.

Do not reinsert the coronary dilatation catheter into the coil dispenser after procedural use.

The safety and effectiveness of these devices have not been established, or is unknown, in vascular regions other than those specifically indicated.*

The safety and effectiveness of this PTCA balloon catheter for the treatment of in-stent stenosis (ISR) has not been established.

The safety and effectiveness of these devices has not been established in the pediatric population.*

With 4.5 mm and 5.0 mm balloon dilatation catheters, some increased resistance may be noted upon insertion or withdrawal into or out of the guiding catheter. Choosing a larger guiding catheter size may minimize this.

Carcinogenic / mutagenic / reproductively toxic (CMR) / endocrine disrupting chemical (EDC). One or more components of this device contain the following substance defined as CMR 1B in a concentration above 0.1% weight by weight: Cobalt; Chemical Abstracts Service (CAS) No. 7440-48-4; EC No. 231-158-0.*

Current scientific evidence supports that medical devices manufactured with stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

ADVERSE EFFECTS

Possible adverse effects include, but are not limited to, the following:

- Allergic reaction or hypersensitivity to latex, contrast agent, anesthesia, device material, and drug reactions to anticoagulation, or antiplatelet drugs
- Vascular access complications which may require transfusion of vessel repair including:
 - Catheter site reactions
 - Bleeding (ecchymosis, oozing, hematoma, hemorrhage, retroperitoneal hemorrhage)
 - Arteriovenous fistula, pseudoaneurysm, aneurysm, dissection, perforation / rupture
 - Embolism (air, tissue, plaque, thrombotic material, or device)
 - Peripheral nerve injury
 - Peripheral ischemia
- Coronary artery or bypass graft complications which may require additional intervention, including:
 - Total occlusion or abrupt closure
 - Arteriovenous fistula, pseudoaneurysm, aneurysm, dissection, perforation / rupture
 - Embolism (air, tissue, plaque, thrombotic material, or device)
 - Thrombosis
 - Stenosis or restenosis
- Pericardial complications which may require additional intervention such as cardiac tamponade, pericardial effusion
- Cardiac arrhythmias (including conduction disorders, atrial and ventricular arrhythmias)
- Cardiac ischemic conditions (including myocardial ischemia, myocardial infarction [including acute], coronary artery spasm and unstable or stable angina pectoris)
- Stroke / cerebrovascular accident (CVA) and transient ischemic attack (TIA)
- System organ failures:
 - Cardio-respiratory arrest
 - Cardiac failure
 - Cardiopulmonary failure (including pulmonary edema)
 - Renal insufficiency
- Blood cell disorders (including heparin induced thrombocytopenia)
- Hypertension / hypotension
- Infection
- Nausea and vomiting
- Palpitation, dizziness, and syncope
- Chest pain
- Fever
- Pain
- Death

REPORTING OF SERIOUS INCIDENTS

Required statement regarding reporting of serious incidents for users and / or patients located in the European Union: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the European Member State in which the user and / or patient is established.*

For other countries outside of the European Member State, the user and / or patient should report any serious incident that has occurred in relation to the medical device to the manufacturer and the authority having jurisdiction in their locale.

MATERIALS REQUIRED

Single Use, Sterile Items (Do not sterilize or reuse)

- Sterile heparinized normal saline
- Guiding Catheter in the appropriate size and configuration to select the coronary artery (diameter not to be less than the minimum recommended guide diameter; see device label)
- Hemostatic valve(s)
- 60% contrast medium diluted 1:1 with normal saline
- 20 cc Luer-lock syringe (optional)
- Appropriately sized guide wire (diameter not to exceed the maximum guide wire for the dilatation catheter; see device label)
- Guide wire introducer
- Guide wire torque device
- Inflation device

PREPARATION FOR USE

The dilatation catheter is coated with HYDROCOAT™ Hydrophilic Coating, which is activated when wet. The total coating length is approximately 42 cm, which includes the distal catheter length measuring approximately 25 cm. Avoid wiping the device with dry gauze, or using alcohol antiseptic solutions, or other solvents to pre-treat the device, as this may impact the coating performance.

Prior to use examine the device carefully for defects or damage. Examine the dilatation catheter for bends, kinks, or other damage. Do not use if the device is damaged or defective.

Prepare equipment to be used following manufacturer's instructions or standard procedure.

Complete the following steps to prepare the coronary dilatation catheter for use:

- Carefully remove the catheter from its protective tubing for preparation. When using a rapid exchange (RX) catheter, do not bend or kink the hypotube during removal.
- Submerge the balloon in sterile heparinized normal saline during balloon preparation to activate the coating.
- Remove the protective mandrel from the distal tip of the dilatation catheter.
- Flush the semi-compliant RX Coronary Dilatation Catheter:

Using a syringe filled with heparinized normal saline, flush the guide wire lumen from the distal end. Flush solution should be seen coming out of the guide wire exit notch located approximately 25 cm proximal to the balloon. Follow this procedure for subsequent flushing.
- Slide the protective sheath off the balloon by grasping the catheter just proximal to the balloon (at the proximal balloon bond site), and with the other hand, grasp the protective balloon sheath and gently remove distally. If unusual resistance is felt during device mandrel and sheath removal, do not use this device and replace with another. Follow device returns procedure for the unused device.
- Prepare an inflation device with the recommended contrast medium according to the manufacturer's instructions.
- Evacuate air from the balloon segment using the following procedure:
 - Fill a 20 cc syringe or the inflation device with approximately 4 cc of the recommended contrast medium.
 - After attaching the syringe or inflation device to the balloon inflation lumen, orient the dilatation catheter with the distal tip and the balloon pointing in a downward vertical position.
 - Apply negative pressure and aspirate for 15 seconds. Slowly release the pressure to neutral, allowing contrast to fill the shaft of the dilatation catheter.
 - Disconnect the syringe or inflation device from the inflation port of the dilatation catheter.
 - Remove all air from the syringe or inflation device barrel. Reconnect the syringe or inflation device to the inflation port of the dilatation catheter, while maintaining negative pressure on the balloon until air no longer returns to the device.
 - Slowly release the device pressure to neutral.
 - Disconnect the 20 cc syringe (if used) and connect the inflation device to the inflation port of the dilatation catheter without introducing air into the system.

CAUTION: All air must be removed from the balloon and displaced with contrast prior to inserting into the body (repeat steps 7a through 7g, if necessary); otherwise, complications may occur.

DIRECTIONS FOR USE

- Insert a guide wire through the hemostatic valve, following the manufacturer's instructions. Advance the guide wire carefully into and through the guiding catheter. When complete, withdraw the guide wire introducer, if used.
- Attach a torque device to the guide wire, if desired. Under fluoroscopy, advance the guide wire to the desired vessel, then across the stenosis.
- Backload the distal tip of the dilatation catheter onto the guide wire, ensuring that the guide wire exits the notch located approximately 25 cm proximal to the balloon.

Note: When backloading the dilatation catheter onto the guide wire, the dilatation catheter should be supported. In advancing the dilatation catheter into the guiding catheter, one's hand should support the dilatation catheter and firmly grasp the proximal shaft. Shaft diameter differences should be taken into consideration when opening and tightening the hemostatic valve and upon withdrawal of the dilatation catheter.
- Advance the dilatation catheter over the guide wire until it approaches the hemostatic valve. Open the hemostatic valve. Insert the dilatation catheter, while maintaining guide wire position, and tighten the hemostatic valve. To facilitate insertion, the balloon must be fully deflated to negative pressure.
 - Tighten the hemostatic valve to create a seal around the dilatation catheter without inhibiting movement of the dilatation catheter. This will allow continuous recording of proximal coronary artery pressure.

Note: It is important that the hemostatic valve be closed tightly enough to prevent blood leakage around the dilatation catheter shaft, yet not so tight that it restricts the flow of contrast into and out of the balloon or restricts guide wire movement.
 - Advance the dilatation catheter until the appropriate proximal marker aligns with the hemostatic valve hub. This indicates that the dilatation catheter tip has reached the guiding catheter tip.
- Advance the dilatation catheter over the guide wire and into the stenosis or stent. Inflate the balloon to a very low pressure (1 atm, 100 kPa or 15 psi) to confirm that the balloon is correctly positioned.

Note: When using the dual wire technique, a dual hemostatic valve should be used and care taken to avoid entanglement, when introducing, torquing, and removing one or both wires. Guide wires should not be rotated more than 180 degrees in either direction during the dual wire procedure. It is recommended that one wire be completely withdrawn from the patient before removing additional equipment.

Note: Post-deployment stent expansion testing was performed on the bench with the MULTI-LINK VISION™ Coronary Stent System and MULTI-LINK ULTRA™ Coronary Stent System stents. All stents should be deployed in accordance with the manufacturer's indications and instructions for use.

- Inflate the balloon (not to exceed 10 total inflations in a stent or 20 total inflations without a stent) to perform PTCA or post-implant dilatation of a stent per standard procedure. Refer to labeling for rated burst pressure and compliance information. Maintain negative pressure on the balloon between inflations.

- Withdraw the deflated dilatation catheter and guide wire from the guiding catheter through the hemostatic valve. Tighten the hemostatic valve.

Note: After the deflated balloon dilatation catheter is withdrawn, it should be wiped clean with gauze soaked with sterile, heparinized normal saline and stored. Prior to reinsertion, the balloon should be submerged in sterile, heparinized normal saline to reactivate the coating.

EXCHANGE PROCEDURE TECHNIQUE

The semi-compliant RX Coronary Dilatation Catheter has been specifically designed for rapid, single operator balloon exchanges. To perform a dilatation catheter exchange:

- Loosen the hemostatic valve.
- Hold the guide wire and hemostatic valve in one hand, while grasping the balloon shaft in the other hand.
- Maintain guide wire position in the coronary artery by holding the wire stationary and begin pulling the dilatation catheter out of the guiding catheter, while monitoring the wire position under fluoroscopy.
- Withdraw the deflated dilatation catheter until the notch in the guide wire lumen is reached (marker indicates notch). Carefully inch the flexible, distal portion of the dilatation catheter out of the rotating hemostatic valve, while maintaining the guide wire's position across the lesion.
- Slide the distal tip of the dilatation catheter out of the hemostatic valve and tighten onto the guide wire to hold it securely in place. Completely remove the dilatation catheter from the guide wire.
- Prepare the next dilatation catheter to be used, as previously described in the **PREPARATION FOR USE** section.
- Backload another dilatation catheter onto the guide wire, as previously described under the **DIRECTIONS FOR USE** section, **step 3**, and continue the procedure accordingly.

DISPOSAL INSTRUCTIONS*

After use, this device, its accessories and packaging should be appropriately classified for disposal, e.g., biohazard, sharps, non-hazardous waste, etc., and carefully disposed of in compliance with facility procedures and applicable laws and regulations.*

REFERENCES

The physician should consult recent literature on current medical practice on balloon dilatation, such as that published by the American College of Cardiology and the American Heart Association (ACC/AHA).

A Summary of Safety and Clinical Performance (SSCP) for this device is available at <https://ec.europa.eu/tools/eudamed>. Search for the device using the Basic UDI-DI 8717648CDC0001LW. This is the SSCP location after the launch of the European Database on Medical Devices/EUDAMED.*

Reference Abbott website for patent markings: www.abbott.com/patents

[™] Indicates a trademark of the Abbott group of companies.

© 2025 Abbott. All Rights Reserved.

Graphical Symbols for Medical Device Labeling

 Manufacturer	 Guiding catheter
 Authorized representative	 French size
 UK Responsible Person	 Rapid exchange
 Importer	 Consult instructions for use or consult electronic instructions for use
 CAUTION: Federal law restricts this device to sale by or on the order of a physician	 Do not re-use
 Use-by date	 Do not resterilize
 Date of manufacture	 Do not use if package is damaged and consult instructions for use
 Batch code	 Single sterile barrier system
 Catalogue number	 Keep away from sunlight
 Unique device identifier	 Keep dry
 Sterilized using ethylene oxide	 Outer diameter
 Non-pyrogenic	 Contains hazardous substances (CAS registry number)
 Packaging unit	 Open here
 Medical device	 Country of manufacture

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This Instructions for Use (IFU) can also be found on the Abbott eLabeling website at vascular.ellfu.abbott. *

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DESCRIPTION

This Instructions for Use (IFU) is applicable to the TREK™ Coronary Dilatation Catheter and MINI TREK™ Coronary Dilatation Catheter, semi-compliant RX Coronary Dilatation Catheters (CDC) which are medical devices intended for use in percutaneous transluminal coronary angioplasty (PTCA) to treat patients with coronary artery disease. These medical devices have an integrated shaft system and a balloon near the distal tip. The shaft has a combination of single lumen and dual lumen tubing. One lumen is used for inflation of the balloon with contrast medium. The second lumen, in the distal shaft, permits the use of a guide wire to facilitate advancement of the dilatation catheter to and through the stenosis to be dilated. The outside diameter of the balloon and distal 25 cm of these dilatation catheters are coated with HYDROCOAT™ Hydrophilic Coating, which is activated when wet. The guide wire lumen is coated with a hydrophobic coating. These coatings help improve lubricity and ensure the device performs as intended.

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This single-use device cannot be reused on another patient as it is not designed to perform as intended after the first usage. Changes in mechanical, physical, and / or chemical characteristics introduced under conditions of repeated use, cleaning, and / or resterilization may compromise the integrity of the design and / or materials,

leading to contamination due to narrow gaps and / or spaces and diminished safety and / or performance of the device. Absence of original labeling may lead to misuse and eliminate traceability. Absence of original packaging may lead to device damage, loss of sterility, and risk of injury to the patient and / or user.

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- a) Balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis, for the purpose of improving myocardial perfusion
- b) Balloon dilatation of a coronary artery occlusion, for the purpose of restoring coronary flow in patients with ST-segment elevation myocardial infarction
- c) Balloon dilatation of a stent after implantation (balloon models 2.00 mm – 5.00 mm only)
- d) Balloon dilatation of *de novo* chronic total coronary occlusions (CTO) (MINI TREK™ RX Coronary Dilatation Catheter only)

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Physicians should review warnings, precautions, and adverse effects for residual risks and possible side effects or hypersensitivity risks that may be relevant when counseling patients about use of this device and should consider the following:

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WARNINGS

The outside diameter of the balloon and distal 25 cm of these dilatation catheters are coated with HYDROCOAT™ Hydrophilic Coating, which is activated when wet, and the guide wire lumen is coated with MICROGLIDE™ Hydrophobic Coating. Please refer to **PREPARATIONS FOR USE** section of these instructions for further information on how to prepare and use these devices to ensure they perform as intended. Failure to abide by the warnings in this labeling might result in damage to the device coating, which may necessitate intervention or result in serious adverse events.

For single use only. DO NOT resterilize and / or reuse, as this can compromise device performance and increase the risk of cross contamination due to inappropriate reprocessing.*

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Persons with known history of allergies to any of the components of this device listed below may suffer an allergic reaction to this coronary dilatation catheter. Prior to its use on the patient, the patient should be counseled on the materials contained in the device, and a thorough history of allergies must be discussed. This device contains: silicone coating, polyethylene oxide coating, polyamide, polyether block amide (PEBAX), polyethylene and stainless steel.*

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The safety and effectiveness of this PTCA balloon catheter for the treatment of in-stent restenosis (ISR) has not been established.

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- Allergic reaction or hypersensitivity to latex, contrast agent, anesthesia, device material, and drug reactions to anticoagulation, or antiplatelet drugs
- Vascular access complications which may require transfusion of vessel repair including:
 - o Catheter site reactions
 - o Bleeding (ecchymosis, oozing, hematoma, hemorrhage, retroperitoneal hemorrhage)
 - o Arteriovenous fistula, pseudoaneurysm, aneurysm, dissection, perforation / rupture
 - o Embolism (air, tissue, plaque, thrombotic material, or device)
 - o Peripheral nerve injury
 - o Peripheral ischemia
- Coronary artery or bypass graft complications which may require additional intervention, including:
 - o Total occlusion or abrupt closure
 - o Arteriovenous fistula, pseudoaneurysm, aneurysm, dissection, perforation / rupture
 - o Embolism (air, tissue, plaque, thrombotic material, or device)
 - o Thrombosis
 - o Stenosis or restenosis
- Pericardial complications which may require additional intervention such as cardiac tamponade, pericardial effusion
- Cardiac arrhythmias (including conduction disorders, atrial and ventricular arrhythmias)
- Cardiac ischemic conditions (including myocardial ischemia, myocardial infarction [including acute], coronary artery spasm and unstable or stable angina pectoris)
- Stroke / cerebrovascular accident (CVA) and transient ischemic attack (TIA)
- System organ failures:
 - o Cardio-respiratory arrest
 - o Cardiac failure
 - o Cardiopulmonary failure (including pulmonary edema)
 - o Renal insufficiency
- Blood cell disorders (including heparin induced thrombocytopenia)
- Hypertension / hypotension
- Infection
- Nausea and vomiting
- Palpitation, dizziness, and syncope
- Chest pain
- Fever
- Pain
- Death

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MATERIALS REQUIRED

Single Use, Sterile Items (Do not sterilize or reuse)

- Sterile heparinized normal saline
- Guiding Catheter in the appropriate size and configuration to select the coronary artery (diameter not to be less than the minimum recommended guide diameter; see device label)
- Hemostatic valve(s)

- 60% contrast medium diluted 1:1 with normal saline
- 20 cc Luer-lock syringe (optional)
- Appropriately sized guide wire (diameter not to exceed the maximum guide wire for the dilatation catheter; see device label)
- Guide wire introducer
- Guide wire torque device
- Inflation device

PREPARATION FOR USE

This device is coated with a hydrophilic coating at the distal 25 cm of the device. Avoid wiping the device with dry gauze, or using alcohol antiseptic solutions, or other solvents to pre-treat the device, as this may impact the coating performance.

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1. Carefully remove the catheter from its protective tubing for preparation. When using a rapid exchange (RX) catheter, do not bend or kink the hypotube during removal.
2. Submerge the balloon in sterile heparinized normal saline during balloon preparation to activate the coating.
3. Remove the protective mandrel from the distal tip of the dilatation catheter.
4. Flush the semi-compliant RX Coronary Dilatation Catheter:
Using a syringe filled with heparinized normal saline, flush the guide wire lumen from the distal end. Flush solution should be seen coming out of the guide wire exit notch located approximately 25 cm proximal to the balloon. Follow this procedure for subsequent flushing.
5. Slide the protective sheath off the balloon by grasping the catheter just proximal to the balloon (at the proximal balloon bond site), and with the other hand, grasp the protective balloon sheath and gently remove distally. If unusual resistance is felt during device mandrel and sheath removal, do not use this device and replace with another. Follow device returns procedure for the unused device.
6. Prepare an inflation device with the recommended contrast medium according to the manufacturer's instructions.
7. Evacuate air from the balloon segment using the following procedure:
 - a) Fill a 20 cc syringe or the inflation device with approximately 4 cc of the recommended contrast medium.
 - b) After attaching the syringe or inflation device to the balloon inflation lumen, orient the dilatation catheter with the distal tip and the balloon pointing in a downward vertical position.
 - c) Apply negative pressure and aspirate for 15 seconds. Slowly release the pressure to neutral, allowing contrast to fill the shaft of the dilatation catheter.
 - d) Disconnect the syringe or inflation device from the inflation port of the dilatation catheter.
 - e) Remove all air from the syringe or inflation device barrel. Reconnect the syringe or inflation device to the inflation port of the dilatation catheter. Maintain negative pressure on the balloon until air no longer returns to the device.
 - f) Slowly release the device pressure to neutral.
 - g) Disconnect the 20 cc syringe (if used) and connect the inflation device to the inflation port of the dilatation catheter without introducing air into the system.

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DIRECTIONS FOR USE

1. Insert a guide wire through the hemostatic valve, following the manufacturer's instructions. Advance the guide wire carefully into and through the guiding catheter. When complete, withdraw the guide wire introducer, if used.
2. Attach a torque device to the guide wire, if desired. Under fluoroscopy, advance the guide wire to the desired vessel, then across the stenosis.
3. Backload the distal tip of the dilatation catheter onto the guide wire, ensuring that the guide wire exits the notch located approximately 25 cm proximal to the balloon.

Note: When backloading the dilatation catheter onto the guide wire, the dilatation catheter should be supported. In advancing the dilatation catheter into the guiding catheter, one's hand should support the dilatation catheter and firmly grasp the proximal shaft. Shaft diameter differences should be taken into consideration when opening and tightening the hemostatic valve and upon withdrawal of the dilatation catheter.

4. Advance the dilatation catheter over the guide wire until it approaches the hemostatic valve. Open the hemostatic valve. Insert the dilatation catheter, while maintaining guide wire position, and tighten the hemostatic valve. To facilitate insertion, the balloon must be fully deflated to negative pressure.
 - a) Tighten the hemostatic valve to create a seal around the dilatation catheter without inhibiting movement of the dilatation catheter. This will allow continuous recording of proximal coronary artery pressure.
Note: It is important that the hemostatic valve be closed tightly enough to prevent blood leakage around the dilatation catheter shaft, yet not so tight that it restricts the flow of contrast into and out of the balloon or restricts guide wire movement.
 - b) Advance the dilatation catheter until the appropriate proximal marker aligns with the hemostatic valve hub. This indicates that the dilatation catheter tip has reached the guiding catheter tip.
5. Advance the dilatation catheter over the guide wire and into the stenosis or stent. Inflate the balloon to a very low pressure (1 atm, 100 kPa or 15 psi) to confirm that the balloon is correctly positioned.
Note: When using the dual wire technique, a dual hemostatic valve should be used and care taken to avoid entanglement, when introducing, torquing, and removing one or both wires. Guide wires should not be rotated more than 180 degrees in either direction during the dual wire procedure. It is recommended that one wire be completely withdrawn from the patient before removing additional equipment.
Note: Post-deployment stent expansion testing was performed on the bench with the MULTI-LINK VISION™ Coronary Stent System and MULTI-LINK ULTRA™ Coronary Stent System stents. All stents should be deployed in accordance with the manufacturer's indications and instructions for use.
6. Inflate the balloon (not to exceed 10 total inflations in a stent or 20 total inflations without a stent) to perform PTCA or post-implant dilatation of a stent per standard procedure. Refer to labeling for rated burst pressure and compliance information. Maintain negative pressure on the balloon between inflations.
7. Withdraw the deflated dilatation catheter and guide wire from the guiding catheter through the hemostatic valve. Tighten the hemostatic valve.
Note: After the deflated balloon dilatation catheter is withdrawn, it should be wiped clean with gauze soaked with sterile, heparinized normal saline and stored. Prior to reinsertion, the balloon should be submerged in sterile, heparinized normal saline to reactivate the coating.

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1. Loosen the hemostatic valve.
2. Hold the guide wire and hemostatic valve in one hand, while grasping the balloon shaft in the other hand.
3. Maintain guide wire position in the coronary artery by holding the wire stationary and begin pulling the dilatation catheter out of the guiding catheter, while monitoring the wire position under fluoroscopy.
4. Withdraw the deflated dilatation catheter until the notch in the guide wire lumen is reached (marker indicates notch). Carefully inch the flexible, distal portion of the dilatation catheter out of the rotating hemostatic valve, while maintaining the guide wire's position across the lesion.
5. Slide the distal tip of the dilatation catheter out of the hemostatic valve and tighten onto the guide wire to hold it securely in place. Completely remove the dilatation catheter from the guide wire.
6. Prepare the next dilatation catheter to be used, as previously described in the **PREPARATION FOR USE** section.
7. Backload another dilatation catheter onto the guide wire, as previously described under the **DIRECTIONS FOR USE** section, **step 3**, and continue the procedure accordingly.

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Reference Abbott website for patent markings: www.abott.com/patents

™ Indicates a trademark of the Abbott group of companies.

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© 2024 Abbott. All Rights Reserved.