

检验报告单

Testing Report

产品名称 Name		丝绸胶带 Silk tape	规格 Specification	Silk tape 1.25cm x 5m, 1 roll/box, 720 rolls/carton	
CO No.		JD20160928	Mfg. Date	05/12/2016	
LOT No.		2016-1	EXP.Date	04/12/2019	
检验依据 Test standard		Zhejinjie-20150021	检验型式 Test type	出厂检验 Delivery Inspection	
分类 No.	检验项目 Test item	标准条款 standard clause	标准要求 Standards		检验结果 Result
A	性能 Function	4.4	医用胶带的剥离强度应不低于1.0N/cm。 Adhesive power \geq 1.0N/cm		3.85N/cm 合格 PASS
		4.5	医用胶带的持粘性应不大于2.5mm。 The viscous \leq 2.5mm		0mm 合格 PASS
B	外观 Appearance	4.1.1	医用胶带应卷齐、平整，表面洁净。 Medical tape should be clean surface, winding tidy, can not have obvious loose or dislocation.		符合 conform 合格 PASS
		4.1.2	医用胶带的胶浆应涂覆均匀，无脱胶、漏胶和背面渗胶。 Medical tape should be glue even, no degum, leak adhesive and permeability glue on the back.		符合 conform 合格 PASS
		4.1.3	胶带解卷时不得有分层现象。 tape is torn with no stratification phenomenon		符合 conform 合格 PASS
	组成 Compose	4.2	医用胶带以丝绸为基材，表面涂布胶粘剂制成。 Medical tape fabric is silk material, the surface is coating with adhesive.		符合 conform 合格 PASS
尺寸 Size	4.2	宽 \pm 0.2cm	Width: \pm 0.2cm	符合 conform	合格 PASS
		长 \pm 0.5m	Long: \pm 0.5m	符合 conform	合格 PASS
结论 Verdict		合格 PASS			

检验 Test person: Yang Yan

复核 Validation: Hu xiaojiang

检验报告单

Testing Report

产品名称 Name		丝绸胶带 Silk tape	规格 Specification	Silk tape 1.25cm x 5m, 24 rolls/box, 720 rolls/carton	
CO No.		JD20160928	Mfg. Date	05/12/2016	
LOT No.		2016-3	EXP.Date	04/12/2019	
检验依据 Test standard		Zhejinjie-20150021	检验型式 Test type	出厂检验 Delivery Inspection	
分类 No.	检验项目 Test item	标准条款 standard clause	标准要求 Standards	检验结果 Result	
A	性能 Function	4.4	医用胶带的剥离强度应不低于1.0N/cm。 Adhesive power \geq 1.0N/cm	3.85N/cm	合格 PASS
		4.5	医用胶带的持粘性应不大于2.5mm。 The viscous \leq 2.5mm	0mm	合格 PASS
B	外观 Appearance	4.1.1	医用胶带应卷齐、平整，表面洁净。 Medical tape should be clean surface, winding tidy, can not have obvious loose or dislocation.	符合 conform	合格 PASS
		4.1.2	医用胶带的胶浆应涂覆均匀，无脱胶、漏胶和背面渗胶。 Medical tape should be glue even, no degum, leak adhesive and permeability glue on the back.	符合 conform	合格 PASS
		4.1.3	胶带解卷时不得有分层现象。 tape is torn with no stratification phenomenon	符合 conform	合格 PASS
	组成 Compose	4.2	医用胶带以丝绸为基材，表面涂布胶粘剂制成。 Medical tape fabric is silk material, the surface is coating with adhesive.	符合 conform	合格 PASS
尺寸 Size	4.2	宽 \pm 0.2cm	Width: \pm 0.2cm	符合 conform	合格 PASS
		长 \pm 0.5m	Long: \pm 0.5m	符合 conform	合格 PASS
结论 Verdict		合格 PASS			

检验 Test person: Hu xiaojiang

复核 Validation: 

检验报告单

Testing Report

产品名称 Name		丝绸胶带 Silk tape	规格 Specification	Silk tape 2.5cm x 5m, 1 roll/box,360rolls/carton	
CO No.		JD20160928	Mfg. Date	05/12/2016	
LOT No.		2016-2	EXP.Date	04/12/2019	
检验依据 Test standard		Zhejinjie-20150021	检验型式 Test type	出厂检验 Delivery Inspection	
分类 No.	检验项目 Test item	标准条款 standard clause	标准要求 Standards		检验结果 Result
A	性能 Function	4.4	医用胶带的剥离强度应不低于1.0N/cm。 Adhesive power \geq 1.0N/cm		3.85N/cm 合格 PASS
		4.5	医用胶带的持粘性应不大于2.5mm。 The viscous \leq 2.5mm		0mm 合格 PASS
B	外观 Appearance	4.1.1	医用胶带应卷齐、平整，表面洁净。 Medical tape should be clean surface, winding tidy, can not have obvious loose or dislocation.		符合 conform 合格 PASS
		4.1.2	医用胶带的胶浆应涂覆均匀，无脱胶、漏胶和背面渗胶。 Medical tape should be glue even, no degum, leak adhesive and permeability glue on the back.		符合 conform 合格 PASS
		4.1.3	胶带解卷时不得有分层现象。 tape is torn with no stratification phenomenon		符合 conform 合格 PASS
	组成 Compose	4.2	医用胶带以丝绸为基材，表面涂布胶粘剂制成。 Medical tape fabric is silk material, the surface is coating with adhesive.		符合 conform 合格 PASS
	尺寸 Size	4.2	宽 \pm 0.2cm	Width: \pm 0.2cm	conform 合格 PASS
	长 \pm 0.5m		Long: \pm 0.5m	conform 合格 PASS	
结论 Verdict		合格 PASS			

检验 Test person: Hu xiaojiang

复核 Validation:

检验报告单

Testing Report

产品名称 Name		丝绸胶带 Silk tape	规格 Specification	Silk tape 2.5cm x 5m, 12 rolls/box, 360 rolls/carton	
CO No.		JD20160928	Mfg. Date	05/12/2016	
LOT No.		2016-4	EXP.Date	04/12/2019	
检验依据 Test standard		Zhejinjie-20150021	检验型式 Test type	出厂检验 Delivery Inspection	
分类 No.	检验项目 Test item	标准条款 standard clause	标准要求 Standards		检验结果 Result
A	性能 Function	4.4	医用胶带的剥离强度应不低于1.0N/cm。 Adhesive power \geq 1.0N/cm		3.85N/cm 合格 PASS
		4.5	医用胶带的持粘性应不大于2.5mm。 The viscous \leq 2.5mm		0mm 合格 PASS
B	外观 Appearance	4.1.1	医用胶带应卷齐、平整，表面洁净。 Medical tape should be clean surface, winding tidy, can not have obvious loose or dislocation.		conform 合格 PASS
		4.1.2	医用胶带的胶浆应涂覆均匀，无脱胶、漏胶和背面渗胶。 Medical tape should be glue even, no degum, leak adhesive and permeability glue on the back.		conform 合格 PASS
		4.1.3	胶带解卷时不得有分层现象。 tape is torn with no stratification phenomenon		conform 合格 PASS
	组成 Compose	4.2	医用胶带以丝绸为基材，表面涂布胶粘剂制成。 Medical tape fabric is silk material, the surface is coating with adhesive.		conform 合格 PASS
	尺寸 Size	4.2	宽 \pm 0.2cm	Width: \pm 0.2cm	conform 合格 PASS
	长 \pm 0.5m		Long: \pm 0.5m	conform 合格 PASS	
结论 Verdict		合格 PASS			

检验 Test person: Hu xiaojiang

复核 Validation: 