

DỰ THẢO NỘI DUNG NHÃN PHỤ





Smith & Nephew Rapid Rhino Nasal Catheter Ant / Posterior Epistaxis
Dụng cụ cầm máu cam khoang mũi trước / sau Rapid Rhino của Smith & Nephew

Số lưu hành: XXXX

Chủ sở hữu số lưu hành: Văn phòng đại diện Smith & Nephew Pte Ltd tại tp HCM

Địa chỉ: Phòng 1335, lầu 13, 39 Lê Duẩn, phường Bến Nghé, Quận 1, tp HCM


Xuất xứ: Costa Rica

Số lô: xem ; Hạn sử dụng xem  trên nhãn gốc

Hướng dẫn sử dụng: xem tờ HDSD đi kèm

Sản phẩm chỉ dùng 1 lần

4 THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF ARTHRO CARE CORP. ANY REPRODUCTION IN PART OR WHOLE WITHOUT THE WRITTEN PERMISSION OF ARTHRO CARE CORP. IS PROHIBITED.



REV		ECN No.	DESCRIPTION	DWG	APPROVED BY	REL DATE
G		55835	PER ECN	AMC	A. CERVENY	23/NOV/15

REF RR 900



LOT 1234567

RAPID RHINO[®]
Nasal Catheter Posterior Epistaxis



2014-01-01

QTY: (1)

Rx Only

STERILE EO



(01)00817470000504(17)140101(10)1234567

NOTES: (UNLESS OTHERWISE SPECIFIED)

1. PRINT LABEL AS SHOWN
2. PRINTING MUST BE CLEAR AND LEGIBLE
3. USE GENERIC LABEL STOCK P/N RR 712


ArthroCare Corporation
7000 West William Cannon Drive
Austin, TX 78735 USA
Phone: (800) 343-5717
www.arthrocare.com

Smith & Nephew
York Science Park
Heslington, York YO10 5DF
United Kingdom
Tel: +44 (0) 1926 482400
ec.rep@smith-nephew.com

P/N 21954-52 Rev. G

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE:		CAD GENERATED DRAWING. DO NOT MANUALLY UPDATE		PARTS LIST	
FRACTIONS	DECIMALS	ANGLES	APPROVALS	DATE	ArthroCare Corporation 7000 West William Cannon Drive Austin, TX 78735
±	.XXx	0.01 ± 2°	DRAWN	10/DEC/13	
		.XXXx	ANALYST	01/JAN/14	LABEL ARTWORK, INNER, RR 900
MATERIAL	SEE NOTES	RESP ENG			
FINISH	SEE NOTES	MFG ENG			SIZE PART TYPE REL LEVEL DWG. NO. REV B 510 2 21954-52 G
DO NOT SCALE DRAWING	QUAL ENG				SCALE: 1:1 CAD FILE: SHEET: 1 OF 1

THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF ARTHROCARE CORP. ANY REPRODUCTION IN PART OR WHOLE WITHOUT THE WRITTEN PERMISSION OF ARTHROCARE CORP. IS PROHIBITED.



REVISONS		DWG	APPROVED BY	REL DATE
REV	ECN No.	DESCRIPTION		
H	60265	PER ECN	AMC	A. CERVENY
				29/JUL/16

REF RR 900
LOT 1234567



RAPID RHINO°
Nasal Catheter Posterior Epistaxis



2014-01-01 QTY: (10)
Rx Only STERILE EO



REF RR 900
LOT 1234567

QTY: (10) 2014-01-01

RAPID RHINO°
Nasal Catheter Posterior Epistaxis

ArthroCare Corporation
7000 West William Cannon Drive
Austin, TX 78735 USA
Phone: (800) 343-5717
www.arthrocare.com

* Trademark of Smith & Nephew. This product may be covered by one or more US Patents. See smith-nephew.com/patents for details.

EC REP
Smith & Nephew
York Science Park
Heslington, York YO10 5DF
United Kingdom
Tel: +44 (0) 1926 482400
ec.rep@smith-nephew.com

- Consult instructions for use
- Do not re-use
- Do not re-sterilize
- Do not use if package is damaged
- Authorized Representative in the European Community

- Keep away from sunlight
- Keep dry
- Use-by date
- Sterilized using ethylene oxide
- Manufacturer

Made in Costa Rica

P/N 21952-52 Rev. H



NOTES: (UNLESS OTHERWISE SPECIFIED)

1. PRINT LABEL AS SHOWN
2. PRINTING MUST BE CLEAR AND LEGIBLE
3. USE GENERIC LABEL STOCK P/N 90021

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE:		CAD GENERATED DRAWING, DO NOT MANUALLY UPDATE		PARTS LIST		
FRACTIONS	DECIMALS	APPROVALS	DATE	ArthroCare Corporation 7000 West William Cannon Drive Austin, TX 78735 LABEL ARTWORK, OUTER, RR 900		
±	.XX± 0.01 ± 2"	DRAWN	M. COLE			28/MAR/14
	.XXX± 0.005	ANALYST	F. MORENO			28/MAR/14
MATERIAL	SEE NOTES	RESP ENG		SIZE	PART TYPE	
FINISH	SEE NOTES	MFG ENG		B	510	
DO NOT SCALE DRAWING		QUAL ENG		REL LEVEL	2	
				DWG. NO.	21952-52	
				SCALE: 1:1	CAD FILE:	
					REV H	
					SHEET: 1 OF 1	