

Comprehensive product portfolio



Treatment centers



Instruments



Laser



Hygiene



3D and Panoramic x-rays



Intraoral x-rays



Intraoral sensors



Phosphor plate scanner

FONA
CLOSER TO YOU



FONA, the most exciting new brand in dental

FONA is a fast-growing, global and dynamic brand, continuously developing its comprehensive product portfolio according to international quality standards. Ensuring excellent price/performance ratio, FONA provides a full range of valuable, reliable and user focused equipment for everyday dentistry. With manufacturing sites in three continents and local support through a rapidly expanding international sales and service team that understands your needs, FONA is truly CLOSER TO YOU.

Copyright FONA Dental, s.r.o., All rights reserved, Subject to technical changes and errors in the text, BR0005, v2, ENG, 20130510

FONA
CLOSER TO YOU

Headquarters:
FONA Dental s.r.o.
Ševčenkova 34
SK - 85101 Bratislava
Slovak Republic
info@fonadental.com
www.fonadental.com

Dealer:

DURABLE EFFICIENCY
Rotary Instruments

FONA Rotary Instruments

DURABLE EFFICIENCY

The comprehensive FONA Instrument range provides all you need for everyday patient treatments. Choose from quality high-speed instruments, air motors, contra-angles and straight handpieces. Reliability and durability every day!



ERGONOMIC DESIGN

Ergonomically shaped and light weight.

TORQUE WHEN YOU NEED IT


Our top quality motors provide exactly the power that you need.

TOP QUALITY

Quality and durability
Made in Germany.




FONA 8080 high-speed handpieces



TECHNICAL SPECIFICATIONS	8080 STANDARD HEAD	8080 MINI HEAD
Bearings	ceramic	ceramic
Light voltage [V]	3,6	3,6
Illumination [Lux]	25.000	25.000
Connection	quick coupling	quick coupling
Power [W]	20	14
Idling speed [rpm]	420.000	420.000
Head diameter / height [mm]	11,4 / 14,8	10,2 / 13,9
Weight [g]	57	55
Pressure [bar]	3,0	2,7
Recommended burr size [ISO]	314/315	313

FONA 8080 high-speed handpieces provide maximum power without compromising in comfort and precision. They combine top technical quality Made in Germany with an attractive price.

Top of the line FONA 8080 high-speed handpieces deliver power up to 20 Watts while ceramic bearings ensure quiet, vibration-free operation. Balanced construction, low weight, hand specific contours and clear light support optimal treatment results. Smooth surface is designed for sure-grip operation and ensures excellent hygiene for dentist and patient.




- Ergonomic design contoured to the shape of the hand
- Well balanced and vibration free
- Smooth surface for sure-grip operation and optimal disinfection
- Anti-retracation system reduces cross contamination risk and prolongs service life
- Precise and powerful chuck system for firm burr fixation

FONA 8060 and 6080 high-speed handpieces



Anti-retraction system



Variety of quick coupling (8060) and fixed (6080) connections



Quality Made in Germany

TECHNICAL SPECIFICATIONS	8060 STANDARD HEAD	8060 MINI HEAD	6080 STANDARD HEAD	6080 MINI HEAD
Bearings	Ceramic	Ceramic	Ceramic	Ceramic
Idling speed [rpm]	400.000	420.000	400.000	420.000
Power [W]	17	14	17	14
Head diameter/height [mm]	11,4 / 14,8	10,2 / 13,9	11,4 / 14,8	10,2 / 13,9
Weight [g]	56	54	53	52
Pressure [bar]	2,7	2,7	2,3	2,3
Recommended burr size	314 / 315	313	314 / 315	313

FONA 8060 and 6080 high-speed handpieces adapt to all daily treatment needs in your dental practice.

FONA 8060's quick coupling connection allows you to change instruments fast with no hassle. FONA 6080 allows you to use both Borden and Midwest connections. Experience vibration-free performance with best results!




- Immediate instrument changes with quick coupling
- Reliable quality and durability
- Push-button system for quick burr change
- Changeable spray ring
- Mini head for better accessibility

QUALITY MADE IN GERMANY
Lasting investment, quality you can rely on.

CERAMIC BEARINGS
Reliable and vibration free performance.


GET CONNECTED
Choose from quick coupling or fixed Midwest or Borden connections.



FONA 6060 high-speed handpieces



Push-button system



Choose from Midwest or Borden fixed connections



Light weight

TECHNICAL SPECIFICATIONS	6060 STANDARD HEAD
Bearings	Ceramic
Idling speed [rpm]	400.000
Power [W]	17
Head diameter/height [mm]	11,4 / 14,8
Weight [g]	43
Pressure [bar]	2,3
Recommended burr size	314 / 315

FONA 6060 provides you with very competitively priced quality instruments, designed for reliability and with German made key parts.

Discover the firm, ergonomic grip and great performance of our 6060 turbines. Make 6060 the backbone of your practice's instrumentation.




- Top quality ceramic bearings
- Vibration free operation
- Push-button system
- Less weight yet same performance
- Reliable chuck system for patient safety

FONA Air motors, straight and contra-angle handpieces



Reliable burr lock mechanism



Choose from Midwest or Borden fixed connections



Durable chrome plated body

TECHNICAL SPECIFICATIONS	40E (STRAIGHT)	40EH (CONTRA)
Bearings	Steel	Steel
Transmission Ratio	1:1	1:1
Idling speed [rpm]	40.000	40.000
Weight [g]	80	70

TECHNICAL SPECIFICATIONS	AIR MOTOR (EM OR EB)
Rotation Speed Range [rpm]	21.000 – 26.000
Torque [N/cm]	2,5
Pressure [bar]	2,5 - 3,0
Weight [g]	82

FONA contra-angle and straight handpieces with regular ISO motor connection meet your daily needs.

FONA 40E and 40EH handpieces are bound to convince you with their price/performance ratio and user friendliness. Get connected with our powerful new Air Motor with either Midwest or Borden connection!



- Push-button burr release system
- External cooling and spray
- Ergonomically shaped
- Safe operation
- Stable torque

FIRM GRIP
Ergonomically shaped for perfect handling.

PUSH BUTTON SYSTEM
Quickly change your burrs at the push of a button.

TORQUE WHEN YOU NEED IT
Powerful air motor with external cooling via Borden and Midwest connection.