

PHARMACARE LIMITED

Ladkrabang Industrial Estate
120 Moo 4, Chalongkrung Road,
Lumplatew, Ladkrabang,
Bangkok 10520, Thailand
Tel : +66 (0) 2326 0990 Fax : +66 (0) 2326 0353, +66 (0) 2739 6112
E-mail : contact@pharmacarelimited.com
www.pharmacarelimited.com



QUALITY ASSURANCE DEPARTMENT

CERTIFICATE OF ANALYSIS/

CERTIFICATE OF STERILIZATION

Cert. No. E 0432/2016

PRODUCT NAME : TIGERPLAST® BLUE METAL DETECTABLE PLASTIC STRIP 100'S,
22 mm x 72 mm
PRODUCT CODE : 45141-000-01
PACKING SIZE : 100 STRIPS/BOX
ORDER NO. : CONTRACT NO.PHAR 001/2016
INVOICE NO. : IVE1600111
BATCH NO. : 16033047
STERILIZATION NO. : 3047
MFG. DATE : 20/04/16
EXP. DATE : 20/04/21
QUANTITY : 360 Boxes

16033047	-	-	-
3047	-	-	-
20/04/16	-	-	-
20/04/21	-	-	-
360 Boxes	-	-	-

SPECIFICATIONS

RESULTS

DESCRIPTION : Each dressing consists of :
1. A blue plasticized PVC film (N) spread evenly with transparent acrylic adhesive mass.
2. A white absorbent pad covered with non-adherent polyethylene net and aluminium detectable backing placed centrally on the adhesive surface.
The adhesive side of the dressing is covered with two strips of non-adhesive protector.
Each dressing is wrapped with printed cold - seal paper.
RESULTS : CONFORM

DRESSING SIZE : 72 x 22 mm. RESULTS : CONFORM

PAD SIZE : 30 x 10 mm. RESULTS : CONFORM

STERILIZATION : The above batch has been sterilized by ethylene oxide in accordance with the requirements of the B.P. 88 Appendix XVIII and ISO 11135 (2007).

REMARKS : - Recommended storage condition : store in a cool dry place, away from heat.
- BY "SEA FREIGHT"

Q.C. ANALYST : Suparnida

Q.A. MANAGER : Phuaporn

DATE : 20-06-16

PHARMACARE LIMITED

Ladkrabang Industrial Estate
120 Moo 4, Chalongkrung Road,
Lumplatew, Ladkrabang,
Bangkok 10520, Thailand

Tel: +66 (0) 2326 0990 Fax : +66 (0) 2326 0353, +66 (0) 2739 6112
E-mail : contact@pharmacarelimited.com
www.pharmacarelimited.com

QUALITY ASSURANCE DEPARTMENT

CERTIFICATE OF ANALYSIS/

CERTIFICATE OF STERILIZATION

Cert. No. E 0317/2017



PRODUCT NAME : TIGERPLAST® PREMIUM ELASTIC FABRIC ADHESIVE
DRESSING STRIP 100'S
PRODUCT CODE : 45133-001-00
PACKING SIZE : 100 STRIPS/BOX
ORDER NO. : CONTRACT NO. PHAR 001/2017
INVOICE NO. : IVE1700084
BATCH NO. : 17023189
STERILIZATION NO. : 3189
MFG. DATE : 10/02/17
EXP. DATE : 10/02/2020
QUANTITY : 504 boxes

17023189	-	-	-
3189	-	-	-
10/02/17	-	-	-
10/02/2020	-	-	-
504 boxes	-	-	-

SPECIFICATIONS

RESULTS


DESCRIPTION : Each dressing consists of :
1. A flesh-colour crosswise elastic fabric (Cotton/Polyamide), spread evenly with zinc oxide rubber based adhesive mass .
2. A white absorbent pad covered with non-adherent polyethylene net, placed centrally on the adhesive surface .
The adhesive side of the dressing is covered with two strips of release liner .
Each dressing is wrapped with heat - seal paper .
CONFORM


DRESSING SIZE : 19 x 72 mm. CONFORM

PAD SIZE : 10 x 30 mm. CONFORM

STERILIZATION : The above batch has been sterilized by ethylene oxide in accordance with the requirements of the B.P. 88 Appendix XVIII and ISO 11135 (2007). CONFORM

REMARKS : BY "SEA FREIGHT"

Q.C. ANALYST : 

Q.A. MANAGER :  DATE : 9-5-17

PHARMACARE LIMITED

Ladkrabang Industrial Estate
120 Moo 4, Chalongkrung Road,
Lumplatew, Ladkrabang,
Bangkok 10520, Thailand
Tel: +66 (0) 2326 0990 Fax : +66 (0) 2326 0353, +66 (0) 2739 6112
E-mail : contact@pharmacarelimited.com
www.pharmacarelimited.com



QUALITY ASSURANCE DEPARTMENT

CERTIFICATE OF ANALYSIS/ CERTIFICATE OF STERILIZATION

Cert. No. E 0319/2017

PRODUCT NAME : TIGERPLAST® CLEAR PLASTIC STRIPS
PRODUCT CODE : 45135-001-01
PACKING SIZE : 100 STRIPS/BOX
ORDER NO. : CONTRACT NO.PHAR 001/2017
INVOICE NO. : IVE1700084
BATCH NO. : 17033197
STERILIZATION NO. : 3197
MFG. DATE : 21/03/17
EXP. DATE : 21/03/2022
QUANTITY : 480 Boxes

17033197	-	-	-
3197	-	-	-
21/03/17	-	-	-
21/03/2022	-	-	-
480 Boxes	-	-	-

SPECIFICATIONS

RESULTS


DESCRIPTION : Each dressing consists of :
1. A clear plasticised PE film spread evenly with transparent surgical adhesive mass.
2. A white absorbent pad covered with non-adherent polyethylene net, placed centrally on the adhesive surface.
The adhesive side of the dressing is covered with two strips of release liner.
Each dressing is wrapped with cold - seal paper.
CONFORM


DRESSING SIZE : 19 x 72 mm. CONFORM

PAD SIZE : 12.5 x 25 mm. CONFORM

STERILIZATION : The above batch has been sterilized by ethylene oxide in accordance with the requirements of the B.P. 88 Appendix XVIII and ISO 11135 (2007). CONFORM

REMARKS : BY "SEA FREIGHT"

Q.C. ANALYST : 

Q.A. MANAGER : 

DATE : 9-5-17

PHARMACARE LIMITED

Ladkrabang Industrial Estate
120 Moo 4, Chalongkrung Road,
Lumplatew, Ladkrabang,
Bangkok 10520, Thailand
Tel: +66 (0) 2326 0990 Fax : +66 (0) 2326 0353, +66 (0) 2739 6112
E-mail : contact@pharmacarelimited.com
www.pharmacarelimited.com



QUALITY ASSURANCE DEPARTMENT

CERTIFICATE OF ANALYSIS/

CERTIFICATE OF STERILIZATION

Cert. No. E 0315/2017

PRODUCT NAME : TIGERPLAST® ANIMAL KINGDOM PLASTIC STRIPS
PRODUCT CODE : 45138-001-00
PACKING SIZE : 100 STRIPS/BOX
ORDER NO. : CONTRACT NO. PHAR 001/2017
INVOICE NO. : IVE1700084
BATCH NO. : 17013187
STERILIZATION NO. : 3187
MFG. DATE : 16/02/17
EXP. DATE : 16/02/2022
QUANTITY : 504 Boxes

17013187	-	-	-
3187	-	-	-
16/02/17	-	-	-
16/02/2022	-	-	-
504 Boxes	-	-	-

SPECIFICATIONS

RESULTS


DESCRIPTION : Each dressing consists of :
1. A character printed PVC film spread evenly with transparent surgical adhesive mass and perforated .
2. A white absorbent pad covered with non-adherent polyethylene net, placed centrally on the adhesive surface.
The adhesive side of the dressing is covered with two strips of release liner.
Each dressing is wrapped with cold - seal paper.


DRESSING SIZE : 19 x 56 mm. CONFORM

PAD SIZE : 12.5 x 20 mm. CONFORM

STERILIZATION : The above batch has been sterilized by ethylene oxide in accordance with the requirements of the B.P. 88 Appendix XVIII and ISO 11135 (2007). CONFORM

REMARKS : BY "SEA FREIGHT"

Q.C. ANALYST : 

-Q.A. MANAGER : 

DATE : 9-5-17