



European Communities Council Directive 93/42/EEC

**ANNEX VII
DECLARATION OF CONFORMITY**

The undersigned, Global Surgical Corporation's authorizing person, declares that the products named and listed on this certificate meet the provisions of this council directive which apply to them, and the CE mark may be affixed.

Technical File: 120-000-011

Device Name: A-Series Surgical Operating Microscope System and A-Series Microscopes

Manufacturer: Global Surgical Corporation
3610 Tree Court Industrial Blvd.
St. Louis, MO 63122, USA

Authorized Representative: Emergo Europe
Molenstraat 15
2513 BH The Hague
The Netherlands

Model(s): See attached list

Intended Use: The Surgical Operating Microscope is a system of modular components intended to improve ergonomics, vision, and lighting during dental and medical procedures/examinations. If the options are chosen, the microscope system may also aid in documentation. Components may include, but are not limited to, support system, binoculars, objective lens, light source, beamsplitters, and video/photography systems. There are four support system options to choose from to meet the clinician's needs: a floor stand model, a floor mounted model, a wall mounted model and a ceiling mounted model.

Device Class: Class I

EC certificate: N/A – Self Certified – Annex IX, Rule 12

Sterility: Device is sold Non-Sterile

Measuring: Device does not perform a measuring function.

UMDNS Code: 12539 – Microscopes, Light, Operating

Testing Standards: IEC 60601-1:2005; IEC 62366:2007; IEC 60601-1-6:2010;
IEC 60601-2-57:2011

Name and Title:
Kevin Hicks
VP of Operations

Signature and Date:  12/19/2016

A-Series Microscope Part Numbers and Descriptions

SUPPORT SYSTEMS:

S.NO	PART NUMBER	DESCRIPTION
1.	M A730F	FLOOR STAND, MOBILE, LOCKING CASTERS, A SERIES
2.	M A730FW	FLOOR STAND COUNTERBALANCE WEIGHT, A SERIES
3.	M A730FM	FLOOR MOUNT, TO-THE-FLOOR, SINGLE MOUNTING PLATE, COLUMN, A SERIES
4.	M A730FMT	FLOOR MOUNT, THROUGH-THE-FLOOR, DUAL MOUNTING PLATES, COLUMN, A SERIES
5.	M A730W	WALL MOUNT, HIGH OR LOW, MOUNTING PLATE, 16" STUDS, SHAFT BRACKET, PLATE COVER, A SERIES
6.	M A730W-KAVO	WALL MOUNT SUPPORT SYSTEM – KAVO, A SERIES
7.	M A730-C8	CEILING COLUMN FOR 8 FOOT MOUNTING SURFACE (NOT DROP CEILING HEIGHT), A SERIES
8.	M A730-C9	CEILING COLUMN FOR 9 FOOT MOUNTING SURFACE (NOT DROP CEILING HEIGHT), A SERIES
9.	M A730-C10	CEILING COLUMN FOR 10 FOOT MOUNTING SURFACE (NOT DROP CEILING HEIGHT), A SERIES
10.	M A730C	CEILING MOUNT, MOUNTING PLATE, STANDARD DROP CEILING, A SERIES
11.	M A730CMF	CEILING MOUNT, MOUNTING PLATE, PLATE COVER, FLUSH-TO-CEILING, CONCRETE OR DRYWALL, A SERIES
12.	M A730-HA	HORIZONTAL ARM, SPRING ARM, EXTENSION ARM (FOR USE WITH ALL FLOORSTANDS, FLOOR MOUNTS, AND LOW WALL MOUNT), A SERIES
13.	M A730-45A	45-DEGREE ARM, SPRING ARM, EXTENSION ARM (FOR USE WITH HIGH WALL AND CEILING MOUNT), A SERIES
14.	M A730-717	MCDONALD ARM, 17" SUPPORT EXTENSION ARM FOR MICROSCOPE SYSTEM

MICROSCOPE BODY:

15.	M A1003	3-STEP A-SERIES MICROSCOPE, AXIS SYSTEM, ANGLE ARM, AND COUPLER, A SERIES
16.	M A1003C	3-STEP A-SERIES MICROSCOPE, AXIS SYSTEM, ANGLE ARM, AND COUNTERBALANCE COUPLER, A SERIES
17.	M A1004	4-STEP A-SERIES MICROSCOPE, AXIS SYSTEM, ANGLE ARM, AND COUPLER, A SERIES
18.	M A1004C	4-STEP A-SERIES MICROSCOPE, AXIS SYSTEM, ANGLE ARM, AND COUNTERBALANCE COUPLER, A SERIES
19.	M A1006	6-STEP A-SERIES MICROSCOPE, AXIS SYSTEM, ANGLE ARM, AND COUPLER, A SERIES
20.	M A1006C	6-STEP A-SERIES MICROSCOPE, AXIS SYSTEM, ANGLE ARM, AND COUNTERBALANCE COUPLER, A SERIES

LIGHT SOURCE:

21.	M A794XA-730	XENON LIGHTSOURCE, A SERIES
22.	M 794XAL	LAMP/HEAT SINK FOR ULTRASOURCE II VERSION A #SYC0049
23.	M A801-LED	LED LIGHT SOURCE FOR THE A-SERIES MICROSCOPE, GREEN AND AMBER (UV) FILTER, A SERIES

BEAM SPLITTER:

24.	M A1061-D50	DUAL PORT BEAMSPLITTER, 50/50 AND 50/50 PORTS, DOVE TAIL MOUNT, A SERIES
25.	M A1061-DVA50	DUAL PORT, VIRTUAL BEAMSPLITTER 95/5 AND 50/50 PORTS, DOVE TAIL MOUNT, A SERIES
26.	M A1061-S50	SINGLE PORT, 50/50 BEAMSPLITTER, DOVE TAIL MOUNT, A SERIES
27.	M A1061-SVA	SINGLE PORT, VIRTUAL BEAMSPLITTER, 95/5 PORT, DOVE TAIL MOUNT, A SERIES

OBJECTIVE LENS:

28.	M A1028-200	200MM OBJECTIVE LENS, 20MM FINE FOCUS ADJUSTMENT, A SERIES
29.	M A1028-250	250MM OBJECTIVE LENS, 20MM FINE FOCUS ADJUSTMENT, A SERIES
30.	M A1028-300	300MM OBJECTIVE LENS, 20MM FINE FOCUS ADJUSTMENT, A SERIES
31.	M 1028ML	MULTI-FOCAL LENS, 200-350MM OPTICAL RANGE
32.	M 1028G-150	OBJECTIVE LENS, 150MM FINE FOCUS (20MM RANGE)
33.	M 1028G-200	OBJECTIVE LENS, 200MM, FINE FOCUS (200MM RANGE)
34.	M 1028G-250	OBJECTIVE LENS, 250MM, FINE FOCUS (20MM RANGE)
35.	M 1028G-300	OBJECTIVE LENS, 300MM, FINE FOCUS
36.	M 1030	OBJECTIVE LENS F/250MM
37.	M 1030A	OBJECTIVE LENS F/250MM/LASER
38.	M 1031	OBJECTIVE LENS F/200MM
39.	M 1031A	OBJECTIVE LENS F/200MM/LASER
40.	M 1034	OBJECTIVE LENS F/300MM
41.	M 1034A	OBJECTIVE LENS F/300MM LASER
42.	M 1035	OBJECTIVE LENS F/350MM
43.	M 1035A	OBJECTIVE LENS F/350MM/LASER
44.	M 1036	OBJECTIVE LENS F/400MM
45.	M 1036A	OBJECTIVE LENS F/400MM/LASER
46.	M 1037	OBJECTIVE LENS F/225MM

LASER FILTER:

47.	M 1047LFM	LASER FILTER MODULE
48.	M CUSTOM	CUSTOM LASER FILTER SLIDE

OPTICAL ACCESSORY:

49.	M A1017	CARR BINOCULAR EXTENDER, 45 DEGREES, A SERIES
50.	M A1019	BINOCULAR ROTATION RING, A SERIES
51.	M A1023C	BINOCULAR CO-OBSERVATION TUBE, 3-AXIS ROTATION, A SERIES
52.	M A1028CLST-3	CLEAR SIGHT LENS COVERS, 3-PACK, A SERIES
53.	M A512	DUAL IRIS DIAPHRAGM, A SERIES
54.	M A517HD1080	HD VIDEO CAMERA, 1080P WITH HDMI CABLE, A SERIES
55.	M A528STHDC	CAMERA ADAPTER, HD VIDEO, A SERIES
56.	M A548F	MONITOR MOUNT, FLOORSTAND, A SERIES
57.	M AFSD	SAMPLE LENS COVER FOR 200-350MM OBJECTIVE LENSES, A SERIES
58.	M APK	PAINT KIT, 1/2OZ BOTTLE, BRUSH, HIGH GLOSS WHITE, A SERIES
59.	M 1039G	HIGH FOLDABLE RUBBER EYECUP
60.	M 1039GL	EYECUP, LONG SILICONE SNAP-ON FOR G-SERIES BINOCULARS
61.	M802GSSV	FIBEROPTIC CABLE, SS WITH VINYL SHEATH, 6MM X 30"
62.	M803GSSV	FIBEROPTIC CABLE, SS WITH VINYL SHEATH, 6MM X 141"
63.	M804GSSV	FIBEROPTIC CABLE, SS WITH VINYL SHEATH, 6MM X 89"
64.	M805GSSV	FIBEROPTIC CABLE, SS WITH VINYL SHEATH, 6MM X 105"
65.	M806GSSV	FIBEROPTIC CABLE, SS WITH VINYL SHEATH, 6MM X 123"

ADAPTER:

66.	M A1050A-725	ADAPTER, A SERIES OPTICS TO FIT G SERIES SUPPORT SYSTEM, FOR USE WITH M A801-LED
67.	M A1050A-725X	ADAPTER, A SERIES OPTICS TO FIT G SERIES SUPPORT SYSTEM, FOR USE WITH XENON LIGHTSOURCE
68.	M A1050G-730	ADAPTER, G SERIES OPTICS TO FIT A SERIES SUPPORT SYSTEM, LED ONLY
69.	M A1050G-730X	ADAPTER, G SERIES OPTICS TO FIT A SERIES SUPPORT SYSTEM, XENON ONLY
70.	M A730-GCK	CEILING MOUNT CONVERSION KIT, G TO A SERIES
71.	M 528J45C	ADAPTER, VIDEO CAMERA, 45MM M517K2 3-CHIP
72.	M 528X-A620	CAMERA COUPLER, CANON A620

73.	M 528XSPC	CAMERA RING COUPLER, 37MM THREAD F/ SONY PC100 AND 110
74.	M 528X-43MM	CAMERA RING COUPLER, 43MM THREAD
75.	M 528RF1	RING FLASH ADAPTER
76.	M 528RF2	RING FLASH ADAPTER
77.	M 528J60C	SINGLE VIDEO ADAPTER, F/ M 1061, M 1061VAV, M 1061VA50 AND C-MOUNT READY CAMERAS, F:60MM
78.	M 528ST45C	STRAIGHT MOUNT VIDEO ADAPTER, F:45MM FOR C-MOUNT HD CAMERA
79.	M 528XSLR-EOS	T-MOUNT ADAPTER FOR CANON EF CAMERAS
80.	M 528XSLR-F	T-MOUNT ADAPTER FOR NIKON F CAMERAS
81.	M 528XSLR-NEX	T-MOUNT TO E-MOUNT ADAPTER FOR SONY NEX5 CAMERA
82.	M 528J60E1.0	VIDEO ADAPTER F/ M1061, VAV, VA50 & ENDOSCOPE CAMERA
83.	M 528X	X-MOUNT ADAPTER WITHOUT CAMERA RING COUPLER
84.	M 528XSLR	X-MOUNT SLR ADAPTER WITHOUT CAMERA RING COUPLER

BINOCULARS:

85.	M A1020-10	BINOCULARS, STRAIGHT, 10X, A SERIES
86.	M A1020-10R	BINOCULARS, STRAIGHT, 10X W/RETICLE, A SERIES
87.	M A1020-10RL	BINOCULARS, STRAIGHT, 10X W/RETICLE (L), A SERIES
88.	M A1020-12	BINOCULARS, STRAIGHT, 12.5X, A SERIES
89.	M A1020-12R	BINOCULARS, STRAIGHT, 12.5X W/RETICLE, A SERIES
90.	M A1021-10	BINOCULARS, INCLINED, 10X, A SERIES
91.	M A1021-10R	BINOCULARS, INCLINED, 10X W/RETICLE, A SERIES
92.	M A1021-10RL	BINOCULARS, INCLINED, 10X W/RETICLE (L), A SERIES
93.	M A1021-12	BINOCULARS, INCLINED, 12.5X, A SERIES
94.	M A1021-12R	BINOCULARS, INCLINED, 12.5X W/RETICLE, A SERIES
95.	M A1022-10	BINOCULARS, INCLINABLE, 10X, A SERIES
96.	M A1022-10R	BINOCULARS, INCLINABLE, 10X W/RETICLE, A SERIES
97.	M A1022-10RL	BINOCULARS, INCLINABLE, 10X W/RETICLE (L), A SERIES
98.	M A1022-12	BINOCULARS, INCLINABLE, 12.5X, A SERIES
99.	M A1022-12R	BINOCULARS, INCLINABLE, 12.5X W/RETICLE, A SERIES

ACCESSORY:

100.	M 557-32	WALL MOUNT BOARD 32" FOR METAL STUD WALLS
101.	M 1071	MICROSCOPE DRAPE STANDARD 20/CS #8178
102.	M 1071L	MICROSCOPE DRAPE F/LASER 10/CS #8173
103.	M 1071XL	MICROSCOPE DRAPE EX LONG 10/CS #8171
104.	M 1072	MICROSCOPE DRAPE NARROW 20/CS
105.	M 1074A	SCOPE DRAPE BARRIER BAG NON-STERILE 300/ROLL
106.	M 1089	STEP BLOCK SPACER 25MM
107.	M 1090	STEP BLOCK SPACER 50MM
108.	M 1091	STEP BLOCK SPACER 75MM
109.	M 547FS22	MONITOR, 22" FLAT LCD W/HDMI INPUT
110.	M DCW	DUST COVER FOR MICROSCOPE
111.	M 800PS	POWER SUPPLY FOR 800 SERIES LED LIGHT