

## **Incubation tray Instructions for use**

| ORDER NO.       | PRODUCT                              | FORMAT      |
|-----------------|--------------------------------------|-------------|
| ZD 9892-3044    | Incubation tray (44 channels, black) | 30 psc/box  |
| ZD 9895-20030-1 | Incubation tray (30 channels, black) | 200 pcs/box |

### **Extended use**

Incubation trays, 1x44/1x30 black are intended for the processing of immunoblots of EUROIMMUN (EUROLINE, EUROLINE-WB, Westernblot) that is a periodic replacement of equipment which follows immunoblot method. They must only be used according to the respective instructions for use of the test kits and by healthcare professionals.

### **Warning:**

- The incubation trays are disposable items and must not be reused.
- The product must only be used by healthcare professionals in suitable laboratory rooms.
- Do not use the tray if the packaging is damaged.
- Before using the product, read the instructions for use carefully. Only use the valid version
- The EUROIMMUN products must not be mixed with or replaced by products from other manufacturers.
- Observe Good Laboratory Practice (GLP) and safety guidelines.

**Waste disposal:** collected and handled as the recommendation of local official regulations.