

## **Incubation tray allergy Instructions for use**

ORDER NO.	PRODUCT	FORMAT
ZD 9897-20030-1	Incubation tray (30 channels, black, volume-reduced for allergy)	200 pcs/box
ZD 9897-12044-1	Incubation tray (44 channels, black, volume-reduced for allergy)	120 pcs/box

### **Extended use:**

Incubation trays allergy, 1x44/1x30 black are intended for the processing of immunoblots of EUROIMMUN (EUROLINE, EUROLINE-WB, Westernblot) that is a periodic replacement of equipment which follows immunoblot method. They must only be used according to the respective instructions for use of the test kits and by healthcare professionals.

### **Warning:**

- The incubation trays are disposable items and must not be reused.
- The product must only be used by healthcare professionals in suitable laboratory rooms.
- Do not use the tray if the packaging is damaged.
- Before using the product, read the instructions for use carefully. Only use the valid version
- The EUROIMMUN products must not be mixed with or replaced by products from other manufacturers.
- Observe Good Laboratory Practice (GLP) and safety guidelines.

**Compatible devices:** EUROBlotMaster 44 and EUROBlot One

**Waste disposal:** collected and handled as the recommendation of local official regulations.