



HUBEI GUARDDENT MEDICAL TECHNOLOGY CO., LTD

Instruction For Use Of Dental Surgical Blade

1. Product description

Dental surgical blades are precision cutting instruments specifically designed for oral and maxillofacial surgery, periodontal surgery, and implant surgery.

2. Intended Used

Installed on the handle of a dental surgical knife, used for cutting, flipping, and trimming of dental soft tissue during surgery

3. Type and Specification

Carbon steel, Stainless steel mater. #10, #11, #12, #15, #15C, #18, #20, #21, #22, #23, #24

4. Warning & Precautions

Dental surgical blades are provided sterile.

5. Packaging and Storage

The products are valid for five years. Please use within the validity period of the product. They should be stored in well ventilated, dry and clean place, avoid direct sunlight.

HUBEI GUARDDENT MEDICAL TECHNOLOGY CO..LTD

Add.: Room 24-29th, 2/F, No. 11 Block, No.1 Guiyuan

Temple South Road, Hanyang, District wuhan City, China

