



HUBEI GUARDDENT MEDICAL TECHNOLOGY CO., LTD

Instruction For Use Of Absorbent Paper Points

1. Product description

Absorbent paper points are auxiliary consumables used in dental root canal therapy for drying the root canal and delivering medicaments. Made from high-absorb paper fibers, they are processed to be sterile and free from optical brighteners or additives, ensuring good biocompatibility.

2. Intended Used

A single-use, conical paper instrument used in dentistry, primarily for absorbing liquids, keeping the operative field dry, collecting samples in endodontic therapy and clinical sampling.

3. Type

15#-40#,45#-80.

4. Warning & Precautions

Absorbent paper points are provided non-sterile.

5. Packaging and Storage

The products are valid for five years. Please use within the validity period of the product. They should be stored in well ventilated, dry and clean place, avoid direct sunlight.

HUBEI GUARDDENT MEDICAL TECHNOLOGY CO., LTD
Add.: Room 24-29th, 2/F, No. 11 Block, No. 1 Gufuyan
Temple South Road, Hanyang, District wuhan City, China

