



Softlines &
Hardlines

Report No. : EFUZ26041508-CG-01
Date : April 21, 2026
Page : 1 of 5

OVERALL RATING	
PASS	X
FAIL	-
DATA	-

CÔNG TY TNHH SẢN XUẤT VÀ THƯƠNG MẠI AN KHANG DENTAL LAB
77/41/16-18, tổ 23, KP 35, Phường Tam Hiệp, Tỉnh Đồng Nai, Việt Nam

TEST REPORT

SAMPLING : Sample Not Drawn By Eurofins Consumer Product Testing Vietnam Co. Ltd
SAMPLE DESCRIPTION : Hàm tháo lắp
Color : /
Material Composition : /
BUYER NAME : /
DELIVERY CONDITION : Sample Received In Good Condition
RECEIVED DATE : April 10, 2026
TESTING PERIOD : April 10, 2026 to April 15, 2026



To be Continued

RESULT SUMMARY

The submitted samples were tested under the following requirements requested by the applicant, please kindly find the attached page(s) for details:

No.	TEST REQUESTED	RESULT			REMARKS
		PASS	FAIL	DATA	
1	Phthalates Content	X			

Pass – Meets Performance Standard

Fail – Fails to Meet Performance Standard

Data – See Test Result

*****For further details, please refer to the following page(s)*****

Signed For And On Behalf Of
Eurofins CPT Vietnam Co. Ltd.




Henry Nguyen
General Manager SL & AN

Results obtained refers to the samples, products or materials received in the Laboratory, as described in above related sample description, and tested in conditions shown in present report. Eurofins Consumer Product Testing Vietnam Co., Ltd. ensures that this job has been performed according to our Quality System and complying contract and legal conditions. Reproduction of this document is only valid if it is done completely and under the written permission of Eurofins Consumer Product Testing Vietnam Co., Ltd. If you have any comments or complaints, please response by sending an email to CSANSL@cpt.eurofinsasia.com and referring to this report number.

To be Continued

SAMPLE PHOTO



Sample A: Hàm tháo lắp

484

TYT

CON

TES

NAI

HỒ

To be Continued

COMPONENT LIST

Component No.	Material	Color	Sample Location
Sample A: Hàm tháo lắp			
A1	Plastic	Translucent pink	Translucent pink plastic (Denture)
A2	Plastic	Off-white	Off-white plastic (Teeth)
A3	Plastic	Orange	Orange plastic (Padding)

TEST RESULT

1. Phthalates Content

Test Method: CPSC-CH-C1001-09.4:2018

Tested item(s)	CAS No.	Reporting Limit (mg/kg)	Limit (mg/kg)	Test result(s) (mg/kg)		
				A1	A2	A3
Diethylhexylphthalate (DEHP)	117-81-7	50	1000	ND	ND	ND
Benzylbutylphthalate (BBP)	85-68-7	50	1000	ND	ND	ND
Dibutylphthalate (DBP)	84-74-2	50	500	ND	ND	ND
Diisononylphthalate (DINP)	28553-12-0	50	1000	ND	ND	ND
Diisodecyl phthalate (DIDP)	68515-49-1	50	1000	ND	ND	ND
CONCLUSION				PASS	PASS	PASS

Note:
ND = not detected, the result is less than reporting limit
mg/kg = milligram/kilogram

NOTE:

The test results are considered as conform to specification based on the general consideration simple acceptance as stated in ISO/IEC GUIDE 98-4:2012



END OF THE REPORT